# Alphalink Cem

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Instructions
<p><b>Description</b></p> <p>Alphalink Cem is a dual-curing fixing composite. The material is fluorescent and radiopaque.</p>

**Indications**
**Alphalink Cem is indicated for**

- the adhesive fixation of indirect restorations such as veneers, inlays, onlays, crowns and bridges. The material is suited to fixing restorations made from glass ceramics, zirconium dioxide, composite, zirconia-reinforced composite, high-precious metal, non-precious metal, titanium.
- the adhesive fixation of root pins.
- building up stumps (as core build up material).

##### Contraindication

- The use of **Alphalink Cem** is contraindicated
- if a patient has known allergies against one or several of the product's components
- if it is impossible to clean and dry the treatment site as necessary or if the prescribed technique cannot be used.

##### Working range

Processing and curing times are dependent on the ambient temperature. As soon as **Alphalink Cem** has been removed from the automix cartridge, the following times are to be considered:

	Ambient temperature ca. 21 °C (ca. 69 °F)	Intraoral 37 °C (98.6 °F)
Processing range	3 - 4 min.	ca. 2 min.
Curing time incl. processing time	ca. 7- 9 min.	ca. 3- 4 min.

<b>Light polymerization</b>		
Light intensity	Optical wavelength range	Exposure time per surface
> 500 mW/cm²	350 - 500 nm	20 sec.

**Note:** After excessive material has been used, please polymerize all edges again for 20 sec.

##### Mixing ratio

With the automix cartridge and the mixing tip, **Alphalink Cem** is always available in the optimal mixing ratio of 1:1.

- Composition** **Alphalink Cem**
- Monomer matrix: diurethane dimethacrylate, 1,4-butane dioldimethacrylate
- Total filler content: 75% by weight, anorganic fillers

##### Application advice

- The material should have reached room temperature before pressing, as the material can be harder to press out and mix at refrigerator temperatures.
- Process **Alphalink Cem** immediately after pressing it out of the automix cartridge and set in the restoration quickly.

##### Use

##### Removal of the temporary restoration

Remove the temporary restoration, clean cavity or stump, remove any residue of the cement that has fixed the temporary restoration and rinse well with water. Subsequently, dry with oil-free compressed air.

##### Try-in of the restoration

- Try in the restoration, check fit and contact points.
- Clean the surfaces to be glued thoroughly with alcohol (e. g. isopropanol). Subsequently, dry with oil-free compressed air.
- Please note:** PMMA materials must **NOT** be cleaned with alcohol. Clean such objects e. g. with a disinfectant based upon quaternary ammonium salts.
- Moisture ingress and contamination must be avoided.
- In case of delicate ceramic restorations, the occlusion test should be carried out very carefully, as before the final fixation of the restoration, there is always a danger of fracturing it.
- If necessary, corrections are carried out with fine diamonds at a medium revolution speed. Any surfaces that were ground must be polished again.
- It is essential for the adhesive fixation of a restoration to dry the working area well (preferably with cotton rolls or parotid gland patches). We recommend the use of a kofferdam.

##### A) Fixation of indirect restorations

##### A1) Pretreatment of the natural tooth

- Apply the etching gel (e. g. Capo Etch, Schütz Dental) to the preparation surface. Etch enamel for 30 sec. and/or dentin for 15 sec.
- Afterwards, rinse off etching gel well with water.
- Wet the surface with a mixture of CapoBond and CapoBond Activator (1:1). Alternatively, a different adhesive system can be used. Please adhere to the manufacturer's instructions.
- We recommend to apply the adhesive with a microbrush.
- Blow off water and solvent residue with a gentle airstream.
- The mixture of CapoBond and CapoBond Activator (1:1) is applied again in a second layer and blown off.

##### A2) Pretreatment of the restoration

- Generally follow the instructions of the restoration material's manufacturer.

Pretreatment of	Procedure
<ul style="list-style-type: none"><li>Metals</li></ul>	<ul style="list-style-type: none"><li>Sandblast the inside surfaces of the restoration with 30 - 50 µm aluminium oxide at a maximum pressure of 2 bar</li> <li>Clean with alcohol and dry the restoration</li> <li>Apply SEBOND Smart* evenly with a thin brush, afterwards let it dry for ca. 1 min.</li></ul>
<ul style="list-style-type: none"><li>Zirconium dioxide</li> <li>Aluminium oxide</li></ul>	<ul style="list-style-type: none"><li>Sandblast the inside surfaces of the restoration with 30 - 50 µm aluminium oxide at a maximum pressure of 2 bar</li> <li>Clean with alcohol and dry the restoration</li> <li>Apply SEBOND* Implant evenly with a thin brush, afterwards let it dry for ca. 1 min.</li> <li><b>Note:</b> For an optimal bond, do not clean the surfaces with phosphoric acid!!!!</li></ul>
<ul style="list-style-type: none"><li>Lithium disilicate glass ceramics e. g. IPS e.max Press, IPS e.max CAD</li></ul>	<ul style="list-style-type: none"><li>Etch with 5% hydrofluoric acid for 20 sec. or in accordance to the glass ceramics e. g. manufacturer's instructions.</li> <li>Thoroughly rinse the restoration with water and dry with oil-free IPS e.max CAD compressed air.</li> <li>Apply SEBOND Implant* evenly with a thin brush, afterwards let it dry for ca. 1 min.</li></ul>
<ul style="list-style-type: none"><li>Composite, e. g. Lava Ultimate</li></ul>	<ul style="list-style-type: none"><li>Sandblast the inside surfaces of the restoration with 30 - 50 µm aluminium oxide at a maximum pressure of 2 bar</li> <li>Clean with alcohol and dry the restoration</li> <li>Apply CapoBond* evenly with a thin brush, blow-dry with a gentle airstream, but do <b>not</b> polymerize.</li></ul>
<ul style="list-style-type: none"><li>PMMA materials, e. g. Tizian Zirconia Reinforced Composite (Schütz Dental)</li></ul>	<ul style="list-style-type: none"><li>Sandblast the inside surfaces of the restoration with 30 - 50 µm aluminium oxide at a maximum pressure of 2 bar</li> <li>Clean with a suitable agent, e. g. disinfectant based on quaterny ammonium salts and dry the restorations</li> <li>Apply CapoBond* evenly with a thin brush, blow-dry with a gentle airstream, but do not polymerize.</li> <li><b>Note:</b> Avoid contact with alcohol as it can lead to cracks.</li></ul>

\*Alternatively, use a different suitable product. In this case, please adhere to the manufacturer's instructions.

##### A3) Application of Alphalink Cem

- Set a new disposable automixing tip onto the dual cartridge. Please dispose of the first bit of mixture coming out of the mixing tip, do not use it for fixation.
- Press out **Alphalink Cem** from the automix cartridge and apply the desired amount directly onto the restoration. Spread evenly.
- As the fixing material will cure inside the mixing tip, it can be used as a cap for the automix cartridge until its next use to protect the contents.

##### A4) Inserting the restoration and removal of excess material

- Bring restoration into position with slight pressure. Fix and hold.
- Remove any excess fixing composite immediately with a disposable brush, a cotton pellet, dental floss or an implant scaler. Pay special attention to the immediate removal of excess material in hard-to-reach areas (approximal, gingiva margin).
- In case of composite of glass ceramic restoration, light-polymerize for 20 sec. from each tooth surface after removal of excess material. Chemical curing is finished after 6 min.
- basically also light-polymerize all edges for 20 sec. after removal of excess material.
- To avoid the formation of an oxygen inhibition layer, the margins of the restoration can be covered with a layer of glycerin gel/air blocker immediately afterwards. This layer can be rinsed off with water after the polymerization has been completed.

##### B) Fixation of root pins

- Prepare the root canal as usual.
- Clean the root canal with a 3% hydrogen peroxide solution (H<sub>2</sub>O<sub>2</sub>) or with a 2,5 - 5,25% sodium hypochlorite solution (NaO-Cl).
- Immediately rinse with water and dry with paper tips.
- Select an appropriate root pin and try in.
- Remove the pin and clean with alcohol.
- Apply a mixture of CapoBond and CapoBond Activator (1:1) with a thin microbrush to the root canal onto the occlusal surface of the stump. Brush in for 15 sec. Remove an excess with paper tips.
- Prepare the root pin according to the manufacturer's instructions and wet it with an appropriate amount of **Alphalink Cem**.
- After setting it in, light polymerize for 40 sec. Afterwards, start the restoration.

##### C) Building up a stump

- If the coronal surface has not yet been etched as described under **B**, apply Capo Etch to the dentin and etch for 15 sec.
- Thoroughly rinse off the etching gel with water, remove all acid residue.
- Extract water and carefully dry the surfaces with paper tips.
- Do not blow dry with air to keep the dentin moist.
- We recommend to apply the adhesive with a microbrush.
- Blow off water and solvent residue with a gentle airstream.
- Apply CapoBond without activator and light polymerize for 2 x 20 sec.
- Set a new disposable automixing tip onto the dual cartridge and apply an appropriate amount of **Alphalink Cem** onto the root pin.
- If a stump matrix is used, fill it and set over the pin.
- Remove excess material and light polymerize for 40 sec.
- Chemical curing is complete after 6 min.
- Finish and polish. An impression can be taken immediately afterwards.

##### Warning

Unpolymerized Alphalink Cem paste contains dibenzoyl peroxide which can cause allergic reactions.

##### Contraindications:

In case that irritations occur or if a patient has known allergies against or hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product. Unpolymerized resin may lead to allergic skin reactions. Therefore, the user should consider protective measures.

##### Interactions with other substances

*Phenolic substances* (like e. g. eugenol) inhibit the polymerization. Therefore, do not use any liners that contain such substances (e. g. zinc oxide-eugenol-cements).

*Oxidative disinfectants* (e. g. hydrogen peroxide) can interact with the polymerization system and thus negatively influence the polymerization. Therefore, do not use an oxidative disinfectant on the automix cartridge. The automix cartridge can be disinfected by e.g. wiping down with medical alcohol.

##### Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e.g. allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. In order to avoid a possible pulp reaction, always prepare a liner in cases with exposed dentin (e. g. a compound containing calcium hydroxide).

##### Storage

- Alphalink Cem** must be stored in a cool place. We recommend refrigerator storage at 3 - 9 °C (37 °F - 48 °F). Leave the mixing tip on the automix cartridge to close it.
- Sebond Implant, Sebond Smart and Capo Bond should be stored at 10 - 25 °C (50 °F - 77 °F). Tightly close the bottle immediately after removal of the primer to keep the volatile components from evaporating.

##### Shelf life

- The maximum shelf life is printed on the label on each primary packaging.
- Do not use the materials after the expiration date.

##### Guarantee

The product was developed for use in dentistry. It must be used exactly in accordance with the instructions. The manufacturer will not be held liable for any damages caused by application of the product for purposes other than those described in these instructions or by improper use. Additionally, it is the user's personal responsibility to check whether the product is suitable for the intended purpose, especially if this purpose has not been described in these instructions.

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