NON-STICK red

revolutionary
ergonomic
design

Deutsch

English

Español

Français

Italiano

Bipolare Pinzetten NON-STICK red

Bipolar Forceps NON-STICK red

Pinzas bipolares NON-STICK red

Pinces bipolaires NON-STICK red

Pinze bipolari NON-STICK red

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Instructions for Use

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Mode d'emploi

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CE 0297

REF 90-989-58-30

Revision No. 3

Date of Release: 2015-05

KLS martin
GROUP

Deutsch.....	5
English	18
Español	31
Français	46
Italiano	59

Symbolerklärung

Symbol explanation

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Légende

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Contents

1	Product Liability and Warranty	19
1.1	General Information	19
1.2	Scope of Delivery	19
1.3	Intended Use.....	19
1.4	Warranty	20
1.5	User's Inspection	20
1.6	Hotline.....	20
2	Notices concerning this Document	21
2.1	Symbols Used in this Document	21
3	Use	22
3.1	Storage and Transportation	22
3.1.1	Screen Basket System	22
3.2	Initial Use	23
3.3	Application	23
4	Cleaning, Maintenance, Disinfection and Sterilization	25
4.1	Preparation for Cleaning	25
4.2	Machine Cleaning and Disinfection.....	26
4.2.1	Validated Program Sequence	26
4.2.2	Material Resistance.....	27
4.3	Sterilization	27
4.3.1	Before sterilization	27
4.3.2	Performing the Sterilization	27
4.3.3	Additional Guidelines	29
5	Disposal	29
6	Reprocessing Information acc. to DIN EN ISO 17664	29

1 Product Liability and Warranty

1.1 General Information

We thank you for having decided to buy a KLS Martin product. This product carries the CE-mark, which means that it satisfies the essential requirements laid down in the EC Directive concerning medical devices.

We are the manufacturer of this product:

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · D-78532 Tuttlingen · Germany

Postfach 60 · D-78501 Tuttlingen · Germany

Tel. +49 7461 706-0 · Fax +49 7461 706-193

info@klsmartin.com · www.klsmartin.com

1.2 Scope of Delivery

- Bipolar forceps NON-STICK **r e d**

Item no. 80-982-12-04 to 80-999-24-04

1.3 Intended Use

The KLS Martin bipolar forceps NON-STICK **r e d** are intended for coagulating, dissecting, grasping and manipulating and tissue. A suitable bipolar cable must be used to connect these products to the appropriate output of an electro-surgical unit that uses bipolar coagulation current.

The maximum output voltage of the generator may not exceed 500 Vp. The appropriate output power setting depends on the specific application, the type of forceps used and the selected coagulation current. In addition, however, the relevant instructions provided in the user manual of the electro-surgical unit must be observed and followed as well.

Suitable bipolar cables:

Products with item numbers 80-287-33-04 / 53-04, 80-286-30-04 / 50-04,

80-287-89-04 / 90-04, 80-287-70-04, 80-291-40-04, 80-286-40-04, 80-293-40-04.

Danger of injuries due to insufficient qualification of the operator!

Safe usage of the bipolar forceps NON-STICK **r e d** requires that the operator is familiar with the methods and application of electro-surgery.

Electrosurgical instruments may only be used by qualified individuals specially trained or instructed in their use!

1.4 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer.

Any warranty exceeding the above provisions shall require a contractual form and shall exclude component-related vandalism, software updates and consumables.

Important Notices

The product may only be repaired by Gebrüder Martin or a qualified person or firm expressly authorized by Gebrüder Martin to perform such work.

If the repair is carried out by a person or firm specially authorized by Gebrüder Martin, the operator of the product is required to obtain from the repairer a certificate with details about the nature and scope of the repair work done. This certificate must show the date of the repair and the details of the person or firm carrying out the work and must be signed.

In all cases where a party other than the product manufacturer performed the work, repaired products must be additionally marked with the repairer's ID label.

Improper interventions or alterations performed by third parties during the period of limitation shall void any and all warranty claims. Unauthorized actions performed on the product shall invalidate any liability claims vis-à-vis Gebrüder Martin.

1.5 User's Inspection

Immediately upon receipt, the goods must be checked for completeness and potential damage in transit. Notice of any such damage must be given immediately.

1.6 Hotline

- Should you have any questions on how to handle the product or use it for clinical applications, please do not hesitate to contact the Product Management:
Tel: +49 7461 706-243
Fax: +49 7461 706-190
- Should you have any technical questions, please do not hesitate to turn to our Martin Service Center:
Tel: +49 7461 706-343
Fax: +49 7461 706-203
E-Mail: msc@klsmartin.com
- Should you have any questions concerning maintenance contracts or training courses, please contact our Technical Service Manager:
Tel: +49 7461 706-332
E-mail: service@klsmartin.com

NOTICE

Each package and the product come with a specific lot number (LOT). Please always indicate this number together with the item number in case of complaint.

2 Notices concerning this Document

Non-observance of these Operating Instructions can lead to serious or even lethal patient injury!

Be sure to read, understand and follow the instructions given below!

- Every user is required to read these Instructions completely and follow them carefully.
- In particular, be sure to heed all cautions, warnings and danger notices.
- Keep this document accessible to users at all times.
- In addition, the Instructions for Use relating to the electrosurgical unit must be observed as well. The safety notes there in must always be observed!

2.1 Symbols Used in this Document

Throughout this document, important information (such as general or safety-related notices) is marked with the following symbols and signal words:

Danger of death or serious injury!

Indicates a situation which, if not avoided, **could** result in death or serious injury!

Danger of minor injury!

Indicates a situation which, if not avoided, could result in minor or moderate injury!

NOTICE

Risk of material damage!


Indicates a situation which, if not avoided, could lead to material damage (loss of time, data loss, device / machine failure, etc.)!

3 Use

3.1 Storage and Transportation

- The forceps must be stored in a clean, cool and dry room.
- Protect against mechanical damage!
- Store and transport in safe containers / packing. Use certified sterile packing (e. g. according to EN 868, ISO 11607) for sterilization, transport and storage. We recommend the KLS Martin screen basket system (see table below).
- Handle with extreme care, do not throw or drop!
- Do not expose to direct sunlight.
- In case of return shipments only send clean and disinfected products in sterile packing.

3.1.1 Screen Basket System

Insert the forceps in the screen basket (1) according to their size. For this purpose, use two fixing elements (3) with retaining clip (4).

An optional cover (2) is available for all screen basket sizes.

Designation	Quantity	Item no.
Screen basket (1), 243 x 255 x 53 mm, for forceps with a maximum length of 23 cm	1	55-805-25-01
Cover (2), 243 x 255 mm	1	55-805-28-01
Screen basket (1), 410 x 255 x 53 mm, for forceps with a minimum length of 23 cm		55-805-42-01
Cover (2), 410 x 255 mm		55-805-45-01
Fixing elements		
Clip (4) for fixing element (PU=10 pcs)	1	55-806-02-04
Fixing element (3, PU=6 pcs)	2	55-806-07-04

3.2 Initial Use

▲WARNING

Danger of infection due to non-sterile handling!

KLS Martin bipolar forceps NON-STICK **r e d** are **NON-STERILE** on delivery!

Non-sterile application and improper sterilization can lead to serious health hazards for patients.

Bipolar forceps NON-STICK **r e d** must always be cleaned and sterilized before using them for the first time and before any further use. See chapter 4 "Cleaning, Maintenance, Disinfection and Sterilization", page 25.

3.3 Application

Prior to each use, the forceps must be subjected to a visual inspection and a functional test and must be cleaned and sterilized. Before using the bipolar forceps check if they are clean, in sound condition and if the insulation is intact. Damaged forceps must be withdrawn from service at once!

Similar to silver, the polished noble-metal tips of the NON-STICK **r e d** forceps can tarnish. This discoloration of the working tips is normal and no reason for concern. It has no negative effect on the functioning of the instruments. The original shine can be re-established by wiping the tips with a silver polishing cloth.

The coagulation current can be activated either with the footswitch or automatically (i. e. via the electro-surgical unit automatic activation function (if available)).

▲WARNING

Peril to life or danger of serious physical injury due to prion infection!

Destroy bipolar forceps after potential contact with prions (risk of CJD contamination!) and do not use again under any circumstances!

▲WARNING

Risk of electric shock and / or burns due to the automatic activation function of the electro-surgical device!

The triggering of an automatic system or unintentional triggering of an activation element can cause uncontrolled activation of HF current.

- Bipolar forceps and cables must only be connected to the electro-surgical device when it is deactivated or in standby mode!
-

⚠ WARNING**Risk of burns by placing electrosurgical instruments on the patient!**

An unintentional HF activation when the function “bipolar automatic operation” is active can occur if the electrosurgical device is placed on the patient and comes into contact with tissue.

- In the intervals between use do not place temporarily unused bipolar forceps on patient!
 - Always make sure the forceps do not accidentally come into contact with tissue, as is done in the case of low-resistance materials (e. g. liquids)!
 - Coagulate only if the working tips of the bipolar forceps are visible. Avoid touching other metallic instruments during coagulation!
-

⚠ WARNING**Risk of explosion due to highly flammable gases!**

Materials saturated with oxygen, such as cotton or waste, can catch fire caused by the sparks that occur during the intentional use of the electrosurgical device.

There is the additional risk of an accumulation of flammable liquids under the patient or in bodily cavities.

- Do not use these products in the presence of flammable or explosive substances!
-

⚠ CAUTION**Improper handling and care as well as non-intended use can lead to premature wear and / or pose a risk for patients and users!**

Therefore, please note:

- These forceps are not intended for use with bipolar cutting current or parameter settings used for bipolar vessel sealing!
- The bipolar forceps NON-STICK **r e d** feature high-quality working tips that must be treated with special care and protected from damage.
- To clean the working tips of the NON-STICK **r e d** forceps during an intervention, only a moist swab may be used (e. g. saturated with physiological salt solution)!

No user repairs are permitted! If the functioning of the bipolar forceps is impaired, you can contact the Martin Service Center, see chapter 1.6 “Hotline”, page 20.

4 Cleaning, Maintenance, Disinfection and Sterilization

4.1 Preparation for Cleaning

The operator or the person processing the product is responsible for the result of the cleaning, disinfection and sterilization performed on the instruments used.

Make sure that

- only product-specific and sufficiently validated cleaning, disinfection and sterilization methods are used,
- the devices used (disinfector, sterilizer) are maintained and checked in regular intervals and that validated parameters are used for every cycle,
- national rules and regulations as well as the hygiene regulations of the physician's practice or the hospital are observed (especially regarding the deactivation of prions).

When choosing the cleaning agent, make sure that

- it is suitable for the cleaning of instruments made of metal and synthetic materials,
- a suitable disinfection agent with proven efficacy is used (e. g. VAH / DGHM or FDA-approved or with CE label) and that it is compatible with the cleaning agent used, and
- that the chemicals used are compatible with the instruments (see chapter 4.2.2 "Material Resistance", page 27).

Disconnect bipolar forceps from connecting cable and clean thoroughly under running water or with an aldehyde-free disinfectant directly after use (within a maximum of 2 hours) to avoid that particles and dried secretions stick to the instruments. Blood and tissue residue should not dry on the instruments.

Due to the better results regarding efficacy and repeatability, machine cleaning should be preferred to manual cleaning.

NOTICE

Risk of damage to the forceps!

- Bipolar forceps NON-STICK **r e d** should not be immersed in hydrogen peroxide (H₂O₂)!
- Do not spread the bipolar forceps NON-STICK **r e d**!
- To avoid scratching the surfaces of the tips and the insulation, do not use abrasive cleaning tools such as metal brushes or steel wool.
- In order to avoid bending, do not apply mechanical load to the tips.

4.2 Machine Cleaning and Disinfection

NOTICE

Risk of damage to the bipolar forceps!

When choosing the disinfectant, make sure that

- the disinfectant has a proven efficacy (e. g. DGHM or FDA-approved or CE label according to DIN EN ISO 15883),
- the program used is suitable for the instrument and includes a sufficient number of washing cycles,
- only sterile or germ-free (10 germs/ml max.) and endotoxin-free (0.25 endotoxin units/ml) water (e. g. purified water / highly purified water) is used for the final rinse,
- the air used for drying is filtered, and
- the disinfectant is maintained and checked on a regular basis.

Observe the criteria mentioned in chapter 4.1 "Preparation for Cleaning", page 25, when choosing the cleaning system that is used.

Insert instruments in the disinfectant. Make sure the instruments do not touch.

4.2.1 Validated Program Sequence

Pre-rinsing

Rinsing temperature: 10 ± 2 °C (50 ± 36 °F)

Soaking time: 1 min

Cleaning

Cleaning temperature: 70 ± 2 °C (158 ± 36 °F)

Soaking time: 5 min

Cleaning detergent: deconex ALKA One-x, concentration 0,5 % (5 ml/l)

Post-rinsing

Rinsing temperature: 10 ± 2 °C (50 ± 36 °F)

Soaking time: 1 min

Thermal disinfection

Disinfection temperature: 90 ± 2 °C (194 ± 36 °F)

Soaking time 5 min

Inspect and pack the instruments immediately after removal from the disinfectant, if necessary after additional post-drying at a clean place.

4.2.2 Material Resistance

Verify that the substances listed below are not among the components of the cleaning or disinfection detergent:

- Organic, mineral, and oxidizing acids (minimum admitted pH-value 5.5)
- Strong lye (maximum admitted pH-value 12)
- Organic solvents (e. g. acetone, ether, alcohol, white gas)
- Oxidizing agents (e. g. peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatized, halogenated hydrocarbons

4.3 Sterilization

4.3.1 Before sterilization

- Perform a visual check for intact insulation, cleanliness and potential damage.
- Similar to silver, the polished stainless steel tips of the NON-STICK **r e d** forceps can tarnish. This discoloration of the working tips is normal and no reason for concern. It has no negative effect on the functioning of the instruments. The original shine can be re-established by wiping the tips with a silver polishing cloth.
- The forceps must be cleaned, disinfected, rinsed with distilled water and finally dried carefully.

4.3.2 Performing the Sterilization

Peril to life or danger of serious physical injury due to prion infection!

Destroy forceps after potential contact with prions (risk of CJD contamination!) and do not use again under any circumstances!

Danger of infection due to non-sterile handling!

Improper sterilization and non-sterile application can lead to serious health hazards for patients.

Sterilization must be carried out according to a validated steam sterilization process, according to the applicable norms EN 285 2009, DIN EN 13060 and ISO 17665-1:2006.

NOTICE

Risk of damage to the bipolar forceps!

- Do not sterilize with hot air.
- Do not use the "Sterrad" plasma sterilization method (STERRAD® is a trademark of Johnson & Johnson Inc.).
- Do not use "flash" autoclaving.
- Do not use dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, or plasma sterilization.

- Use certified sterile packing (e. g. according to EN 68, ISO 1607) or sterile containers for the sterilization. We recommend the KLS Martin screen basket system, see section 3.1.1 "Screen Basket System", page 22.
- Sterilization at high temperature or a long sterilization time reduces the lifespan of the instruments.
- Ensure sufficient drying.
- As regards sterilizer handling, loading and parameter setting, please follow the recommendations of the manufacturer of your machine.
- The operator must ensure that the sterility is maintained after the sterilization process.

Only use the sterilization methods listed below. Do not use other sterilization procedures.

Method	Exposure Temperature	Exposure Time	Drying Time
Fractionated vacuum (dynamic air removal)	min. 132 °C (270 °F) max. 138 °C (280 °F)	min. 3 min max. 20 min	min. 30 min
Gravity ²	min. 132 °C (270 °F)	min. 40 min	min. 30 min

² The less efficient gravitation procedure may only be applied if the fractionated vacuum procedure is not available

- Steam sterilizer according to EN 13060 and EN 285
- Validation according to EN SO 7665 (previously EN 54 / ANSI AMI SO 11134), valid IQ / OQ (commissioning) and product specific performance qualification [PQ]
- Make sure not to exceed a sterilization temperature of 138 °C (280 °F); plus tolerance according to EN ISO 17665 [previously EN 554 / ANSI AAMI ISO 11134]

NOTICE: The fundamental suitability of the instruments for effective steam sterilization by the fractionated vacuum procedure and the gravity procedure was demonstrated by an independent accredited test laboratory (Test report 111739-10). Typical conditions of use in hospitals and doctor's practices as well as the specified procedure were considered.

4.3.3 Additional Guidelines

Cleaning, maintenance and sterilization instructions are available in the Red Brochure entitled "Proper Maintenance of Instruments", published by the "Arbeitskreis Instrumentenaufbereitung" (Instrument Processing Working Group) and to be obtained from Gebrüder Martin GmbH & Co. KG on request.

Useful information on instrument processing is also available on our website at:
<http://kls-martin.com/phnet/>

5 Disposal

Frequent reprocessing has only a minor influence on the service life of surgical instruments, which is mainly a function of use-related wear and tear or damage. Please dispose of your instruments in accordance with relevant local regulations, or have them properly recycled, once they have reached the end of their life cycle.

The national regulations for waste disposal must be complied with!

6 Reprocessing Information acc. to DIN EN ISO 17664

	<p>Danger of prion infection and risk of damage of the NON-STICK red bipolar forceps due to improper use!</p> <ul style="list-style-type: none"> • Destroy forceps after potential contact with prions (risk of CJD contamination!) and do not use again under any circumstances! • Do not use the "Sterrad" plasma sterilization method! • Do not sterilize with hot air! • Do not use "flash" autoclaving! • Do not immerse in hydrogen peroxide (H₂O₂)! • To avoid scratching polished working tips and insulation, do not use abrasive cleaning tools!
Reprocessing limitations	Damage such as bent tips, loss of insulation due to cracks or splintered plastic coating and wear.

INSTRUCTIONS

Workplace	Central Sterile Supply Department (CSSD), instrument processing department
Storage and transport	Store forceps in a clean, cool and dry room. Protect against mechanical damage! Use protective caps for forceps with sensitive tips. Handle with extreme care, do not throw or drop. Store and transport in safe containers / packing.
Cleaning preparation	Disconnect forceps from connecting cable and clean thoroughly under running water or with a aldehyde-free disinfectant directly after use (within a maximum of 2 hours) to avoid that particles and dried secretions stick to the instruments.

Cleaning: manual	<p>A validated manual procedure is not indicated, the machine process should be preferred. In case the mechanical procedure is not available, make sure that when choosing the cleaning system used</p> <ul style="list-style-type: none"> • it is suitable for the cleaning of instruments made of metal and synthetic materials, • a suitable disinfection agent with proven efficacy is used (e. g. VAH / DGHM or FDA-approved or with CE label) and that it is compatible with the cleaning agent used, and • that the chemicals used are compatible with the instruments. 			
Cleaning: Mechanical procedure	<p>Validated program sequence: Pre-cleaning: Rinsing temperature: 10 °C ± 2 °C (50 °F±36 °F), insertion time: 1 min Cleaning: Cleaning temperature: 70 °C ± 2 °C (158 °F±36 °F), insertion time: 5 min, cleaning agent: deconex ALKA One-x, concentration 0.5 % (5 ml/l) Final rinse: Rinsing temperature 10 °C ± 2 °C (50 °F±36 °F), insertion time: 1 min</p>			
Disinfection	<p>Thermal disinfection for mechanical treatment 90 °C±2 °C, (194 °F±36 °F) insertion time: 5 min Immersion in disinfectant solution. Observe manufacturer's information on concentration, exposure time and temperature. Observe material compatibility (see chapter 4.2.2, page 27), 90 °C (194 °F), min. 5 minutes.</p>			
Maintenance, checks and function testing	<p>Check forceps for damage prior to each use. Withdraw damaged forceps from service immediately. Similar to silver, the polished stainless steel tips of the NON-STICK red forceps can tarnish. However, this has no negative effect on the functioning of the instruments. The original shine can be re-established by wiping the tips with a silver polishing cloth.</p>			
Packaging	<p>Use certified sterile packing (e. g. according to EN 868, ISO 11607) or sterile containers.</p>			
Sterilization	<p>Prior to sterilization, forceps must be cleaned, disinfected, rinsed with distilled water and finally dried carefully.</p>			
	Method	Exposure Temperature	Exposure Time	Drying Time
	Fractionated vacuum (dynamic air removal)	min. 132 °C (270 °F) max. 138 °C (280 °F)	min. 3 min max. 20 min	min. 30 min
	Gravity	min. 132 °C (270 °F)	min. 40 min	min. 30 min
<p>As regards sterilizer handling, loading and parameter setting, please follow the recommendations of the manufacturer of your machine.</p>				
<p>Last revised 02 / 2012</p>				

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- Non sterilizzare le pinze bipolari NON-STICK **red** con sterilizzatori a gas plasma STERRAD! (STERRAD® is a trademark of Johnson & Johnson Inc.).
Martin Nederland/Marned B.V.
KLS Martin UK Ltd.