## Composite Wetting Resin

### Description and Indications:


Composite Wetting Resin is a 45% filled, radiopaque light cured wetting resin. It may be used during the placement of composite if the composite has become dry during incremental layering, or if the oxygen inhibition layer has been removed or disturbed. Composite Wetting Resin helps to improve the glide of the instrument during sculpting and contouring.

### Procedure:


1. Remove Luer lock caps and attach the Inspiral® Brush tip to syringe.
2. Apply a thin layer to cured composite increment as needed. Air or vacuum thin (vacuum excess resin using 18-20 gauge Endo-Eze® tip on Luer Vacuum adapter).
3. If using Composite Wetting Resin for instrument glide, apply a small quantity to instrument or surface of soft composite.

### Precautions:

1. Any resin can be sensitizing to skin or can cause allergic reaction if individual is already sensitized. Therefore, wash area of skin exposed to resin completely to remove the material. If dermatitis, rash or any other allergic reaction is observed, see a physician.
2. When dispensing intraorally, never force syringe plungers. If there is excess resistance to flow, replace tip with a new one and verify flow before proceeding.
3. Use only recommended tips with each specific material. Tips are disposable.
4. We recommend using Ultradent® Syringe Covers sealed with the Ultradent® Impulse Sealer for a protective barrier to prevent cross-contamination.

**EN** Apply a thin layer to cured composite increment as needed. Air or vacuum thin (vacuum excess resin using 18-20 gauge Endo-Eze® tip on Luer Vacuum adapter).

If using Composite Wetting Resin for instrument glide, apply a small quantity to instrument or surface of soft composite.

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