# Detaseal® bite

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# Detaseal® bite

**GB**  
Instructions for use

**Bite registration material, A-silicone based, addition curing fast setting – high final hardness**


Fig. 1

## 1. Mixing and Dosing

Extrude material from the cartridge with Automix2-system gun (Fig. 1). Place cartridge into mixing gun and remove closure cap by rotation. Eject a small amount of material (bleeding) until equal flow is being extruded from both orifices (Fig. 2). Insert mixing cannula (rose) into guide grooves on cartridge and cannula and rotate clockwise to a stop (Fig. 3). Dispenser is now ready to mix silicone in any amount needed. Squeeze out the material with uniform pressure. Store cartridge with mixing cannula attached for subsequent use.


Fig. 2

## 2. Instruction for use

Extrude **Detaseal® bite** from cartridge with mixing cannula or with attached contouring tip onto teeth (Fig. 4). At terminal occlusion application from vestibular is possible. Let registration cure in terminal occlusion and remove from the mouth after 1 min. 30 sec.


Fig. 3

## 3. Disinfection

After removal from the mouth rinse registrate under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.

## Important hints

- Remove any retraction astringents or mouth wash solutions prior to the registration – could affect setting reaction.
- Do not intermix with condensation curing silicones.
- Avoid staining of clothing. Cured impression materials cannot be removed.
- Remove all residues of registration material from the patient's mouth.
- Please follow the instructions on the safety data sheet!


Fig. 4

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to get tested the allergenic potential before the application of the material.

Only for dental use by qualified personnel.

## Indications for use:

- Bite registration with crown and bridge impressions and for orthodontic indications, key for intraoral registration

## Technical Data:

DIN 13903

- **Mixed volume:** 50 ml (cartridge)
- **Mixing ratio:** 1 : 1
- **Colour code:** base: yellowish green catalyst: white
- **Mixing time:** omitted (Automix2-system)
- **Working time:** approx. 45 sec.\*
- **Setting time:** approx. 2 min. 15 sec.\*
- **Time in mouth:** approx. 1 min. 30 sec.
- **Final hardness:** approx. 90 Shore A
- **Linear dimensional change:** approx. 0.2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

## Ordering information:

- Detaseal® bite Automix2**
- Standard packing **02139**
- 2 x 50 ml cartridges
- 8 mixing cannulas, rose
- 8 contouring tips
- Ecopacking **02020**
- 8 x 50 ml cartridges
- 32 mixing cannulas, rose
- 32 contouring tips
- Mixing cannulas **02344**
- 6 mm, 48 pcs.
- Mixing cannulas Automix2
- rose, 48 pcs. **02707**
- Contouring tips **02611**
- 12 pcs.

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