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FANTESTIC® Z CORE DC

CE 0482

**Hersteller:**R-dental Dentalerzeugnisse GmbH • Winterhuder Weg 88 • 22085 Hamburg  
T 040-30707073-0 • F 040-30707073-73 • I www.r-dental.com**FANTESTIC® Z CORE DC (KwikMix®, 5 ml)**

REF FCS4065 (Farbe A2 Universal)

REF FCS4067 (Farbe weiß-opak)

**FANTESTIC® Z CORE DC (AutoMix, 25 ml)**

REF FCK4065 (Farbe A2 Universal)

REF FCK4067 (Farbe weiß-opak)

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REF FCK4067 (shade white-opaque)

**Description and indications**

FANTESTIC® Z CORE DC is a flowable, self-curing, high radiopaque core build-up composite (microhybrid composite) with light-cure option (dual-cure, DC).

FANTESTIC® Z CORE DC is indicated for adhesive core build-ups of vital and non-vital teeth, fillings by insufficient tooth structure, for adhesive luting of glass fibre reinforced composites (e. g. pre-silanized glass-fibre posts, like FANTESTIC® POST) and for adhesive luting ("cementing") of inlays, onlays, crowns and bridges. For all indications the use of a dual-cure adhesive system is required before application of the composite.

**Composition and characteristics**

FANTESTIC® Z CORE DC is based on multi- and difunctional methacrylates and inorganic filler particles with a size of 0.02-10 µm. The filler content is about 64 % by weight and 48 % by volume.

FANTESTIC® Z CORE DC exhibits a short setting time with low heat development. The dual-cure option enables the user to carry out luting ("cementations") and core build-ups in cases where light-curing cannot guarantee a polymerization.

FANTESTIC® Z CORE DC shows an excellent standability, consistency and wetting of tooth structure ("wettability").

FANTESTIC® Z CORE DC is delivered in commercial 5 ml-double-chamber syringes (KwikkMix®) and 25 ml-cartridges (AutoMix).

FANTESTIC® CORE DC meets the requirements of DIN EN ISO 4049, type 2, class 3 (Dentistry - Polymer based dental restorative material).

**Contraindications**

The application of FANTESTIC® Z CORE DC is contraindicated, if a dry working area or the recommended application technique are not possible and to any of the ingredients of FANTESTIC® Z CORE DC.

**Side effects**

In rare cases FANTESTIC® Z CORE DC may cause a sensitizing reaction at patients with a hypersensitivity to any of the ingredients. In such cases FANTESTIC® Z CORE DC should not be used further on.

In a direct contact with the pulp, irritations are possible. Deep cavities which are close to pulp areas should be covered with a thin layer of a calcium hydroxide liner (e. g. Cal de Luxe®, REF CDL2503) or a bioactive, light-cure pulp capping material with MTA-filler (ReViCal®, REF RVC1002) according to the instructions for use.

**Incompatibilities with other materials**

Do not store FANTESTIC® Z CORE DC in proximity of eugenol containing materials nor let it come in contact with eugenol containing materials. The proximity of unsealed products containing eugenol is also harmful. Further on unpolymerized FANTESTIC® Z CORE DC should not come into contact with eugenol containing materials.

**Recommended procedures****Preparation of the 5 ml-double-chamber syringe (KwikkMix®)**

Before first use screw off the original cap and waste it. Before screwing on a new mixing cannula the double-chamber syringe should be tested if both pastes will be extruded equally. Then the mixing cannula is put on and locked by turning it 90° clockwise.

FANTESTIC® Z CORE DC (KwikkMix®) will be mixed light protected and automatically using the MIXING CANNULA MCN 2 (REF MIK1041) + INTRA-ORAL-TIP MCN 2 (REF IOT2041) of the syringe system.

Instead of the INTRA-ORAL-TIP MCN 2 a fine root-canal application tip is available (ROOT-CANAL-TIP, REF RCT2041), with which the material could be applied directly and efficiently.

Further on instead of MIXING CANNULA MCN 2 + INTRA-ORAL-TIP MCN 2 a NEEDLE MIXING CANNULA (REF MIK1041N) with an integrated, bendable application needle is available.

For this type of double-chamber syringe a mixing pistol for a comfortable application is available (MIXING PISTOL N 5, REF MIP1041).

Attention: Store the used double-chamber syringe with the fixed used mixing cannula in the dark!

**Preparation of the FANTESTIC® Z CORE DC 25 ml-cartridge (AutoMix)**

Before first use screw off the original cap and waste it. Before screwing on a new mixing cannula the cartridge should be tested if both pastes will be extruded equally. Then the mixing cannula is put on and locked by turning it 90° clockwise.

FANTESTIC® Z CORE DC (AutoMix) will be mixed light protected and automatically using the yellow MIXING CANNULA LN (REF MIK1099) + INTRA-ORAL-TIP L (REF IOT2001) or using the blue MIXING CANNULA MLN (REF MIK2099) + INTRA-ORAL-TIP L (REF IOT2001).

Further on instead of MIXING CANNULA LN + INTRA-ORAL-TIP L a yellow NEEDLE MIXING CANNULA LNN (REF MIK1099N) or a blue NEEDLE MIXING CANNULA MLNN (REF MIK2099N) with an integrated, bendable application needle is available.

For this 25 ml-cartridge a commercial mixing pistol (MIXING PISTOL KN, REF MIP1083) is required.

Attention: Store the used double-chamber syringe with the fixed used mixing cannula in the dark!

**A. Adhesive luting ("cementing") of root canal posts****A 1. Isolation**

A dry area and use of rubber dam (REF KDL31XX, XX = depends on strength) or latex-free dam (REF KDL3500) is recommended.

**A 2. Preparation of the root canal posts**

Before an adhesive luting root canal posts must be silanized (e. g. silanizing with FANTESTIC® CerBond, REF FCF4100 or ZM-Bond®, REF ZMF4200). Pay attention to the manufacturer's instructions.

**A 3. Preparation of the root canal**

Before starting the preparation, the tooth must be cleaned of residues and so on. The lege artis prepared root canals are cleaned (e. g. with sodium hypochlorite solution) and rinsed. Remove excess adhesive from the root canal using paper points.

**The following adhesives are recommended for conditioning:****A 4.a Application of the self-etching adhesive FANTESTIC® UNIVERSAL BOND according to the instructions for use**

FANTESTIC® UNIVERSAL BOND (REF FUB6000) can be used as an adhesive for the self-etch technique, for the selective enamel etching or for the total-etch (Etch &amp; Rinse) technique.

**A 4.b Total-Etch (Etch & Rinse) technique and application of a dual-cure adhesive, like FANTESTIC® FLOWSIVE or FANTESTIC® FLOWSIVE 2 each together with FANTESTIC® FLOWSIVE ACTIVATOR (Etch & Rinse) according to the instructions for use**

Apply FANTESTIC® ETCH (REF EGS3000) precisely into the root canal for 15 s. Then rinse for about 20 s with water and dry with oil- and water-free air. The root canal dentin should be maintained in a slightly moist condition.

Equivalent amounts of FANTESTIC® FLOWSIVE (REF FFF3100) respectively FANTESTIC® FLOWSIVE 2 (REF FFF3105) each together with FANTESTIC® FLOWSIVE ACTIVATOR (REF FAF3150) will be mixed for dual-cure light protected for about 15 s.

Attention: Do not interchange the lids of the bottles.

**A 4.c Application of the self-etching, dual-cure adhesive system FANTESTIC® FLOWSIVE SE (Self-Etch) according to the instructions for use**

FANTESTIC® FLOWSIVE SE (REF FFF5155) Liquid A + Liquid B will be mixed light protected 1:1 within the enclosed FANTESTIC® Mixing moulds (REF FFF5180) for about 10 s.

Attention: Do not use component Liquid A or Liquid B alone and do not interchange the lids of the bottles of Liquid A and B.

The homogeneous mixture at this point is self etching (SE) for dentin and enamel.

Then the mixture will be applied with a brush into the slightly wet root canal and will be agitate intensively for about 15 s. Repeat this step 1-2 times. Remove excess adhesive using paper points. Then aerate the solvent for about 15 s in a weak oil and waterfree air stream. The mixture should be light-cured (if locally possible) for about 20 s to get an optimal chemical adhesion.

**A 5. Adhesive luting ("cementing") of root canal posts**

Pay attention to the manufacturer's instructions for use of the posts.

For definite adhesive luting ("cementing") of glass fibre posts or fibre reinforced composite (FANTESTIC® POST, REF FPS2X1X, x depends on geometry and size) the self-curing or dual-cure mode is valid.

Before every application a small amount of material (2-3 mm, size like rice corn or pea) should be pressed out and should be disposed of immediately.

FANTESTIC® Z CORE DC will be applied in the self-curing modus within the working time of about 1 min 30 s onto the post (it can be silanized using FANTESTIC® CerBond, REF FCF4100 or ZM-Bond®, REF ZMF4200), and then the post will be screwed into the root canal and held by slight pressure of fingers until the material is set.

Attention: FANTESTIC® Z CORE DC should not be luted directly in the root canal because of the applied adhesive a premature polymerization could occur.

The self-curing adhesive will auto-cure within about 3 min 30 s after working time.

For a fast stabilizing of the post the coronal area can be light-cured for about 20 s with a polymerization lamp (wavelength range 400-500 nm) with a light intensity of at least 1.000 mW/cm².

The filled MIXING CANNULA functions as a seal until the next application. Do not re-use the original cartridge cap.

Then the definite core build-up can be processed.

## B. Adhesive core build-up

### B 1. Isolation

A dry area and use of rubber dam (REF KDL31XX, XX = depends on strength) or latex-free dam (REF KDL3500) is recommended.

### B 2. Preparation of the cavity

It is prepared lege artis and old restorations and caries will be removed. Pay attention to the manufacturer's instructions for use of the posts.

### B 3. Pulp protection

Deep cavities which are close to pulp areas should be covered with a thin layer of a calcium hydroxide liner (e. g. Cal de Luxe®, REF CDL2503) or a bioactive, light-cure pulp capping material with MTA-filler (ReViCal®, REF RVC1002) according to the instructions for use.

The following adhesives are recommended for conditioning:

**B 4.a Application of the self-etching adhesive FANTESTIC® UNIVERSAL BOND according to the instructions for use**  
FANTESTIC® UNIVERSAL BOND (REF FUB6000) can be used as an adhesive for the self-etch technique, for the selective enamel etching or for the total-etch (Etch & Rinse) technique.

**B 4.b Total-Etch technique and application of a dual-cure adhesive, like FANTESTIC® FLOWSIVE or FANTESTIC® FLOWSIVE 2 each together with FANTESTIC® FLOWSIVE ACTIVATOR (Etch & Rinse) according to the instructions for use**

Apply FANTESTIC® ETCH (REF EGS3000) precisely into the root canal for 15 s. Then rinse for about 20 s with water and dry with oil- and water-free air. The root canal dentin should be maintained in a slightly moist condition.

Equivalent amounts of FANTESTIC® FLOWSIVE (REF FFF3100) respectively FANTESTIC® FLOWSIVE 2 (REF FFF3105) each together with FANTESTIC® FLOWSIVE ACTIVATOR (REF FAF3150) will be mixed for dual-cure light protected for about 15 s.

Attention: Do not interchange the lids of the bottles.

**B 4.c Application of the self-etching, dual-cure adhesive FANTESTIC® FLOWSIVE SE (Self-Etch) according to the instructions for use**

FANTESTIC® FLOWSIVE SE (REF FFF5155) Liquid A + Liquid B will be mixed light protected 1:1 within the enclosed FANTESTIC® Mixing moulds (REF FFF5180) for about 10 s.

Attention: Do not use component Liquid A or Liquid B alone and do not interchange the lids of the bottles of Liquid A and B.

The homogeneous mixture at this point is self-etching (SE) for dentin and enamel.

Then the mixture will be applied with a brush into the slightly wet root canal and will be agitated intensively for about 15 s. Repeat this step 1-2 times. Remove excess adhesive using paper points. Then aerate the solvent for about 15 s in a weak oil and waterfree air stream. The mixture should be light-cured (if locally possible) for about 20 s to get an optimal chemical adhesion.

A sticky oxygen inhibition layer will remain and must stay absolutely dry and clean until the application of FANTESTIC® Z CORE DC.

### B 5. Adhesive core build-up

Before application a small amount of material (2-3 mm, size like rice corn to pea) should be pressed out and should be disposed of immediately.

FANTESTIC® Z CORE DC will be applied directly in thin layers within a working time of about 1 min 30 s (by using a matrix if necessary and apply it from bottom upwards to prevent voids).

To facilitate the placement of FANTESTIC® Z CORE DC a matrix band should be placed around the prepared tooth.

FANTESTIC® Z CORE DC may be contoured by using a composite instrument. A flat-ended interproximal carver (Heidemann spatula) is recommended.

The self-cure system will auto-cure intraorally within about 3 min 30 s after working time.

Then the core build-up will be polymerized for 40 s with a polymerization lamp (wavelength range 400-500 nm) with a light intensity of at least 1.000 mW/cm<sup>2</sup>. With this technique an optimum of physical properties will be obtained.

An explorer can be used to test that the FANTESTIC® Z CORE DC has completely set. Remove the matrix not earlier than the material has set.

The final core build-up preparation can be carried out by using crown preparation burs.

### Precautions to application, product life and/or storage

Do not adjust the viscosity by using any resin.

Avoid contact of the etching gel with soft tissues, skin or eyes. If there is an accidental contact, flush immediately with copious amounts of water.

Pay attention to the Material Safety Data Sheets (MSDS >> [www.r-dental.com](http://www.r-dental.com)).

The ambient light of the dental lamp may start polymerization of the composite.

Do not use FANTESTIC® Z CORE DC in case of allergic reactions. Areas close to the pulp should be protected.

Avoid contact of FANTESTIC® Z CORE DC with soft tissues, skin or eyes, especially by anyone having known resin allergies. Commercial medical gloves do not protect against methacrylates in every case.

Do not store FANTESTIC® Z CORE DC in proximity of eugenol containing materials nor let it come in contact with eugenol containing materials. The proximity of unsealed products containing eugenol is also harmful. Further on unpolymerized FANTESTIC® Z CORE DC should not come into contact with eugenol containing materials.

Before first use screw off the original cap of the double-chamber syringe/cartridge and dispose it. Before screwing on a new mixing cannula the double-chamber syringe/cartridge should be tested to ensure both pastes will be extruded equally.

Before each application a small amount of material (2-3 mm, size like rice corn to pea) should be pressed out and should be disposed of immediately.

The filled MIXING CANNULA functions as a seal until the next application. Do not re-use the original double-chamber-syringe/cartridge cap.

Shelf life see label. Do not use after expiry date.

Store unopened material until the first use in the refrigerator. After the first use store FANTESTIC® Z CORE DC below 68 °F and store it horizontally, that means lying.

Opened FANTESTIC® CORE DC should be used within 3 months, see card box.

Avoid exposure to sunlight. Keep out of the reach of children!

FANTESTIC, Cal de Luxe, ZM-Bond and KwikkMix are registered trademarks of R-dental Dentalerzeugnisse GmbH.

Caution: For dental use only.

### Warranty

R-dental Dentalerzeugnisse GmbH warrants this product will be free from defects in material and manufacture. R-dental Dentalerzeugnisse GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. The user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and R-dental Dentalerzeugnisse GmbH's sole obligation shall be repair or replacement of the product.

### Limitation of Liability

Except where prohibited by law, R-dental Dentalerzeugnisse GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

© 2021 R-dental Dentalerzeugnisse GmbH. All rights reserved.