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Fleximeter®-Strips

BK 250, BK 251, BK 252, BK 253, BK 254

Bausch
WE MAKE OCCLUSION VISIBLE

Hersteller

Dr. Jean Bausch GmbH & Co. KG

Oskar-Schindler-Str. 4

50769 Köln - Germany

Tel: +49 221 70936-0

Fax: +49 221 70936-66

E-Mail: info@bauschdental.de

Internet: www.bauschdental.de

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Internet: www.bauschdental.de

Date of issue: 2019-12-05

1 General description and purpose of the medical device

Fleximeter®-Strips are flat, precisely cut silicone tabs in various colors with a disc- and rod-shaped, tapered end and defined thicknesses of 1.0mm, 1.5mm and 2.0mm. Thickness precisions: 1.0mm ± 0.2mm, 1.5mm ± 0.2mm, 2.0mm ± 0.2mm

For professional use only.

Indication

Fleximeter®-Strips are used for specification of the preparation distance or the preparation height when preparing teeth for restorations (e.g. crowns, bridges and double crowns), and to determine the material thickness needed to adjust bite guard splints. Fleximeter®-Strips, along with the liquid Arti-Spot® or Arti-Spray® occlusion spray, can be used in combination to mark contact points. Dentists can use the Fleximeter®-Strips separately or in pairs, in the patient's mouth as well as on models. Fleximeter®-Strips can also be used in the laboratory.

Contraindication

No known contraindications.

Side Effects

Possible side effects may be specific allergic reactions. However, given the product's history of more than ten years, there has not been a verified documented report of allergic reactions. Side effects or interactions may occur if the products are used with new, unknown products or materials.

2 Notes

- Before each use, Fleximeter®-Strips must be cleaned, disinfected and sterilized.
- Fleximeter®-Strips (instruments) can be steam-sterilized (humid heat, 134°C).
- Fleximeter®-Strips can be used in combination with the products Arti-Spot® or Arti-Spray®.
- Before disposing of Fleximeter®-Strips, they either have to be sterilized or disposed of with the normal contaminated practice waste.

3 Package content

- Fleximeter®-Strips (non-sterile – Type: according to product label)
- Instructions for use

4 Preparation

The instrument must be cleaned, disinfected and sterilized before each use. This also applies in particular to the initial application after delivery, since the instrument is delivered non-sterile. Please refer to the detailed requirements for reprocessing, described in section 6 Reprocessing.

5 Application

- Fleximeter®-Strips are used taking into account standard dental hygiene procedures.

- To accurately determine occlusal measurements, Fleximeter®-Strips are inserted between the affected teeth of the upper and lower jaw.
- Fleximeter®-Strips can be used to change the occlusal spacing.
- To determine the raising of the bite, Fleximeter®-Strips are inserted between the affected teeth of the upper and lower jaw.
- Depending on the respective indication, the entire surface or one end of the Fleximeter®-Strips can be used.
- Fleximeter®-Strips can be coated with a color layer of Arti-Spot® or Arti-Spray®.

Caution:

- If a patient accidentally swallows Fleximeter®-Strips, do not induce vomiting!
- If a patient aspirates Fleximeter®-Strips, immediately take action for removal and call an ambulance if necessary.

6 Reprocessing

In principle, the following should be followed: "Hygiene requirements for the processing of medical devices: Recommendation of the Commission for Hospital Hygiene and Infection Prevention (KRINKO) at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM)," Federal Health Bulletin 2012 • 55:1244-1310.

The instruments must be reprocessed immediately after each treatment (within a maximum of 2 hours):

- Thoroughly clean the instruments by hand, using an appropriate small brush with firm bristles.
- Note: The cleaning should be carried out in a water bath, without any additional cleaning agent, below the water surface, to achieve a sufficient cleaning of the Fleximeter®-Strips and to avoid protein fixation as well as to protect the environment against contamination with splashed water.
- Rinse the instruments with water (at least drinking water quality).
- Place the parts in a practice-customary cleaning and disinfecting bath. Examples:
 - o Becht Premium Concentrate
 - o Dürr Dental ID 213 instrument disinfection
 - o Pluradent Pluline instrument bath
 - o Schülke & Mayr gigasept® instru AF
 - o Note: "List of disinfectants and disinfecting procedures tested and approved by the Robert Koch Institute" or the VAH list of disinfectants.
 - o Note: To prevent protein fixation, formaldehyde-containing cleaning agents and disinfectants may be used only after adequate cleaning.
 - o Note: Strictly follow the instructions provided by the manufacturer of the cleaning agent or disinfectant. In particular, always follow the required concentration and residence time!
 - o Note: In case of automatic cleaning, the manufacturer's instructions have to be followed strictly!
- Final lavage of Fleximeter®-Strips with water (at least drinking water quality, recommendable: demineralized water with microbiological quality equal to drinking

- water)
- Drying.
- Visual check for corrosion, damaged surfaces, chipping, damage to shape and contamination. Damaged instruments shall be discarded (limited number of reconditioning cycles see Section "Reusability"). In case of residual contamination, the entire cleaning procedure with all stages shall be repeated (cleaning, intermittent rinsing, disinfection, final rinsing and drying).
- The instrument and all parts must be free of any residue and dry before further preparation.
- The instrument does not require servicing.
- All parts of the instrument should be packed and sealed in adequately sized, single-use sterilization bags that comply with EN 868-2/IF ISO 11607 (suitable for steam sterilization). Follow the instructions of the manufacturer of the sterilization bags and sealing machines and the current standard requirements.
- Sterilization must be completed in a validated procedure using moist heat in an autoclave in accordance with DIN EN 13060 Type B and DIN EN 285 and ANSI AAMI ST79. Follow the instructions for use of the autoclave manufacturer.
- Sterilize the instrument using moist heat (saturated steam) and a pre-vacuum procedure for **5 minutes at 134°C**.
- After completion of sterilization, the instrument parts must be stored dry and dust-free in the sealed sterilization packaging.
- The recommended storage period for sterile medical products is described in Standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling.

7 Reusability

Frequent reconditioning does not have any effect or limit on the use of the instrument, as the end of the product lifetime is determined by wear and damage due to use.

The use of damaged or soiled instruments is the responsibility of the user.

In case of non-compliance, all liability is excluded.

8 Disposal Instructions

Before disposing of Fleximeter®-Strips, they either have to be sterilized or disposed of with the normal contaminated practice waste.

9 Notification of Incidents

Serious adverse events occurring in connection with the product must be reported to the manufacturer and the competent authority of the Member State.

10 Symbols

- Manufacturer
- Production Date
- Steam sterilizable with saturated steam (autoclaving) at 134°C
- „Follow manufacturers instructions“
- Reference number
- Unique Device Identifier
- Conformity with the relevant EU regulations

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