GC Fuji ORTHO™ LC

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For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Bonding metal brackets and attachments.
- Bonding ceramic brackets.
- Bonding porcelain appliances.
- Band cementation only on extra band retentive in desired.

CONTRAINDICATIONS

- Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylic monomer or methacrylic polymer.
- Avoid use of this product on dentures.
- Cementation of stainless steel crown-retained appliances (e.g., Herbst Appliances).

DIRECTIONS FOR USE

- Etaminek megfelelően készítsék a fogakat:
 - Using a prophyl cup or brush, clean the bonding surfaces of the teeth with plain (non-fluorinated) pumice and water.
 - Rinse thoroughly with water.
 - Using a sponge or a cotton pellet, apply ORTHO CONDITIONER to the bonding surfaces of the teeth for 20 seconds.

Rinse thoroughly. Although successful results have been obtained with no enamel pretreatment, conditioning of the enamel bonding surfaces will increase the adhesive's bond strength. It is therefore recommended, especially for new uses of GC Fuji ORTHO LC.

OPTIONAL

Each the enamel bonding surfaces according to the etchant manufacturer's instructions. Rinse thoroughly. Note: The enamel bonding surfaces must be moist. An overly dry or desiccated enamel surface will adversely affect the bond strength. The liquid and the resin must be in contact with the enamel bonding surfaces for at least 10 seconds. Incorporate the remaining powder and mix thoroughly for an additional 10-15 seconds (total 20-25 seconds).

2. Apply the adhesive:
 - The standard powder to liquid ratio is 3.0g / 1.0g. Level large scope of powder to 20 seconds. For accurate dispensing of powder, tap the bottle gently. Do not shake the bottle. Hold the liquid bottle vertically and squeeze gently. Do not close bottles immediately after use.

For powder or a small mix, use the small scoop and one drop of liquid. Divide the powder into 2 equal parts. Mix the first portion with ALL the liquid and wait for about 10 seconds. Incorporate the remaining powder and mix thoroughly for an additional 10-15 seconds (total 20-25 seconds).

Note:

Working time is approximately 3 minutes from start of mixing at 23°C (73.4°F). Higher temperatures will shorten the working time, and lower temperatures will lengthen it.

Bonding Procedure

- Brackets
 - Coat the bonding surface of the bracket completely with the mixed resin.
 - Position the coated bracket on the tooth.
 - Press the bracket firmly against the enamel surface. Place an orthodontic scaler on the bracket for 20-30 seconds. Press the bracket again. Continue placing additional brackets. Place all brackets in the desired position. The use of a curing light or an adhesive may be necessary. If adhesive dries in a cement, "back" the bracket with a scaler for approximately 10 seconds.

AFTER REMOVING EXCESS ADHESIVE, BE CAREFUL NOT TO DISTURB THE BRACKETS BEFORE THE ADHESIVE IS CURED.

DECREASE THE BOND STRENGTH

1. **Use a curing light (or visible light)**, cure each bracket for 10-20 seconds each from the occlusal, mesial, distal and gingival aspects. This will reduce the bond strength. Use a curing light with a light meter to ensure adequate light output.
 - Insert leveling wire and complete initial procedures. A light cure wire (e.g., OTE NI) or equivalent is recommended for initial bonding.

1. Bonding ceramic brackets: Chemically retentive bases must be etched and silanated. If not etched and silanated, follow the manufacturer's directions for cleaning, anti-etching and silanating. No treatment is required for mechanically retentive bases.

2. Bonding to porcelain restorations: Prepare porcelain surfaces by etching and silanating according to manufacturer's directions. Some use of GC Fuji ORTHO LC.
 - Bonding brackets to amalgam/metal alloy restorations: Lightly roughen the metal surface with a fine diamond point before placing GC Fuji ORTHO LC.

B. Acrylic Appliances

- Use a scaler on the bracket or a rotary instrument, roughen and place relation holes in the internal surfaces of the appliance.
- Prepare the appliance as mentioned in section #1.
- To facilitate removal of the appliance, place petroleum jelly on the occlusal surfaces and a silanating agent on the lingual surfaces.
- Fit the internal surfaces of the appliance with GC Fuji ORTHO LC. Seal into proper position.

C. Removable orthodontic appliances

1. Apply adhesive during light curing in a distal-to-mesial direction, light-cure each side of the appliance for 30 seconds each on the buccal, palatal and lingual surfaces.
2. Appliance activation can begin at the patient's bedtime the same evening.

D. Debonding Procedure

- Brackets
 - Grasp the bracket with a ligature cutter and, while supporting the tooth with fingers from the lingual side, give a sharp clockwise rotation to the tooth. The resistance is met, by again holding in the opposite direction and after desiccating (using an air syringe) the tooth around the bracket. Remove any residual adhesive with a scaler or rotary instrument.

To remove bonded appliances, flex the appliance in multiple sites with hand pressure. If this does not remove the appliance, use a scaler or other appliance, desiccate the adhesive, and remove the appliance.

STORAGE

For optimum performance, store in a cool and dark place (4-25°C / 39-77.0°F).

PACKAGES

1. Introductory Package
 - 5g powder (1) / 8g (6.8ml) liquid (1) with a double-ended powder scoop, a plastic spatula, set of mixing pads (No. 22)
2. Standard Package
 - 40g powder (1) / 8g (6.8ml) liquid (2) with a double-ended powder scoop, a plastic spatula, set of mixing pads (No. 23)
3. Professional Package
 - 40g powder with a double-ended powder scoop (1) / 8g (6.8ml) liquid (1)

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton swabbed in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. In case of contact with or liquid with any other glass ionomer materials.
4. In rare cases the product may cause sensitivity to some persons. If such reactions are experienced, discontinue the use of the product and refer to a physician.
5. Personal protective equipment (PPE) such as gloves, face masks and eye protection should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gcgroup.com>

CLEANING AND DISINFECTING

MULTI-USE DELIVERY SYSTEMS to avoid cross-contamination between patients. The device requires mid-level disinfection. Immediately after use remove device and label for destruction. Discard device if damaged.


DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level regressive healthcare-grade infection control product according to regional/national guidelines.

Undesired effects/Reporting

If you experience any kind of undesired effect, reaction or similar events experienced by you or of this product, including those not listed in this document, please contact your local sales office or contact your nearest office as well as our internal vigilance system: vigilance@gc.dental

In any case you should always refer to the safety of this product.

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Prior to use, carefully read the instructions for use. **EN**

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GC Fuji ORTHO™ LC

LIGHT-CURED ORTHODONTIC BONDING ADHESIVE

For use only by a dental professional in the recommended indications. **RECOMMENDED INDICATIONS**

1. Bonding metal brackets and attachments.
2. Bonding ceramic brackets.
3. Bonding acrylic appliances to teeth with a moistened cotton roll immediately after band cementation only when extra band retention is desired.

CONTRAINDICATIONS

1. Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylic monomer or methacrylic polymer.
2. Bonding polycarbonate brackets.
3. Cementation of stainless steel crown-retained appliances (e.g., Herbst Appliances).

DIRECTIONS FOR USE

1. Enamel Preparation
 - a) Using a prophyl cup or a brush, clean the bonding surfaces of the teeth with plain (non-fluoriated) pumice and water.
 - b) Rinse thoroughly with water.
 - c) Using a sponge or a cotton pellet, apply ORTHO CONDITIONER to the bonding surfaces of the teeth for 20 seconds. Rinse thoroughly. Although successful results have been obtained with no enamel pretreatment, conditioning of the enamel bonding surfaces will increase the adhesive's bond strength. It is therefore recommended, especially for new users of GC Fuji ORTHO LC.

OPTIONAL: Etch the enamel bonding surfaces according to the etchant manufacturer's instructions. Rinse thoroughly.

NOTE: The enamel bonding surfaces must be moist. An overly dry or desiccated enamel surface will adversely affect the bond strength. The optimum level of surface moisture can be obtained by wiping the bonding surfaces of the teeth with a moistened cotton roll immediately prior to bonding procedure. If the bonding surfaces become desiccated during the bonding procedure, rehydrate by wiping them with the moistened cotton roll before bonding with GC Fuji ORTHO LC.

2. Powder and Liquid Dispensing
 - a) The standard powder to liquid ratio is 3.0 g / 1.0 g. 1 level large scoop of powder to 2 drops of liquid.
 - b) For accurate dispensing of powder, tap the bottle gently. Do not shake or invert.
 - c) Hold the liquid bottle vertically and squeeze gently.
 - i) Close bottles immediately after use.

For rebonds or a small mix, use the small scoop and one drop of liquid.

3. Mixing
 - a) Divide the powder into 2 equal parts. Mix the first portion with ALL the liquid and mix for about 10 seconds. Incorporate the remaining powder and mix thoroughly for an additional 10-15 seconds (total 20-25 seconds).

NOTE: The working time is approximately 3 minutes from start of mixing at 23°C (73.4°F). Higher temperature will shorten the working time, and lower temperatures will extend it.

4. Bonding Procedure
 - a) Coat the bonding surface of the bracket completely with the mixed adhesive.
 - b) Position the coated bracket on the tooth.
 - c) Press the bracket firmly against the enamel surface. Using an explorer or a scaler, remove the excess adhesive at this time. Press the bracket again. Continue placing additional brackets. Place all brackets in a quadrant or in the full arch. More than one mix of adhesive may be necessary. If bracket drift is a concern, "tack" the brackets by light curing for approximately 5 seconds.

AFTER REMOVING EXCESS ADHESIVE, BE CAREFUL NOT TO DISTURB THE BRACKETS BEFORE THE ADHESIVE IS CURED. MOVEMENT OF THE BRACKETS AT THIS POINT MIGHT DECREASE THE BOND STRENGTH.

- a) Using a curing light (470nm wavelength), cure each bracket for 10 seconds each from the occlusal, mesial, distal and gingival aspects. It is very important that the curing light be tested periodically with a light meter to ensure adequate light output.
- b) Insert leveling wire and complete initial procedures. A light force wire (e.g., .016 NiTi or equivalent) is recommended at the initial bonding. **Notes:**
 - 1) Bonding ceramic brackets: Chemically retentive bases must be etched and silanated. If not etched and silanated by the manufacturer, etch and silanate before bonding with GC Fuji ORTHO LC. If bracket bases become contaminated, follow the manufacturer's directions for cleaning, and/or etching and silanating. No treatment is required for mechanically-retentive bases.
 - 2) Bonding brackets to porcelain restorations: Prepare porcelain surfaces by etching and silanating according to manufacturer's directions. Bond with GC Fuji ORTHO LC as per above.
 - 3) Bonding brackets to amalgam/metal alloy restorations: Lightly roughen the metal surface with a disc or a fine diamond point before placing GC Fuji ORTHO LC.

- a) Using a microetcher or a rotary instrument, roughen and place retention holes in the internal surfaces of the bracket and petiole.
- b) Prepare the teeth as mentioned in section #1.
- c) To facilitate removal of the appliance, place petroleum jelly on the occlusal surfaces of the teeth.
- d) Fill the internal surfaces of the appliance with GC Fuji ORTHO LC.
- e) Seat into proper position.
- f) Immediately remove excessive adhesive from the periphery of the appliance.
- g) Moving the curing light pin in a distal-to-mesial direction, light-cure each side of the appliance for 30 seconds each on the buccal, palatal and occlusal surfaces.
- h) Appliance activation can begin at the patient's bedtime the same evening.

5. Bonding Procedure
 - A. Brackets
 - i) Grasp the bracket with a ligature cutter and, while supporting the tooth with fingers from the lingual side, give a sharp clockwise rotation to the bracket. If resistance is met, try again by twisting in the opposite direction after desiccating (using an air syringe) the tooth around the bracket. Remove any residual adhesive with a scaler or rotary instrument.
 - B. Appliances
 - i) To remove bonded appliances, flex the appliance in multiple sites with band removing pliers. If this does not remove the appliance, section the appliance, desiccate the adhesive, and remove the appliance.

STORAGE: Recommended for optimal performance, store in a cool and dark place (4-25°C / 39.2-77.0°F).

1. Introductory Package: 15 g powder (1), 8 g (6.8 mL) liquid (1) with a double-ended powder scoop, a plastic spatula, a set of mixing pads (No. Z2)
2. Standard Package: 40 g powder (1), 8 g (6.8 mL) liquid (2) with a double-ended powder scoop, a plastic spatula, a set of mixing pads (No. Z3)

3. Replacement Packages:
 - a) 40 g powder with a double-ended powder scoop
 - b) 8 g (6.8 mL) liquid

4. CAUTION
 1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
 2. In case of contact with eyes, flush immediately with water and seek medical attention.
 3. Do not mix powder or liquid with any other glass ionomer materials.
 4. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
 5. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referred in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or for The Americas <http://www.gcamerica.com>

They can also be obtained from your supplier.

CLEANING AND DISINFECTING


MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERGE: Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects - Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

as well as to our internal vigilance system: vigilance@gc.dental In this way you will contribute to improve the safety of this product.

Last revised: 06/2019

MANUFACTURED BY
GC CORPORATION
781 Hasunuma-cho, Tabatahri-ku, Tokyo 176-8685, Japan
REGULATORY AUTHORITY
GC AMERICA, INC.
3737 West 127th Street, Anso, IL 60863 U.S.A.

REGULATORY AUTHORITY
GC CORPORATION
781 Hasunuma-cho, Tabatahri-ku, Tokyo 176-8685, Japan
REGULATORY AUTHORITY
EU: GC EUROPE N.V.
Huisvestingweg 10, 3810 Lelystad, The Netherlands
3-9001 Leuven, Belgium TEL: +32 37 74 10 00
USA: GC AMERICA, INC.
3737 West 127th Street, Anso, IL 60863 U.S.A.
TEL: +1 708-597-0000
www.gcamerica.com

GC SOUTH AMERICA
Rua Heliópolis, 378 - Jardim - São Paulo, SP, BRAZIL
CEP: 06022-051 - TEL: +55 11-2925-0965
CNPJ: 06.708.717/0001-00
RESP. TEC: Mayara de Sales Ribeiro - CRDOSP 105.982

GC ASIA DENTAL PTE. LTD.
11 Tampines Central 1, #04-01, Singapore, Singapore 529541
TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
170 Bayswater Road, Bayswater, VIC 3009 Australia
TEL: +61 3 9391 8200

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OPTION: Etch the enamel bonding surfaces according to the etchant manufacturer's instructions. Rinse thoroughly.

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