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- Partial ceramic veneers
- Preparation for full ceramic restorations
- All requirements for adhesive bonding must be fulfilled

CONTRAINDICATIONS

Restorations made of Initial LRF BLOCK (GC) are contraindicated for patients diagnosed with excessive masticatory functions, in particular teeth grinders and clenchers. An absolute contraindication applies for restoring devitalized teeth with castor-oil based Initial LRF BLOCK (GC) for patients with hyperfunctions.

End crowns – premolars
Owing to the small adhesive surface and the small root diameters, endo crowns for premolars are contraindicated.

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PREPARATION DESIGN
When designing restorations, the following minimum dimensions should be maintained:

- Minimum wall thickness should be 1,5 mm, 1,0 mm at margins.
- Prepare margins with free chamfer rounded shoulder (Figure 1).
- Minimum thickness of the restoration should be 1.5 mm in pit and fissure areas and 2.0 mm in cusp areas.
- All internal edges and angles should be rounded. Avoid having margins in direct occlusal contact with the opposing tooth.

Cutter-shaped preparation should be avoided by all kind of restorations.

- For veneer, the preparation should be located in the enamel and a minimum thickness of 0,6 mm shall be respected.

MILLING
Please refer to respective instructions of the milling unit. Select the milling program for Initial LRF BLOCK (GC) or feldspar ceramic block materials (Wet grinding is advised). After grinding, check the restoration for discolorations linked to the grinding process

- Protuberances – Cracks or chips.

If the restoration is defective, discontinue the manufacture of the restoration.

FINISHING AND POLISHING
Final polishing
Finish and polish with coarse and fine silicone paste. DIAPOLISHER PASTE (GC) can be used for final polishing to obtain high gloss. Pay attention to margins and contact points when polishing the restoration. The correct speed must be ensured and generation of heat must be avoided. Clean the restoration before cementing.

CAUTION:
When adjusting or polishing the restorative, use a dust collector and wear a dust mask to avoid inhaling dust. Wear protective glasses.

2. Glaze firing using GC Initial LRF Glaze Paste
Glazing the restoration is to be done by using Initial LRF Glaze Paste (GC). Apply the Initial LRF Glaze Paste (GC) with a fine brush on the restoration. Always mix the Initial LRF Glaze Paste (GC) in the jar before use. When needed you can dilute the Initial LRF Glaze Paste (GC) using Initial Glaze Liquid (GC). Dilute the desired amount of the Initial LRF Glaze Paste (GC) on a separate plate. Never dilute in the jar. Initial Inviso powder stains (GC) can be added & mixed directly with the Initial LRF Glaze Paste (GC). There is no need to adjust the firing temperatures.

TIP:
For placement in a ceramic furnace use a firing tray. Depending on the restoration's geometry place the restoration dried on a platinum foil or on a with platinum eroded pin. A small amount of GC Initial Firing Foam can be used in combination with a metal pin (Do not use a max CAD crystallisation Tray – Ivoclar).

Remove the restoration immediately after firing from the furnace and leave to cool down. It is advised to limit the number of repeated firings to 3.

CEMENTATION
1. Preparing Etching protocol
Apply hydrofluoric acid gel (5-9%) for 60 seconds to the inner surfaces of the restoration.

- Wash with water spray or an ultrasonic cleaner and dry.
- Condition the etched surfaces with a silane coupling agent such as CERAMIC PRIMER I (GC) or G-Multi PRIMER (GC) and allow it to dry.

2. Cementing

- Cement with an adhesive resin cement such as G-CEM LinkAge (GC) or self-adhesive resin cement such as G-CEM LinkAce (GC).

NOTE:

- Please use CERAMIC PRIMER I (GC) or G-Multi PRIMER (GC), prior to use G-CEM LinkAge (GC) or G-CEM LinkAce (GC). Please refer to the respective instructions for use.
- In case the preparation is non-retentive, an adhesive resin cement such as G-CEM LinkAge (GC) is preferred.
- Phosphoric acid (35-37%) can be also used for the purpose of cleaning the surface, prior to scrubbing with a microbrush for 10-15 seconds.

CAUTION:
Hydrofluoric acid is for extraoral use only and must not be applied in the oral cavity.

CTE
CTE (25-500°C) [10-6/K]: 22,0 – 22,5 μm/(m°K)

CLASSIFICATION
Ceramic material, Type II / classification 2

SHADES

5 V- shades available in two translucencies:

- High translucency (HT): A1 HT, A2 HT, A3 HT, A3.5 HT, B1 HT
- Low translucency (LT): A1 LT, A2 LT, A3 LT, A3.5 LT, B1 LT

One special shade (BL)

- One transparency shade (BL)

SIZE
Available in 3 sizes: 12, 14, 14L PACKAGE - 5 blocks per box

PACKAGE
5 GC Initial LRF BLOCK per box
Glaze Paste: 1x GC Initial LRF Glaze Paste, Jar, 4g

STORAGE
Recommended for optimal performance, store at room temperature (4-29°C / 39-77°F) away from direct sunlight and high humidity. Shelf life:
GC Initial LRF BLOCK: 10 years from date of manufacture
GC Initial LRF Glaze Paste: 5 years from date of manufacture

CAUTION
1. In rare cases this product may cause sensitivity in some people. If such reactions are experienced, discontinue the use of product and consult a physician.
2. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Manufactured by
KLEMA DENTALPRODUKTE GMBH
KoblenstraÙe 3a
AT-6812 Meiningen
Tel: +43 5522 36837; Fax: +43 5522 36839; www.klema.at

DISTRIBUTED BY
GC CORPORATION
1761 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
www.gc.europa.net

GC EUROPE N.V.
Papestraat 3 • 3420 Aartselaar • Leuven 1240 • Interleuvenlaan 33, B-3001 Leuven, Belgium
TEL: +32 16 74 10 00 www.gc.europa.com

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Prior to use, carefully read the instructions for use.

EN

GC Initial™ LRF BLOCK

LEUCITE REINFORCED GLASS CERAMIC CAD/CAM BLOCK

For use only by a dental professional in the recommended indications.

INDICATIONS

Fabrication of full contour all ceramic indirect restorations: partial and full single anterior and posterior crowns, inlay, onlay, laminated veneer.

- Additional requirements:
- Preparation for full ceramic restorations
- All requirements for adhesive bonding must be fulfilled

CONTRAINDICATIONS

Hyperfunktion
Restorations made of Initial LRF BLOCK (GC) are contraindicated for patients diagnosed with excessive masticatory functions, in particular teeth grinders and clenchers. An absolute contraindication applies for restorations desvitalized teeth with restorations made of Initial LRF BLOCK (GC) for patients with hyperfunctions.


Endo crowns - premolars
Owing to the small adhesive surface and the small root diameters, endo crowns for premolars are contraindicated.

Bridge constructions
Since Initial LRF BLOCK (GC) consist of fine glass ceramic with a limited strength of approx. 200 MPa, this material is not suitable for the fabrication of any type of bridges construction.

PREPARATION DESIGN

When designing restorations, the following minimum dimensions should be maintained:

- Minimum wall thickness should be 1.5 mm, 1,0 mm at margins. Prepare margins with deep chamfer or rounded shoulder (Figure 1).
- Minimum thickness of the restoration should be 1.5 mm in pit and fissure areas and 2,0 mm in cuspal areas.
- All internal edges and angles should be rounded. Avoid having margins in direct occlusal contact with the opposing tooth.
- Gutter-shaped preparation should be avoided by all kind of restorations.
- For veneer, the preparation should be located in the enamel and a minimum thickness of 0,6 mm shall be respected.

MILLING

Please refer to respective instructions of the milling unit. Select the milling program for Initial LRF BLOCK (GC) or Feldspat ceramic block materials (Wet grinding is advised). After grinding, check the restoration for: discolorations linked to the milling process - Protuberances - Cracks or chips. If the restoration is defective, discontinue the manufacture of the restoration.

FINISHING AND POLISHING

1. **Mechanical polishing**
Finish and polish with coarse and fine silicone point. DIAPOLISHER PASTE (GC) can be used for final polishing to obtain high gloss. Pay attention to margins and contact points when polishing the restoration. The correct speed must be ensured and generation of heat must be avoided. Clean the restoration before cementing.

CAUTION:

When adjusting or polishing the restorative, use a dust collector and wear a dust mask to avoid inhaling dust. Wear protective glasses.

2. **Glaze firing using GC Initial LRF Glaze Paste**
Glazing of the restoration is to be done by using Initial LRF Glaze Paste (GC). Apply the Initial LRF Glaze Paste (GC) with a fine brush on the restoration. Always mix the Initial LRF Glaze Paste (GC) in the jar before use. When needed you can dilute the Initial LRF Glaze Paste (GC) using Initial Glaze Liquid (GC). Dilute the desired amount of the Initial LRF Glaze Paste (GC) on a separate plate. Never dilute in the jar. Initial InVivo powder stains (GC) can be added & mixed directly with the Initial LRF Glaze Paste (GC). There is no need to adapt the firing temperatures.

TIP:
For placement in a ceramic firing unit use a firing tray. Depending on the restoration's geometry place the restoration direct on a platinum foil or on a with platinum embedded pin. A small amount of GC Initial Firing Foam can be used in combination with a metal pin (Do not use e.max CAD crystallisation Tray – Ivoclar).

GC INITIAL LRF GLAZE PASTE FIRING PARAMETERS

Pre-heating temperature	Drying time	Raise of temperature	Vacuum	Final temperature	Holding time
450°C	4min	45°C/min	Yes	840-860°C	1min

Remove the restoration immediately after firing from the furnace and leave it to cool down. It is advised to limit the number of repeated firings to 3.

CEMENTATION

- Preparations Etching protocol**
 - Apply hydrofluoric acid gel (5-9%) for 60 seconds to the inner surfaces of the restoration.
 - Wash with water spray or an ultrasonic cleaner and dry.
 - Condition the etched surfaces with a siane coupling agent such as CERAMIC PRIMER II (GC) or G-Multi PRIMER (GC) and allow it to dry.

2. Cementing

Cement with an adhesive resin cement such as G-CEM LinkForce (GC) or self-adhesive resin cement such as G-CEM LinkAce (GC).

NOTE:

- Please use CERAMIC PRIMER I (GC) or G-Multi PRIMER (GC), prior to use G-CEM LinkForce (GC) or G-CEM LinkAce (GC). Please refer to the respective instructions for use.
- In case the preparation is non-retentive, an adhesive resin cement (such as G-CEM LinkForce (GC)) is preferred.
- Phosphoric acid (35-37%) can be also used for the purpose of cleaning the surface, preferably scrubbing with a microbrush for 10-15 seconds.
- CAUTION:** Hydrofluoric acid is for extraoral use only and must not be applied in the oral cavity.

CTE
CTE (25-500°C) [10 °°K]: 22,0 - 22,5 µm/(m°K)

CLASSIFICATION

Ceramic material, Type II / classification 2

SHADES

- 5 V- shades available in two translucencies:
- High Translucency (HT): A1 HT, A2 HT, A3 HT, A3.5 HT, B1 HT
- Low Translucency (LT): A1 LT, A2 LT, A3 LT, A3.5 LT, B1 LT
- One bleach shade (BL)

SIZE

Available in 3 sizes: 12, 14, 14L PACKAGE - 5 blocks per box

PACKAGE
Blocks: 5 GC Initial LRF BLOCK per box
Glaze Paste: 1x GC Initial LRF Glaze Paste, Jar, 4g

STORAGE

Recommended for optimal performance, store at room temperature (4-25°C / 39,2-77,0°F) away from direct sunlight and high humidity.

Shelf life
GC Initial LRF BLOCK: 10 years from date of manufacture
GC Initial LRF Glaze Paste: 5 years from date of manufacture

CAUTION

- In rare cases this product may cause sensitivity in some people. If such reactions are experienced, discontinue the use of product and consult a physician.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or for the Americas <http://www.gcamerica.com>

They can also be obtained from your supplier.

Last revised: 10/2016

Rx Only

CE 0483

IGC

MANUFACTURED BY
KLEMA DENTALPRODUKTE GMBH
Kobalcherstraße 3a
AT-6812 Meiningen
Fon +43 5223 366837, Fax +43 5222 366839, www.klema.at

RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +65 6546 7588
<http://www.gcaustralia.com>

GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
TEL: +81 2-9301-8200

GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium
TEL: +32 16 74 10 00
www.gceurope.com

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