Hera

HERAVEST® M

Ⓓ	Gebrauchsanweisung	3
Ⓖ	Instructions for use	6
Ⓕ	Mode d'emploi	9
Ⓔ	Instrucciones de uso	12
Ⓘ	Istruzioni per l'uso	15

KULZER
MITSUI CHEMICALS GROUP

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Stand: 04.2009

Graphite-free-phosphate-bonded, precision investment material for investing all types of chrome cobalt molybdenum and in precious metal - for use with silicone and gel duplicating materials

Storage

Keep the powder in a cool, dry place.

Do not store the mixing liquid below +5 °C as it is sensitive to low temperatures.

Working temperature

Approx. 22 °C (room temperature)

Mixing ratio:

Silicone-Duplicating

100 g Powder – 18 ml Liquid

160 g Powder – 29 ml Liquid

Gel-Duplicating

100 g Powder – 16 to 18 ml Liquid

160 g Powder – 25 to 29 ml Liquid

If the amount of liquid is reduced, the concentration must be adapted accordingly, e.g. to 50 % for 160 g/25 ml

Liquid concentration for 1000 ml Liquid		
Liquid concentration (examples)		
Concentration	Concentrate (ml)	dist. water (ml)
50 %	500	500
65 %	650	350
75 %	750	250

Concentrations							
Mixing ratio Powder: Liquid	Concentration %	Sachet 1 x 160 g Fabricating 1 Model		Sachet 2 x 160 g (320 g) Fabricating 2 Models		Sachet 3 x 160 g (480 g) Fabricating 3 Models Investing 1 Ring	
		Liquid (ml)	dist. water (ml)	Liquid (ml)	dist. water (ml)	Liquid (ml)	dist. water (ml)
Silicone-Duplicating							
160 g – 29 ml	75	22	7	43	15	65	22
	65	19	10	37	21	56	31
Gel-Duplicating *							
160 g – 29 ml	75	22	7	43	15	65	22
	65	19	10	37	21	56	31
160 g – 27 ml	65	17	10	35	19	53	28
	50	13	14	27	27	40	41
160 g – 25 ml	50	12	13	25	25	37	38

If the amount of liquid is reduced, the concentration must be adapted accordingly, e.g. to 50 % for 160 g/25 ml

* The model must only be duplicated with gel if it has been fully invested in investment

Expansion control (%)				
	Silicone-Duplicating		Gel-Duplicating*	
	Model	Investing	Model	Investing
Cobalt-chrome-Molybdenum (CoCrMo)	75	75	50–75	50–75
Gold casting alloy	65	65	65	50–65
* also refer to the mixing ratio				

Stirring

Use a spatula to stir the investment material by hand until the powder has been wetted thoroughly.

Mixing in vacuum

60 seconds.

Place in vacuum for 15 seconds but without switching the paddles on.

Mix for 60 seconds. Check the vacuum/paddles frequently for correct functioning. Inadequate vacuum results in inaccurate fit and bubbles on the casting.

Working time

Approx. 4 minutes at approx. 22 °C (room temperature)

The working time depends on the room temperature. Warmth curtails the working time.

Fabricating the model in duplicating silicone

Fabricating investment model

Vibrator speed 6000 min⁻¹

Before casting the silicone duplicating mould, spray it with Hera SWE 2000 silicone wetting agent, wait 2–3 minutes and blow it dry. The excellent flow properties facilitate casting non-porous models. Set the vibrator to high frequency and low intensity. Once the duplicating mould has been covered, switch the vibrator off immediately. Vibrating too intensely causes the material to segregate, impairs the precision of fit and results in rough surfaces (bubbles) as well as sedimentation.

Drying time Investment model

30 minutes. Carefully remove the model after 30 minutes. We recommend applying adhesive to the model or heating the working model briefly to approx. 50 °C to ensure that the wax pattern adheres properly.

Fabricating the model in duplicating gel

Fabricating investment model

Vibrator speed 6000 min⁻¹

Carefully remove the plaster model from the gel duplicating mould. Ensure that the mould is not distorted. The excellent flow properties facilitate casting non-porous models. Set the vibrator to high frequency and low intensity. Once the duplicating mould has been covered, switch the vibrator off immediately. Vibrating too intensely causes the material to segregate, impairs the precision of fit and results in rough surfaces (bubbles) as well as sedimentation.

Drying time Investment model

30 minutes / Carefully remove the model after 30 minutes.

Time taken for model to dry

60 minutes at 220 °C

Hardening the model

The model should be a light colour before immersion in the cold hardening liquid for 5–8 seconds. After immersion, return the model to the drying chamber for 2 minutes.

Full flask investing

Vibrator speed 6000 min⁻¹

Set the vibrator to high frequency and low intensity. Switch off the vibrator as soon as the pattern is covered with investment to avoid possible segregation of the investment (sedimentation). The model must only be duplicated with gel if it has been fully invested in investment.

The model must only be duplicated with gel if it has been fully invested in investment.

Suitable ring systems
a) Kulzer metal denture base system with paper sleeve
b) Other alternatives can also be used for casting metal denture bases

Placing the rings in a preheated furnace

Allow the investment mould to set for at least 30 minutes before placing it in the pre-heating furnace. Roughen or grind the top and bottom surfaces of the investment mould before placing it in the furnace.

Preheating

Preheating process for CoCrMo frameworks			
Program level	Temperature (°C)	Heat-rate (°C/min.)	Hold time (min.)
1	e.g. 23 °C room temp. – 270 °C	8	60
2	270 °C – 580 °C	8	60
3	580 °C – 1030 °C	8	40

Preheating process for precious metal alloy frameworks			
Program level	Temperature (°C)	Heat-rate (°C/min.)	Hold time (min.)
1	e.g. 23 °C room temp. – 270 °C	8	60
2	270 °C – 580 °C	8	60
3	580 °C – 700 °C	8	40

⚠ Safety Instructions

- Do not open the furnace during preheating as this would cause aggressive combustion.
- The preheat furnace should only be filled to half of its maximum capacity.
- Investments contain quartz. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer). Recommendation: Wear a protective face mask, type FFP 2 EN 149:2001. Prevent the release of dust when opening the bag and filling the mixing bowl.
- Rinse the empty bag with water before crumpling it up.
- Always moisten any dust before removing it from the workplace.
- After casting, place the completely cooled mould in water until it is fully moistened to prevent dust when devesting.
- Use an extractor with a fine dust filter when sandblasting.

We always attempt to include the latest results for our investment materials (development and application) in the instructions. Excellent precision of fit and surfaces are achieved with Kulzer equipment and materials. The results may vary if other materials (e.g. casting ring liners, waxes, equipment etc.) or different techniques are used.

We appreciate your comments, feedback, and suggestions.
Contact in the countries and more information are **available by the displayed QR code** or at our website **www.kulzer.com**

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Kulzer GmbH