GC INITIAL ZIRCONOXID-DISK CAD/CAM-ROHLINGE AUS ZIRKONOXID

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RECOMMENDED INDICATIONS
GC Initial Zirconia Disk is intended for use by professional dental technicians in the fabrication of custom, all-ceramic restorations. GC Initial Zirconia Disk is recommended for the manufacture of single-tooth and bridgework restorations such as crowns and bridges with one or two pontics. GC Initial Zirconia Disk can be used for anterior as well as posterior restorations.

KONTRAINDICATIONS
In rare cases this product may cause sensitivity in some patients. If such reactions are experienced, discontinue the use of product and consult a physician. Similar to all full ceramic restorations, zirconium dioxide frameworks are not suitable for patients suffering from bruxism. The thickness of the framework should be more than 0.4mm. In cases where space is restricted or the preparation is unsuitable, an alternative form of restoration should be chosen.

When designing crowns, bridges and primary components, the following minimum dimensions should be adhered to:

- Minimum wall thickness for anteriors: 0.4 mm
- Minimum wall thickness for laterals / posteriors : 0.6 mm
- Minimum wall thickness for abutment teeth: 0.6 mm
- Cross-section of connectors for anteriors: 6 mm²
- Minimum section of connectors for laterals/posteriors (for three-element bridges): 9 mm²
- Restorations should be designed in such a way as to support the veneering ceramic at the cusps area, so that these can be veneered with an almost uniform layer of ceramic.

PRECAUTION
1. Disk is pre-sintered. Handle with care.
2. Processing: only by skilled dental technicians.
Use only suitable parameters for processing ceramics.
After milling, check the framework for:

- discolorations
- outbreaks
- cracks

If the framework is deficient, discontinue manufacture of restoration.

RECOMMENDED INDICATIONS
GC Initial Zirconia Disk is intended for use by professional dental technicians in the fabrication of custom, all-ceramic restorations. GC Initial Zirconia Disk is recommended for the manufacture of single-tooth and bridgework restorations such as crowns and bridges with one or two pontics. GC Initial Zirconia Disk can be used for anterior as well as posterior restorations.

KONTRAINDICATIONS
In rare cases this product may cause sensitivity in some patients. If such reactions are experienced, discontinue the use of product and consult a physician. Similar to all full ceramic restorations, zirconium dioxide frameworks are not suitable for patients suffering from bruxism. The thickness of the framework should be more than 0.4mm. In cases where space is restricted or the preparation is unsuitable, an alternative form of restoration should be chosen.

When designing crowns, bridges and primary components, the following minimum dimensions should be adhered to:

- Minimum wall thickness for anteriors: 0.4 mm
- Minimum wall thickness for laterals / posteriors : 0.6 mm
- Minimum wall thickness for abutment teeth: 0.6 mm
- Cross-section of connectors for anteriors: 6 mm²
- Minimum section of connectors for laterals/posteriors (for three-element bridges): 9 mm²
- Restorations should be designed in such a way as to support the veneering ceramic at the cusps area, so that these can be veneered with an almost uniform layer of ceramic.

PRECAUTION
1. Disk is pre-sintered. Handle with care.
2. Processing: only by skilled dental technicians.
Use only suitable parameters for processing ceramics.
After milling, check the framework for:

- discolorations
- outbreaks
- cracks

If the framework is deficient, discontinue manufacture of restoration.

MILLING
Please refer to the instructions of the manufacturer of the milling unit. Use the milling program for zirconia.

COLORING
Use coloring liquids appropriate for green state (semi-sintered) zirconia.

SINTERING
Use a suitable high-temperature oven.

	Heat	Heat	Hold	Cool
Temperature (°C)	1000	1450	1450	1000
Time (hour)	2	4.5	2	1

Note: The parameters for heat-up and cool-down depend on object size as well as the oven control system. Please refer to the oven's operator manual.

FINISHING AND POLISHING
The hard-sintered GC Initial Zirconia Disk should be grinded/contoured further only if absolutely necessary. Only water-cooled diamond tools should be used for finishing and polishing. The local overheating may cause the fracture of the material. Do not apply high pressure.

Polishing of the crown surface is highly recommended to avoid damage to the opposing teeth. For higher esthetic results we recommend veneering the hard-sintered Zirconia substructures with GC Initial ZR-FS. For full contour restorations GC Initial IQ Lustre Pastes NF can be used. Please follow the instructions for use of the ceramic manufacturer.

PACKAGE
1 piece per box, diameter 98.5mm
HT: High Translucency in 4 different thicknesses

- GC Initial Zirconia Disk HT 12mm/14mm/18mm/25mm CIP

ST: Standard Translucency in 4 different thicknesses

- GC Initial Zirconia Disk ST 12mm/14mm/18mm/25mm CIP

CAUTION
1. Processing of GC Initial Zirconia Disk produces dust which can lead to eye and skin irritation, and may cause lung damage. Always ensure that the extraction unit on your milling machine is in good working order. Wear protective goggles, gloves and face mask.
2. Do not use autoclave to sterilize the framework.

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INDICATIONS RECOMMENDADAS
GC Initial Zirconia Disk is intended for use by professional dental technicians in the fabrication of custom, all-ceramic restorations. GC Initial Zirconia Disk is recommended for the manufacture of single-tooth and bridgework restorations such as crowns and bridges with one or two pontics. GC Initial Zirconia Disk can be used for anterior as well as posterior restorations.

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