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KKD Kondisil	V-1 Putty	V-3 low consistency
Chemical nature	Condensation type, elastomeric silicone impression material Type 0 – very high consistency – putty	Condensation type, elastomeric silicone impression material Type 3 – low consistency – light bodied
Meets the requirements of	EN ISO 4823	EN ISO 4823
Usage	for preliminary impressions	corrective impressions
Colour	white	blue
Physical properties		
Shore hardness	immediately: 70 – after 2 hours: 79	immediately: 33 – after 2 hours: 39
Max. strain in compression	1.3%	12.2%
Memory after deformation	98.6%	99.2%
Linear dimensional change	0.2%	1.1%
Storage	frost-free, do not store above 25° C	frost-free, do not store above 25° C
Shelf life	refer to package marking	refer to package marking
Mixing rate	1 measuring spoon Putty (B) (21 g) to 2.5 cm paste accelerator (P) (0,5 g)	8 cm base paste (B) (6 g) to 5 cm paste accelerator (P) (1 g)
Mixing	manual kneading. When using automatic mixing devices follow manufacturer's instructions	Push the thorn of the top cap into the opening of the tube and turn the cap repeatedly. Press the necessary amount of the material onto the mixing pad. When using automatic mixing devices follow manufacturer's instructions
Periods of time:		
Mixing time*	60 seconds	45 – 60 seconds
Total working time*	2 minutes	2 minutes
Time in mouth	3 ½ minutes	3 ½ minutes
Setting time*	5 ½ minutes	5 ½ minutes
Recommended time for pouring the model*	30 minutes to 72 hours	30 minutes to 72 hours
*(related to 23° C room-temperature and 50 % +/- 5 % relative humidity)		

#### Directions for use:

- ➔ excess accelerator and/or higher temperatures (related to 23° C room temperature and 50 % +/- 5 % relative humidity) shorten the working time, resp. setting time; insufficient accelerator and/or lower temperatures prolong working and setting time
- ➔ keep attention to a perfect adhesion. Moisture and saliva residues provide insulation. For that reason the impressions should be cleaned thoroughly and dried carefully before each subsequent impression step to assure perfect chemical bonding between the two materials.
- ➔ as model material we recommend to use, depending on final use, stones and plasters of class III (TopDent Hartgips, KKD Cal) and class IV (KKD Satin Roc, KKD Veradur 03)
- ➔ electroplating is possible with all usual dental copper or silver plating baths
- ➔ store impression under normal room conditions
- ➔ containers should be closed tightly after each use

#### Precautions

- ➔ avoid spots on clothing, since elastomeric impression materials are chemically inert
- ➔ avoid eye contact. If necessary wash eyes with copious amounts of water and seek ophthalmological help
- ➔ in order to avoid impurities we recommend to use vinyl gloves