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## Miramatic carpule canulas

GB Instructions for use

- Pull the plunger out up to the stop and insert the syringe cartridge
- (- Pull the plunger, ensuring it abuts the syringe cartridge's rubber plunger)
- Attach canula to Miramatic-head
- Check marks at syringe head (syringe head has to be locked, both marking lines shape one line)

Important: Only Miraject-canulas must be used in Miramatic syringes.

### Notes for Miramatic-Aspi:

Self aspiration becomes activated by a small counter pin inside the syringe head. The automatic aspiration is active whenever the pressure on the plunger is being released during use. When pressure is applied to the plunger, the cartridge is pushed onto the counter pin inside the syringe head. The counter pin condenses the cylinder cartridge's rubber membrane, resulting in negative pressure, activating the self-aspiration. By releasing the pressure, the rubber returns to its original shape, avoiding further leaking of the injection liquid. You activate the manual aspiration by pressing on the disk with your thumb. If the needle has perforated a blood vessel this becomes visible inside the cartridge due to the negative pressure.

### Manufacturer Information on Reprocessing of Medical Devices according to DIN EN ISO 17664:2004

Manufacturer: Häger & Werken GmbH & Co. KG, Ackerstr. 1, 47269 Duisburg, Germany  
Tel.: +49 (203) 99 269-0, Fax: +49 (203) 29 92 83, www.hagerwerken.de

Product: Miramatic Spritze, Miramatic Aspi (REF 355 463 – 355 469; 355 471 – 355 475; 355 495 – 355 496; 355 499)


**Warning notices:** No special notices

**Limitation of reprocessing:** Frequent reprocessing has only got a small impact on this instrument. The instrument shelf life is mainly determined by wear and damage by use.

**Place of handling:** Clean surface with a disposable cloth or paper towel.

**Storage and transport:** No special requirements.

**Preparations for cleaning:** Dismantle the product and if necessary (stiffness) lubricate the Miramatic head with a drop of instrument oil.

**Cleaning:** Only use suitable cleansers, carefully following the instructions of the respective manufacturer. After cleaning, remove cleanser thoroughly under running water. Do not exceed 93 °C for mechanical cleaning.

**Disinfection:** Only use suitable disinfecting solutions, carefully following the instructions of the respective manufacturer. Remove disinfecting solution thoroughly under running water. Do not exceed 93 °C for mechanical disinfection.

**Maintenance:** No special requirements.

**Packaging:** Standardised packing material for sterilisation may and should be used.

**Sterilisation:** Vapour sterilisation at 134 °C for 5 minutes or vapour sterilisation at 121 °C for 15 minutes.

**Control/Functional check:** Sight check on damages, wear, deformation.

**Storage:** No special requirements.

The a.m. instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine control are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

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Fabrikant: Häger & Werken GmbH & Co. KG, Ackerstr. 1, 47269 Duisburg, Germany