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Mirawhite[®] shine

GB

Mirawhite shine is a gel for gentle tooth whitening. The mode of action of Mirawhite shine bleaching gel is based on Ardox-X[®] technology. Ardox-X[®] is a hydrocarbon oxoborate complex, which effects tooth whitening by an oxidation process. Unlike peroxides, no free radicals are built during the oxidation process.

Gentle and Safe

The new tooth whitening technology of Mirawhite shine combines basic characteristics of tooth whitening without compromise in the result:

- The Mirawhite shine tooth whitening gel is non-abrasive and water and saliva soluble. It causes no damages or sensitivities of teeth and gingiva.
- Unlike products containing peroxide, Mirawhite shine does not soften bondings of fillings, so it is ideal for people with many fillings.
- Mirawhite shine is gentle on the teeth. Active oxygen molecules dissolve stains from the teeth, gently removing discoloration from the tooth surface. The teeth restore their natural whiteness.

Application



Mirawhite shine is applied with the brush twice a day, in the morning and evening after toothbrushing, for a period of at least 7 days.

Simple steps to naturally white teeth:

1. Brush your teeth as usual and let them dry, so that no moisture is left on the tooth surfaces. The dry teeth may not get in touch with the lips or the tongue.
2. Twist the bottom part of the pen until the gel is visible on the brush. Apply a thin layer of the whitening gel tooth by tooth in an up and down motion. If necessary, re-moisten the pen with gel.
3. The mouth should remain open, so the gel can dry. The lips must not get in touch with the teeth for min. 60 seconds. Then rinse with water. In order to obtain an optimal result, do not eat, drink or smoke for 30 minutes after application.

Re-cap the pen after use in order to avoid drying-out.

Best results are achieved after two weeks.

Should some teeth (e.g. canine teeth) feature a poorer whitening effect, these can be re-treated individually. Omit already sufficiently whitened teeth.

Whitening Results

Use the shade guide in front of the mirror to measure the whiteness of your teeth before and after treatment with Mirawhite shine. If used as directed, the first whitening results become visible after 7 days. Whitening by up to 5 shades can be achieved within two weeks.

The final result always depends on the cause of the discoloration. Teeth that appear brown or yellow – mainly as a result of the consumption of nicotine, coffee, tea and red wine, are whitened by an average of 4 shades. Grayish teeth – caused by medicaments – are whitened by an average of 2 – 3 shades.

Precautions

1. Keep out of the reach of children.
2. Use Mirawhite shine only as directed.
3. Avoid eye contact. In case of eye contact, rinse thoroughly with water.
4. Store in a cool, dry place (8 – 21 °C). Protect against heat and direct sunlight.
5. Do not use during pregnancy.

Mirawhite shine gel with Ardox-X[®] has been tested by:

- TUFTS University, Boston, USA
- Loma Linda University, California, USA
- University of Indiana, USA
- St. Bartholomew's and Royal London School of Medicine, UK
- Queens University Belfast, UK
- Dental Academic Centre Amsterdam, NL
- Mubio Research BV, University of Maastricht, NL
- ACTA Amsterdam, NL

For further information please visit www.miradent.de.

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Hager & Werken GmbH & Co. KG