# 3M | RelyX™ Universal

- ⓔⓃ **Resin Cement**
- ⓔⓄ **Befestigungs-Composite**
- ⓔⓇ **Ciment de scellement composite**
- ⓔⓉ **Cemento composito**
- ⓔⓈ **Cemento definitivo de resina**
- ⓔⓅ **Cimento de resina**
- ⓔⓁ **Composietcement**
- ⓔⓁ **Ρητινώδης κονία για μόνιμες συγκολλήσεις**
- ⓔⓋ **Resincement**
- ⓔⓂ **Yhdistelmämuovipohjainen kiinnityssementti**
- ⓔⓁ **Resincement**
- ⓔⓃ **Kompositt-sement**

Instructions for Use  
Gebrauchsinformation  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones de uso  
Instruções de uso

Gebruiksaanwijzing  
Οδηγίες Χρήσης  
Bruksanvisning  
Käyttöohje  
Brugsanvisning  
Bruksanvisning

**en ENGLISH****Product Description**

RelyX™ Universal is a dual-curing universal resin cement that is used for both adhesive and self-adhesive cementation of indirect restorations and posts. RelyX Universal is available in various colors matching the RelyX™ Try-In Pastes.

RelyX Universal resin cement is used in combination with Scotchbond™ Universal Plus Adhesive for adhesive cementation. Scotchbond Universal Plus Adhesive is self-etching. Optionally, the adhesive strength can be enhanced further through additional etching of the tooth structure with Scotchbond™ Universal Etchant.

Standard cases as well as post cementations can be treated securely and efficiently when RelyX Universal is used in self-adhesive mode.

RelyX Universal Resin Cement is available in an Automix Syringe developed by 3M. The Automix Syringe has a built-in rotary valve that opens and closes when the appropriate mixing tip is attached and removed. The mixing tips ensure high quality of the mix and minimum waste of pastes.

RelyX Universal uses Scotchbond Universal Plus Adhesive as a universal restoration primer.

Scotchbond Universal Plus Adhesive contains MDP monomer as well as silane. The MDP monomer bonds to Zirconia and alumina ceramic surfaces as well as to metals. The silane component allows the adhesive to bond to glass ceramics.

 These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products, please refer to the corresponding Instructions for Use.

**Intended Purpose**

Intended purpose: dental luting agent for cementation of indirect restorations.

Intended users: trained dental professionals who have theoretical and practical knowledge regarding the use of dental products.

Clinical benefit: Holds indirect restorations in place with and without adhesive to restore tooth integrity and maintain functionality.

**Indications****Use of RelyX Universal in Self-Adhesive Mode or in Adhesive Mode with Scotchbond Universal Plus for the Following Indications:**

- Final cementation of all-ceramic, composite, or metal inlays and onlays, crowns, and bridges
- Final cementation of posts made of ceramic, glass fiber-reinforced composite or metal, and screws
- Final cementation of all-ceramic, composite, or metal restorations on implant abutments

**Use of RelyX Universal in Adhesive Mode with Scotchbond Universal Plus for the Following Indications:**

- Final cementation of all-ceramic or metal Maryland bridges and 3-unit inlay/onlay bridges
- Final cementation of all-ceramic or composite veneers and occlusal veneers (tabletops)

Scotchbond Universal Etchant is used for etching of enamel and dentin using the acid etching technique.

**Precautionary Measures****For Patients and Dental Personnel**

RelyX Universal base paste and Scotchbond Universal Plus: contact with eyes may cause severe eye damage. Wear eye protection. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Scotchbond Universal Etchant: contact with eyes or skin may cause burns. Wear suitable protective clothing, gloves, and eye/face protection. In case of contact, flush area with plenty of water and seek medical advice.

**For Patients**

Scotchbond Universal Plus and RelyX Universal contain substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of Scotchbond Universal Plus in patients with known acrylate allergies and avoid use of RelyX Universal in patients with known acrylate or perester allergies.

If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use.

**For Dental Personnel**

RelyX Universal and Scotchbond Universal Plus contain substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured paste. If skin contact occurs, wash skin immediately with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If RelyX Universal or Scotchbond Universal Plus contacts the glove, remove and discard glove, wash hands immediately with soap and water, and re-glove. If allergic reaction occurs, seek medical attention as needed.

Scotchbond Universal Plus is highly flammable. Avoid static discharge. Keep away from heat, sparks, open flame, pilot lights, and other sources of ignition.

Please report any serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

3M Safety Data Sheets/Safety Information Sheets can be obtained at [www.3M.com](http://www.3M.com) or from your local subsidiary.

**Tooth Preparation**

- Protection of the pulp: cover areas in close proximity to the pulp by applying spots of a hard-setting, eugenol-free calcium hydroxide material to protect the pulp from irritation. To avoid bite increases during the cementing of the later restoration, the pulp protection should be applied prior to taking an impression for the final restoration.
- Seat temporary restorations using a eugenol-free product (e.g., RelyX™ Temp NE). The use of temporary cements containing eugenol can inhibit the polymerization process of Scotchbond Universal Plus and RelyX Universal during the final cementing.
- Do not use hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) or sodium bicarbonate (NaHCO<sub>3</sub>) as the residues can adversely affect the adhesive strength and setting reaction of Scotchbond Universal Plus and RelyX Universal.
- Do not use substances such as desensitizers, disinfectants, astringents, dentin sealants, rinsing solutions containing EDTA, scanning powder, etc., after the final cleaning. Their residues may have a detrimental effect

on the bonding strength and setting reaction of the adhesive or resin cement.

#### Preparation of the Cavity/Tooth Stump

- ▶ Remove loose preparation debris by spraying with water.
- ▶ Clean the stump or the cavity thoroughly, rinse with a water spray, and lightly air dry in only 2–3 bursts of water-free and oil-free air, or use, e.g., cotton pellets to dry it off. Do not overdry.
  - The cavity should be just dry enough that the surface has a slightly glossy appearance. As with any permanent cement, overdrying can lead to post-operative sensitivity.

#### Root Canal Preparation for Post Cementation

- ▶ Remove the root canal filling, leaving at least 4 mm of root canal filling (guttapercha) apically.
- ▶ Prepare the root canal for the post, rinse with water, and dry with paper tips.
- ▶ The use of a rubber dam during the cementation of posts is recommended.

#### Preparation for Veneers

- ▶ Prepare and clean the tooth as usual for veneers.
- ▶ The tooth must then be etched as described in the section “Etching of the Tooth Structure”.

#### Preparation for Maryland and Inlay/Onlay Bridges

Maryland and inlay/onlay bridges are not recommended for patients with bruxism or periodontitis.

Abutment teeth must have an adequate enamel surface for bonding. They should be healthy or only slightly restored and the periodontal conditions should be good. It is the sole responsibility of the dentist to ensure proper selection of indication and technique. The guidelines of the relevant national professional associations must be observed for such indications.

- ▶ Prepare retentive elements as necessary.
- ▶ The tooth must then be etched as described in the section “Etching of the Tooth Structure”.
- ▶ Use a rubber dam and keep the cavity free of any contamination during cementation.

#### Preparation for Chairside CAD/CAM Restorations

Please follow the Instructions for Use of the restoration material.

#### Trial Fitting of the Restoration

- ▶ Remove the temporary restoration and thoroughly clean the tooth structure of any residue from the temporary cement.
- ▶ Make a trial insertion of the final restoration and check its fit and contact points.
  - If a low viscosity silicone is used for the trial fit, the silicone residue must then be removed completely.
- ▶ Use RelyX Try-In Pastes for the selection of the matching cement shade.

#### Pre-Treatment of Restorations

Avoid any contamination whatsoever of the surfaces to be treated during pre-treatment and until final cementation.

#### Pre-Treatment of Zirconia, Aluminum Oxide, and Metal Surfaces

Please follow the Instructions for Use of the restoration material or abutment. In the absence of specific instructions, we recommend the following procedure:

- ▶ Sandblast surfaces to be cemented. Use 30 µm to 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to obtain a rough, matte surface.

- ▶ Clean the blasted surface with alcohol and dry it with water-free and oil-free air.
- ▶ If sandblasting was done before the trial fitting (e.g., in the dental laboratory), the surfaces must be cleaned after the trial fitting (we recommend using sodium hypochlorite solution, NaOCl 5%), rinsed thoroughly with water, and dried.
- ▶ The surfaces of Zirconia, aluminum oxide, or metal restorations must not be treated with phosphoric acid.

#### Use of a primer

RelyX Universal has very strong self-adhesive bond strength to Zirconia, aluminum oxide, and metal. Optionally, Scotchbond Universal Plus Adhesive can be used as a primer if additional adhesive strength is required. We recommend using Scotchbond Universal Plus as a primer for veneers, occlusal veneers (tabletops), inlay/onlay bridges, and Maryland bridges.

- ▶ Use the disposable applicator to apply Scotchbond Universal Plus to the entire surface to be luted and rub it in for 20 sec. During treatment, rewet the disposable applicator as required.
- ▶ Subsequently direct a gentle stream of air over the liquid for a minimum of 5 sec until it no longer moves and the solvent has evaporated completely.

#### Pre-Treatment of Etchable Glass Ceramic Restorations

Please follow the Instructions for Use of the restoration material. In the absence of specific instructions, we recommend the following procedure:

- ▶ Etch the bonding surface of the glass ceramic restoration with an HF etching gel.
- ▶ Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.
- ▶ Use the disposable applicator to apply Scotchbond Universal Plus Adhesive to the entire bonding surface and rub it in for 20 sec. During treatment, rewet the disposable applicator as required.
- ▶ Subsequently direct a gentle stream of air over the liquid for a minimum of 5 sec until it no longer moves and the solvent has evaporated completely.

Alternatively, a suitable silane (e.g., RelyX™ Ceramic Primer) can be applied in accordance with the Instructions for Use.

#### Pre-Treatment of Composite Restorations

Please follow the Instructions for Use of the restoration material. In the absence of specific instructions, we recommend the following procedure:

- ▶ Sandblast surfaces to be cemented. Use 30 µm to 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to obtain a rough, matte surface.
- ▶ Clean the blasted surface with alcohol and dry it with water-free and oil-free air.
- ▶ If sandblasting was done before the trial fitting (e.g., in the dental laboratory), the adhesive surfaces must be cleaned after the trial fitting (we recommend using sodium hypochlorite solution, NaOCl 5%), rinsed thoroughly with water, and dried.
- ▶ Use the disposable applicator to apply Scotchbond Universal Plus Adhesive to the entire bonding surface and rub it in for 20 sec. During treatment, rewet the disposable applicator as required.
- ▶ Subsequently direct a gentle stream of air over the liquid for a minimum of 5 sec until it no longer moves and the solvent has evaporated completely.

#### Pre-Treatment of Glass Fiber-reinforced Posts

We recommend the following procedure when using RelyX™ Fiber Posts:

- ▶ Clean the RelyX Fiber Post with alcohol and dry it with water-free and oil-free air.
- ▶ When using RelyX Fiber Posts, application of a primer is not required.

Please follow the instructions for use for other glass fiber-reinforced posts. In the absence of specific instructions, we recommend the following procedure:

- ▶ Use the disposable applicator to apply Scotchbond Universal Plus Adhesive to the entire bonding surface and rub it in for 20 sec. During treatment, rewet the disposable applicator as required.
- ▶ Subsequently direct a gentle stream of air over the liquid for a minimum of 5 sec until it no longer moves and the solvent has evaporated completely.

Alternatively, a suitable silane (e.g., RelyX Ceramic Primer) can be applied in accordance with the Instructions for Use.

#### Etching of the Tooth Structure (Optional, Depending on the Type of Material and Indication)

Scotchbond Universal Plus Adhesive is self-etching. Adhesion can be maximized by selective enamel etching or by the etching of enamel and dentin (“Total Etch Procedure”). Uncut enamel must be etched.

Total-Etch is required for veneers, occlusal veneers/tabletops and for Maryland and inlay/onlay bridges.

**Dentin may be etched solely when Scotchbond Universal Plus is used subsequently; do not etch the dentin for the self-adhesive use of RelyX Universal.**

- ▶ We recommend using a rubber dam.
- ▶ Follow the instructions for use for the application of the etching gel.
- ▶ When using Scotchbond Universal Etchant, remove the syringe from the foil pouch, unscrew the cap, and attach an etching gel needle tip.
- ▶ To avoid contamination of the etching gel syringe during treatment, we recommend using the etching gel syringe with a standard dental protective sleeve; see Hygiene and Disinfection.
- ▶ Apply the etching gel to the precise area of the enamel and/or dentin and allow it to react for 15 sec. Then rinse thoroughly for 15 sec with a water spray while using suction.
- ▶ Dry with water-free and oil-free air or with cotton pellets. Do not overdry!
- ▶ Avoid any contamination of the etched surface. If contamination occurs, repeat the etching, rinsing, and drying.
- ▶ Continue immediately with the application of the adhesive.
- ▶ After etching, remove the etching gel needle tip from the syringe and replace the cap. Dispose of the etching gel needle tip.

#### Application of Scotchbond Universal Plus Adhesive to the Tooth Structure (Optional, Depending on the Type of Material and Indication)

RelyX Universal must be used with Scotchbond Universal Plus Adhesive and enamel etching for restorations that require maximum adhesive strength to the enamel (veneers, occlusal veneers/tabletops, Maryland and inlay/onlay bridges). For all other indications, RelyX Universal may be used with or without Scotchbond Universal Plus Adhesive, provided that the instructions for the restorative material do not stipulate otherwise. RelyX Universal may be used solely with Scotchbond Universal Plus Adhesive. Other adhesives have not been approved for use.

- ▶ Pour the required quantity of Scotchbond Universal Plus in a mixing well.
- ▶ Protect the adhesive in the mixing well from light.
- ▶ Carefully close the bottle tightly after dosing. The lid perceptibly snaps into place when pressed.
- ▶ Avoid any contamination with blood, saliva, or sulcus fluid during application.
- ▶ A rubber dam is recommended to keep the treated area sufficiently dry.
- ▶ Use the disposable applicator to apply the adhesive to the entire tooth structure and rub it in for 20 sec. Avoid contact between the adhesive and the oral mucosa.
- ▶ If necessary, rewet the disposable applicator during treatment.
- ▶ Subsequently direct a gentle stream of air over the liquid for a minimum of 5 sec until it no longer moves and the solvent has evaporated completely. In case of application in the root canal, remove any excess with a paper tip and blow a gentle stream of air over it for 5 sec as well.
- ▶ Light curing of Scotchbond Universal Plus is not required; it may optionally be polymerized for 10 sec at a light wavelength range between 400 nm and 500 nm with any commonly available light-curing device.

#### Preparation of the Automix Syringe Before Application of the Cement

##### When Using a New RelyX Universal Automix Syringe

- ▶ Remove the Automix Syringe from the foil pouch immediately before use only and discard the foil pouch.
- ▶ After removal from the foil pouch, RelyX Universal resin cement must be completely used within 6 months and before lapse of the expiration date, so the date of removal from the pouch must be noted in a suitable place.

##### Before Every Application

- ▶ If the Automix Syringe has been stored in a refrigerator, allow it to reach room temperature before use.
- ▶ There is a built-in rotary valve at the front opening of the Automix Syringe that protects the cement. This valve will open, allowing the paste to be extruded, only if the mixing tip has been attached correctly.
- ▶ Do not attach the mixing tip until immediately before use; lock it into place by turning it 90° in the direction of the printed arrow until it stops. This motion opens the Automix Syringe and allows the paste to be extruded. Do not press the plunger unless the mixing tip has been correctly attached.
- ▶ To extrude the paste, apply moderate centric force; do not bend the plunger. Do not apply excessive force. If extrusion force appears high, check to ensure that an unused mixing tip has been attached and correctly locked.
- ▶ Extrude a small quantity of paste onto a pad so that the base and catalyst pastes in the Automix Syringe are precisely mixed in equal amounts. Discard the paste which has been squeezed out.
- ▶ To avoid contamination of the Automix Syringe during treatment, we recommend using the Automix Syringe with a standard dental protective sleeve; see Cleaning and Disinfection.

##### Cementing the Post

Do not use Lentulo spirals to insert the cement in the root canal as this can excessively accelerate setting.

- ▶ For application into the root canal, attach an elongation tip contained in the package (or a RelyX Unicem Applicator Elongation Tip) to the mixing tip until it locks into place. Support the mixing tip to prevent any bending.

- ▶ Insert the elongation tip as deeply as possible into the root canal and apply RelyX Universal, beginning apically. Keep the tip of the elongation tip immersed in the cement and slowly move the elongation tip upwards as the level of the paste rises.
- ▶ Do not remove the elongation tip from the cement until the root canal has been completely filled.
- ▶ This so-called immersion filling should be carried out slowly and take at least 5 sec to minimize the inclusion of air bubbles.
- ▶ Place the post in the root canal filled with cement; apply moderate pressure to hold it in position. We recommend rotating the post slightly during insertion to avoid the inclusion of air bubbles.
- ▶ Remove the excess of cement with a suitable instrument or a cotton pellet.
- ▶ Polymerize the cement through the post when using translucent posts. The exposure time depends on the translucency of the post being used; for RelyX Fiber Post, it is 20 sec.
- ▶ Remove and discard the used mixing tip with the RelyX Elongation Tip or RelyX Unicem Aplicap Elongation Tip after application; the syringe closes itself with the built-in rotary valve.

### Cementing of Veneers, Inlays/Onlays, Crowns, and Bridges

#### Application to Tooth or Restoration

Protect the working area from contamination with water, blood, saliva and sulcus fluid during the application and setting phase.

- ▶ Wet the entire area of the cavity walls, the cavity floor, and, if necessary, the surface of the restoration to be cemented evenly with RelyX Universal.
  - Keep the opening of the mixing tip immersed in the material during the entire application to prevent the inclusion of any air bubbles.
- ▶ The paste requires a certain amount of time to flow through the mixing tip. The flow rate cannot be accelerated by increasing the pressure on the plunger.
- ▶ As soon as the pressure on the plunger decreases, the material flow stops, and the paste begins to set. Do not use force to press out paste that has set as this could cause damage to the mixing tip and the Automix Syringe.
- ▶ For veneers: owing to the relatively short working time of RelyX Universal (self-curing begins after 2 min), we recommend cementing a maximum of 2 veneers at the same time.
- ▶ Seat the restoration and stabilize long enough for the cement to set fully; see Times.
- ▶ Remove and discard the used mixing tip after application; the syringe closes itself with the built-in rotary valve.

#### Removal of Excess

Notes on oxygen inhibition: As is the case with all composite materials, RelyX Universal is subject to oxygen inhibition, i.e., the upper layer (about 50 µm) which is in contact with atmospheric oxygen during the polymerization does not harden. Leaving sufficient excess permits the removal of the uncured layer during finishing/polishing without leaving behind a deficit. Large excess quantities can be removed more easily!

#### Options for Removing Excess

##### Removal of Excess for Best Aesthetic Results

- ▶ Remove the excess before polymerization, e.g., using a sponge pellet, while using a suitable instrument to hold the restoration in place.

- ▶ Cover the margins of the restoration with a glycerin gel to prevent oxygen inhibition.
- ▶ Immediately light-cure or wait for the self-curing to finish (see Polymerization and Shaping).

#### Fast Removal of Excess

- ▶ Use an appropriate instrument (e.g., scaler) to remove excess cement after brief light exposure (approximately 2–3 seconds with a conventional polymerization device); clean the proximal areas in particular without delay. Use a suitable instrument to hold the restoration in position during this procedure.
- ▶ The tack-curing time of 2–3 seconds for RelyX Universal resin cement is significantly longer than for RelyX Unicem or RelyX Ultimate.
- ▶ The 3M Elipar™ S10 and Elipar™ DeepCure-S LED light-curing devices have a built-in tack-cure mode that permits a controlled brief curing of 1 sec.

#### Polymerization and Finishing

- ▶ The working time and setting time of RelyX Universal (see Times) depend on the ambient and oral temperature. The times shown are based on conditions relevant for the dental practice. As is the case with every composite cement, the setting of RelyX Universal slows down significantly at room temperature.
- ▶ RelyX Universal is a dual-curing cement and therefore also sensitive to natural or artificial light. The working time is significantly reduced during application under operating lights!
- ▶ We recommend light-curing the cement through the restoration when doing ceramic and composite work. Select the exposure times appropriate for the number of surfaces (see Times). Recommended light-curing devices: Elipar™ S10 or Elipar™ DeepCure-S LED light-curing devices.
- ▶ Finish any remaining uneven areas on the edges of restorations and polish the marginal area with diamond devices, aluminum oxide coated discs (e.g., Sof-Lex™), and diamond polishing paste.
- ▶ Then check the occlusion.

#### Times

	RelyX Universal min:sec	Scotchbond Universal Plus min:sec
Maximum working time from start of mixing	02:00	–
Light curing in a light wavelength range between 400 nm and 500 nm and a light intensity $\geq 1200 \text{ mW/cm}^2$ *	–	
– single surface, from occlusal	00:10	(optional) 00:10
– any other surface additional	00:10	–
– RelyX Fiber Post/ RelyX Fiber Post 3D posts, from occlusal	00:20	–

	RelyX Universal min:sec	Scotchbond Universal Plus min:sec
Self-curing – setting time after start of mixing	06:00	–

\* If lower light intensity is used, the exposure time must be extended in accordance with the instructions of the light-curing device manufacturer.

#### After Finishing

- ▶ Carefully check the sulcus of the treated teeth and the surrounding areas; remove any cement residues still remaining. This is especially important when cementing to implant abutments.

#### Cleaning and Disinfection of the Automix Syringe

- ▶ Before use of the syringe, put on new, unused disposable gloves during all clinical applications of the Automix Syringe to reduce the risk of cross-contamination.
- ▶ For hygienic reasons, use a new mixing tip and protective sleeve for every patient to avoid risk of infection.
- ▶ Place the protective sleeve over the Automix Syringe with the previously attached mixing tip and use the tip of the mixing tip to carefully punch through the end of the protective sleeve.
- ▶ Avoid contact between the reusable parts (e.g., the body of the syringe) and the patient's mouth.
- ▶ A mixing tip (disposable product) is used to apply the product to the patient and is discarded after being used once on the patient.
- ▶ After using the Automix Syringe in the protective sleeve, remove the sleeve carefully so that there is no contamination from the outer surface of the protective sleeve. Discard the used protective sleeve along with the mixing tip.
- ▶ After removal of the protective sleeve, wipe the Automix Syringe thoroughly for 30 seconds with a ready-to-use CaviWipe™ or an equivalent agent.
- ▶ Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals during use and reprocessing.

Note: If the Automix Syringe/etching gel syringe is not visibly clean, i.e., if, contrary to expectations, there is contamination with blood or saliva, for example, discard the Automix Syringe immediately.

Do not place the Automix Syringe/etching gel syringe into a disinfection bath or a cleaning-disinfection device (washer-disinfector).

#### Cleaning and Disinfection of Scotchbond Universal Plus General

**Use new, unused disposable gloves during all clinical applications of the bottle to eliminate the risk of contamination.**

The dental assistant should dose the material beforehand to avoid any contamination of the bottle from the dentist's gloves during treatment.

Note: If the bottle is not visibly clean, i.e., if, contrary to expectations, there is contamination with blood or saliva, for example, discard the bottle immediately.

Do not put the adhesive bottle into a disinfection bath or a cleaning and disinfecting device (washer-disinfector).

#### Inspection, Maintenance, Testing

Examine the medical device for damage, discoloration, and contamination before every use. Do not use damaged and contaminated devices under any circumstances.

#### Notes

- RelyX Universal and Scotchbond Universal Plus can prematurely polymerize when exposed to natural or artificial lighting. Therefore avoid intensive light exposure during application.
- Disposable applicator, etching gel needle tip, RelyX Universal mixing tip, RelyX Elongation Tip/RelyX Unicem Aplicap Elongation Tip are single-use devices.
- RelyX Universal mixing tip, RelyX Elongation Tip/RelyX Unicem Aplicap Elongation Tip would be damaged when removing set material.
- Please keep the single-use accessory in the original bag to guarantee the traceability of the batch.

#### Disposal

- Dispose of Scotchbond Universal Etchant syringes only when they are completely empty. The residue of the etching gel – if heavily diluted with water or neutralized – can be poured down the drain.
- Dispose of the RelyX Universal contents or container in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste to avoid health risks because of improper handling.

#### Storage and Shelf Life

Store RelyX Universal in the foil pouch at 02–25 °C/36–77 °F.

After removal from the foil pouch, use RelyX Universal within 6 months and before the lapse of the expiration date. Constantly high humidity influences the setting and must be avoided.

Store Scotchbond Universal Plus and Scotchbond Universal Etchant at 02–25 °C/36–77 °F.

Do not use after the expiration date.

#### Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

#### Warranty


3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

#### Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental, or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

## Symbol Glossary

Reference Number & Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.4 Use-by date		Indicates the date after which the medical device is not to be used.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1 5.4.2 Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.
ISO 15223-1 5.4.4 Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive with notified body involvement.
Medical Device		Indicates the item is a medical device.
Flammable		Product is classified as flammable. Source: Globally Harmonized Symbols for Classifying Hazardous Chemicals.

Reference Number & Symbol Title	Symbol	Description of Symbol
Corrosive		Product is classified as corrosive. Source: Globally Harmonized Symbols for Classifying Hazardous Chemicals
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).
Recycle		Indicates packaging is recyclable with mixed paper.
Low-density polyethylene		Indicates plastic component made of low-density polyethylene is recyclable. Official Journal of the EC; Commission Decision (97/129/EC)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

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3M, Elipar, RelyX, Scotchbond and Sof-Lex are trademarks of
Used under license in Canada. © 2020, 3M. All rights reserved.
3M Customer Care/Safety Data Sheet Information: U.S.A. 1-800-634-2249.
3M or 3M Deutschland GmbH.
3M, Elipar, RelyX, Scotchbond et Sof-Lex sont des marques de commerce de 3M ou 3M Deutschland GmbH.