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en INSTRUCTIONS FOR USE

SUTURE MATERIALS of SILK

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MATERIAU DE SUTURE en SOIE

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MATERIALE DA SUTURA in SETA

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MATERIAL DE SUTURA de SEDA

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MATERIAL DE SUTURA, à base de SEDA

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SILK SUTURE MATERIALS● **SERAFLEX®****WARNINGS AND PRECAUTIONS**

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide; for SERAFLEX® in cassette packs: gamma irradiation).

- The product is intended for single use only.

This does not apply to cassette packs.

- The shelf life after opening should be considered for cassette packs.
- Do not resterilize or re-use the product as critical changes in performance and safety are possible.

- Do not use after expiration date.

- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or unintentionally opened prior to use.

- Use the product only for the intended indication.
- Please contact the manufacturer for further information.

DESCRIPTION AND INTENDED PURPOSE

SERAFLEX® is a natural, non-absorbable, surgical suture material derived from the cocoon fibres of the silkworm moth (*Bombyx mori* L.). SERAFLEX® is a braided, silicone-coated suture and is available dyed black with Logwood black (haematin), Colour Index No. 75290 or undyed. SERAFLEX® in the cassette pack is supplied in an isopropanol-containing preserving solution.

The suture materials made of silk are available in various gauge sizes and lengths. They are supplied in combinations with various types of needles made of surgical steel or without needles. They comply with the current edition of the European Pharmacopoeia (EP) monograph for "Sutures, sterile non-absorbable" (fila non resorbilia sterilia)". In deviation from this standard, the maximum thread length in the case of the cassette packs is not 4, but 100 m.

Silk threads are suitable for all wound sutures and ligatures in which non-absorbable thread material is indicated and long-term stability of the suture material is not required. Threads left in tissues for a long time can lose tensile strength. Silk threads can also be used as holding sutures and for marking purposes.

INSTRUCTIONS FOR USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the physician, the surgical technique, and the type and size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security.

In order to prevent needle damage during suturing it is recommended that the needle be grasped at an area between the end of its first third (as seen from its end) and its midpoint. Unneeded suture material may be used for ligature or may be connected with surgical eye needles. For a list of suitable needles and application instructions, please refer to the Instructions for Use of Serag-Wiessner sterile eye needles. The suitability of other manufacturers' products has not been tested. When using other brands, their specific instructions must be observed.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to foreign material. Ligatures are a special type of suture intended to close hollow organs. In the body, silk suture material is slowly and progressively encapsulated by connective tissue.

CONTRAINDICATIONS

The product is not intended for use at the heart or in the central vascular or nervous systems. Because of the loss of tensile strength that occurs in the body, silk threads must not be used when prolonged stability of the suture material is required. The device must not be used in patients with allergies to any of the product components. The use of silk threads on the lungs and tendons is not recommended.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The requirements of the suture material and the techniques required vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. Like all foreign bodies, the device can give rise to calculus formation when in prolonged contact with salt-containing solutions (e.g. in the biliary or urinary tract).

In the case of sutures in infected tissue or tissue predestined for infection, the use of suitable ancillary measures is called for, as the multifilament thread structure can prolong infection. Surgical instruments should be handled with care in order to avoid damage to threads. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled and disposed of with particular care.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of the device:

local irritation, inflammatory reactions to the foreign body; formation of suture fistulas or granulomas; allergic reactions.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label. Store cassette packs protected from UV-radiation.

INSTRUCTIONS FOR REMOVAL FROM THE PACK

Minibox: For easy removal, the thread should be withdrawn from the Minibox not obliquely, but vertically upwards.

Minibox - microsutures in protective sleeve: Remove the sleeve from the minibox attachment; grasp the needle with the needle holder and withdraw the thread from the protective sleeve; if a **short thread** is required, it should be cut through together with the protective sleeve; in the case of **double-armed threads**, two individual threads can, if required, be withdrawn after the loop of thread at the end of the protective sleeve has been cut through.

Cassette pack (also called suture dispenser or outer pack): Remove the tamperproof closure (dispenser cap closure) and use the enclosed label to note the expiry date on the outer pack; observe the following instructions to guarantee aseptic removal for a period of 6 months after opening. Flip up the dispenser cap at the clip and fix it to the back part of the cap. If the outer pack has already been opened, i.e. the dispenser cap closure has been removed or opened, **disinfect** the inner surface of the dispenser **before use** with a commercially available alcohol-based product. Pull out the end of the suture with a sterile instrument. Discard the first 10 cm of the suture. Pull out the desired length of suture, cut it off flush, and close the cap. **After finishing the procedure, but no more than 4 hours later, or if required again disinfect** the inner surface of the dispenser and the cap with a commercially available alcohol-based product and close the cap.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG