sta-seal

DE

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sta-seal

Instructions for use

Precision impression material, silicone based, condensation curing, medium consistency – medium bodied

Fig. 1


Fig. 2

1. Preparation of the impression

Depending on the impression technique select an appropriate tray. For optimal adhesion of the impression to the tray apply a thin layer of **sili Adhesive** spray or **sili Adhesive** fluid onto the tray and let dry for 60 sec. (Fig. 1).

2. Mixing

Extrude **sta-seal** and **cat[®] universal FUTUR** from tubes in equal length onto mixing pad (Fig. 2). Be careful, that the extruded strings show a uniform cross-section. With dosage of liquid catalyst 1 cm **sta-seal** to 1 drop of catalyst. Mixing 30 sec. with spatula until homogeneous mix and uniform colour is achieved.

3. Disinfection

After removal rinse impression under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.

4. Model casting

Casting between 2 hours and 24 hours after removal of the impression. Recommended materials: Dental plasters of class III and IV.

5. Galvanisation

Impressions can be electroplated with commercially available copper and silver plating baths.

6. Cleaning of trays

Remove set material with a blunt instrument. **sili Adhesive** film can be easily removed with **sili Solvent** spray (free of FCHC). Work in well ventilated rooms. Clean and disinfect as usual.

Important working hints

- Overdosing of catalyst or increased temperatures accelerate, underdosing or decreased temperatures retard the setting time.
- Catalyst bottles or tubes should be closed tightly immediately after use.
- Avoid staining of clothing. Hardened impression materials cannot be removed.
- Avoid contact of catalyst with eyes, may cause irritation. If necessary flush eyes with plenty of water and seek medical care immediately.
- Don't leave any residual material in the patient's mouth.
- Please follow the instructions of the safety data sheet!


Avoid contact with eyes and skin!
Wear suitable eye protection and gloves!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to get tested the allergenic potential before the application of the material.

For use by trained specialist staff.

sta-seal / cat[®] universal FUTUR:

Indications for use:

- Universal impression material for edentulous impressions
- Mucostatic and relined impression

Technical Data:

DIN EN ISO 4823 – Type 2

Mixing ratio:

- sta-seal / cat[®] universal FUTUR**
1 cm sta-seal = 1.10 g /
1 cm cat[®] = 0.05 g
- sta-seal / sta-seal cat[®]**
1 cm sta-seal = 1.10 g /
1 drop cat[®] = 0.02 g

Colour code:

- sta-seal:** lime green
- cat[®] universal FUTUR:** blue
- sta-seal cat[®]:** light blue

Mixing time:

approx. 30 sec.

Working time:

1 min. 15 sec.*

Setting time:

approx. 3 min. 45 sec.*

Time in mouth:

approx. 2 min. 30 sec.

Strain in compression:

14.50 %

Recovery from deformation:

99.30 %

Linear dimensional change:

0.45 %

Application:

At 23 °C ± 2 °C / 73 °F ± 4 °F,
50 ± 5 % rel. humidity

Storage:

15 °C / 59 °F
25 °C / 77 °F

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

Ordering information:

sta-seal
Standard packing 02010
160 ml tube

tube 3-pack 02011
3 x 160 ml tubes

cat[®] universal FUTUR
5 x 35 ml paste 02015

sta-seal cat[®] 02018
5 x 10 ml bottles

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 07243/510-0 · Fax: 07243/510-100
www.detax.de · post@detax.de

Made in Germany

05/2016

C-silicones

medium bodied

medium consistency

silone®

normal setting
final hardness approx.
45 Shore A

Universal material for all dentulous impressions, single phase impression technique, for partial-, coppering and inlay impressions

medium bodied

medium consistency

sta-seal

normal setting
final hardness approx.
35 Shore A

Precision impression material for edentulous impressions, atrophic conditions, flabby ridges, mucostatic and relin impressions

heavy bodied

high consistency & plasticity

sta-seal f

ideal final hardness
extended plastic phase
extended setting time

Special impression material for functional impressions with extended setting, functional border moulding, functional impression in one session with si-plast trays

catalyst paste

catalyst liquid

cat^p / cat^f

inert system
non-irritating

cat^p universal paste, silasoft® cat^f, liquid, sta-seal cat^f, liquid

special instrument

accessory

Deta-Cut

sili adhesive

sili solvent

Sterilizable, curved, for cutting the preliminary impression

adhesive for C-silicones

removal of the adhesive and cleaning of trays

Ordering information:

silone®
Standard packing **02008**
160 ml tube

tube 3-pack **02009**
3 x 160 ml tubes

sta-seal
Standard packing **02010**
160 ml tube

tube 3-pack **02011**
3 x 160 ml tube

sta-seal f
Standard packing **02012**
3 x 80 ml dosing syringes
10 ml sta-seal cat^f bottle

Large packing **02013**
9 x 80 ml dosing syringes
3 x 10 ml sta-seal cat^f bottles

Catalyst
cat^p universal FUTUR
5 x 35 ml paste **02015**

sta-seal cat^f **02018**
5 x 10 ml bottles

Deta-Cut **02323**

sili adhesive
spray can, 250 ml **02051**

fluid, bottle 25 ml **02053**

sili solvent
spray can, 250 ml **02052**

Auxiliaries
easy application syringe **02140**
1 pc.

Mixing pad **02031**
large, vellum paper
30 sheets

Mixing pad **02032**
small, vellum paper
30 sheets

Mixing pad **02033**
large, plastic coated
30 sheets

Mixing pad **02034**
small, plastic coated
30 sheets

Mixing spatula **02036**
stainless steel
with wooden handle

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