# roeko

## Stay-put

IMPREGNATED

### Instructions for use

**Indications:**  
Stay-put impregnated consists of a braided polyester cord in which a silver-plated copper wire is braided. The retraction cord is impregnated with aluminium chloride hexahydrate for the temporary retraction and haemostasis of the gingival margin. The mechanical use of the retraction cord causes compression. The compression contracts the treated tissue and causes haemostasis. The aluminium chloride hexahydrate supports the efficacy of the cord during preparation of the gingiva or during impression procedures due to its astringent properties and the haemostasis in the gingival margin. It thereby evokes a slight shrinkage (contraction) of the gingiva, a reduction of secretion from the sulcus and a sealing of smaller vessels which are bleeding by causing coagulation of proteins in the superficial layers of the damaged or non-damaged gingival margin. The result is the presentation of the cervical and subgingival preparation limit as a requirement for dental impressions, for preparations of cavities, for placing silicate, composite or amalgam fillings, for fitting matrices and strips in the gingival sulcus area and for fitting inlays, crowns and bridges.

The impregnated retraction cord is available in size x-fine, fine, medium and thick.

**Composition:**  
Polyester cord with a silver-plated copper wire braided into it and impregnated with aluminium chloride hexahydrate ( $\text{AlCl}_3 \times 6\text{H}_2\text{O}$ , Ph. Eur.).

1 cm cord contains		
0 x-fine:	0.5 mg $\text{AlCl}_3 \times 6\text{H}_2\text{O}$	REF 521 020
1 fine:	0.5 mg $\text{AlCl}_3 \times 6\text{H}_2\text{O}$	REF 521 021
2 medium:	0.75 mg $\text{AlCl}_3 \times 6\text{H}_2\text{O}$	REF 521 022
3 thick:	1.0 mg $\text{AlCl}_3 \times 6\text{H}_2\text{O}$	REF 521 023

**Application:**  
For best results when using Stay-put impregnated cut a piece of cord corresponding to the circumference of the tooth, overlapping is not necessary. After preparing and drying the tooth first place the centre of the cord into the sulcus by using a retraction cord packer. Next work out towards the ends of the cord. During the preparation for the impression a certain retraction of the gingiva is achieved through the mechanical influence of the cord and the astringent effect of the aluminium salts. After approximately 1 – 3 minutes the sulcus is opened and any bleeding which has been caused during preparation work is halted, so that immediately after removal of the retraction cord the impression, in particular to reproduce the cervical preparation limit, can be taken. The area which has been treated with Stay-put impregnated should be rinsed with water.

**Contraindications:**  
Stay-put impregnated with aluminium chloride hexahydrate should not be used in case of acute gingivitis or acute necrotizing ulcerative gingivitis. The use of Stay-put impregnated in the inflamed area should be considered carefully, weighing the benefits against the risks, as the retraction of the gingival margin could cause tissue damage to a certain extent.

**Important interactions:**  
Coagulation can occur if contact is made with alkaline reacting substances or with phosphates. Gelatine preparations will be tanned.

**Note:**  
Stay-put impregnated does not influence the ability to drive a vehicle or attend a machine.

Stay-put impregnated has a bitter taste due to the aluminium chloride hexahydrate.

Avoid touching the cord with electrotomy and laser equipment.

**Warning:**  
Medical devices should be kept out of reach of children.

**Use and storage:**

- for single use only
- the package must be closed immediately after use. Store in a dry, dark place below 25°C.
- Stay-put impregnated should not be sterilized in the autoclave or by dry heat.
- shelf life is 3 years. After opening, the bottle may be used until the expiration date. Do not use after the expiration date shown on the label and packaging.

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