LINK Series®

INSTRUCTIONS

EN

TMS LINK SERIES

TMS LINK Series Pins consist of single and double shear-self-threading pins linked to resilient plastic shanks. No chucks are necessary because the shanks fit directly into any latch type handpiece.

TMS LINK Series shanks aid in the alignment of pins to the drilled channels for smooth and safe insertion.

TMS LINK Series Pins and KODEX® Drills are color coordinated for quick and easy identification.

TMS LINK Series Pins may be installed with any standard or miniature head latch type handpiece at low speeds of approximately 750-1000 RPM.

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TMS LINK SERIES CATALOG INFORMATION

SINGLE-SHEAR COMPLETE KIT CONTENTS: 20 TMS LINK SERIES Pins • 1 KODEX Drill • 1 Universal Hand Driver

DOUBLE-SHEAR COMPLETE KIT CONTENTS: 20 TMS LINK SERIES Pins • 2 KODEX Drills • 1 Universal Hand Driver

SINGLE-SHEAR OR DOUBLE-SHEAR BULK KIT CONTENTS: 50 TMS LINK SERIES Pins

TMS LINK SERIES PIN	COMPLETE KITS	CAT. NO. BULK KITS	CAT. NO. KODEX DRILL	CAT. NO. COLOR CODE	DRILL SIZE
MINUTA SINGLE-SHEAR	L-501	L-502	K-89	PINK /.0135" (.350mm)	
MINKIN SINGLE-SHEAR	L-511	L-512	K-91	RED /.017" (.425mm)	
MINIKIN SINGLE-SHEAR, TITANIUM ALLOY	L-811	L-812	K-91	RED /.017" (.425mm)	
MINIM SINGLE-SHEAR	L-521	L-522	K-92	SILVER /.021" (.525mm)	
MINIM DOUBLE-SHEAR	L-531	L-532	K-92	SILVER /.021" (.525mm)	
REGULAR SINGLE-SHEAR	L-541	L-542	K-96	GOLD /.027" (.675mm)	
REGULAR DOUBLE-SHEAR	L-551	L-552	K-96	GOLD /.027" (.675mm)	

TMS LINK SERIES FEATURES

- TMS LINK Series Pins and corresponding KODEX Drills are color-coordinated for quick and easy identification.
- Resilient plastic shanks conveniently fit any latch-type handpiece. No chucks necessary!

- Plastic shanks align TMS LINK Series Pins automatically to the drilled channel for smooth and safe pin insertion.
- Attractive, see-through plastic flip-top package permits easy access for quick pin removal.

Coltène/Whaledent Inc.
 235 Ascot Parkway
 Cuyahoga Falls, OH 44223 / USA
 Tel. USA & Canada 1 800 221 3046
 +1 330 916 8800
 Fax +1 330 645 8704
 info.us@coltene.com

Coltène/Whaledent AG
 Feldwiesenstrasse 20
 9450 Altstätten / Switzerland
 Tel. +41 (0) 71 757 53 00
 Fax +41 (0) 71 757 53 01
 info.ch@coltene.com

Coltène/Whaledent GmbH + Co. KG
 Raiffeisenstrasse 30
 89129 Langenau / Germany
 Tel. +49 (0) 7345 805 0
 Fax +49 (0) 7345 805 201
 info.de@coltene.com

©02/03/2019

Made in USA

P/N 30002426M

EN**SUGGESTED PIN PLACEMENT PROCEDURE**

1. Determine pin positions and indent channel sites slightly with No. 1/2 or No. 1/4 round bur (Fig. 1).

NOTE: Channel sites should be at least 0.5mm from the dento-enamel junction.

2. Drill pin channel in dentin with KODEX® drill, making no more than two passes. A Reduction Gear Handpiece will optimize channel sizing and pin placement. Engine must be engaged clockwise during insertion and removal of twist drill (Fig. 2).

• For MINIKIN™ (.017") and MINUTA® (.0135") pins; drill channel in one smooth pass; to full depth and out.

• For REGULAR (.027") and MINIM (.021") pins; using a depth limiting drill, drill channel in two passes; first to half depth and out, then to full depth and out.

Inspect KODEX drills after repeated use; blunt tip may cause oversizing of channel.

3. Insert TMS LINK® Series Pin into a latch-type reduction gear handpiece and position pin over channel before activating engine for insertion (Fig. 3).

4. Operate engine at low speed (750-1000 RPM) and apply light downward pressure. The TMS LINK Series Pin will thread into drilled channel and shear when seated. (Fig. 4, placement of two Double Shear Pins) NOTE: The optional universal hand driver may be used to install TMS Link pins manually.

5. TMS Bending Tool may be used to bend pins to match the contours of the final restoration. Build up superstructure with material of choice (Fig. 5).

CLEANING, DISINFECTION AND STERILIZATION

Drills, TMS Bending Tool and Hand Wrench: Prior to each use, remove debris ultrasonically and sterilize by autoclave ONLY between patient use. (To minimize corrosion, place drills in separate bags from other instruments). Immediately air blow dry. A disinfectant may be used prior to use.

Pins: Prior to each use, place pins with sterilization indicator, into an autoclavable bag and seal. Sterilize in an autoclave at 130-135 degrees C for 12-15 minutes, following standard procedures.

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