
# VECTOR® 900

## Gebrauchshinweise/Hints for use

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
**Description:** Special screw with integrated bars for protrusive double plates, comes with key for expansion screw. Screw elements are made from stainless steel DIN 1.4305 and bars from MIM 17.4 PH.

**Indication:** The VECTOR® 900 screw is indicated for protrusive double plates for the treatment of Angle Class II cases (lower jaw retrusion/distal occlusion) in growing patients. Can be used at the same time for bilateral expansion of the upper jaw. For the lower jaw, we recommend a combination with the VECTOR® 100 expansion screw (REF 2492). The angle tester acc. to Adamik (REF 2332) is ideally suited to determine and control the correct insertion position of the screw (60° - 70° in relation the occlusal plane).

**Contraindications:** The stainless steel alloy used is highly corrosion-resistant and suited for oral use. However, in patients hyper sensitive to nickel, the minor nickel percentage may cause allergic reactions. If the allergic reaction is confirmed by the corresponding allergy test, it may be necessary to discontinue treatment.

**Cleaning/care:** CETRON® cleaning powder is best suited for cleaning the screw. CETRON® spray can be used additionally for the care of the screw thread.

### Dimensions and expansion parameters:

W x H x L: 14.0 x 20.0 x 13.0 mm  
Frontal expansion per 1/1 rotation: 0.7 mm  
Frontal expansion per 1/4 rotation: 0.175 mm  
Total frontal expansion: 5.0 mm

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