## VITA AKZENT® LC

Gebrauchsanweisung | Kurzversion
Instructions for use | short version
Mode d’emploi | Version courte
Instrucciones de uso | Versión abreviada
Istruzioni d’uso | -versione breve

## VITA – perfect match. VITA

## VITA

VITA Zahnfabrik · H. Rauter GmbH & Co. KG
Spitalgasse 3 · 79713 Bad Säckingen
Germany · Allemagne · Alemania · Germania
Tel. +49 (0) 7761 / 562-0
Fax +49 (0) 7761 / 562-299
www.vita-zahnfabrik.com
info@vita-zahnfabrik.com

922-12410 – 0421 (X.)TPP – Version (02)

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**Lichthärtende Malfarben**

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**Instructions for use** **EN**

**Light-curing stains**

Please read the instructions for use before using the product.

**Intended use**

VITA AKZENT LC is a light-curing methacrylate-based stain/glaze system for extraoral surface characterization of dental restorations made of hybrid ceramic, resin veneering materi-als, CAD/CAM composites, prefabricated teeth and denture base resins. It can also be used for internal characterization with the layering technique of veneering composites.

**Patient target group**

No restrictions

**Indication**

Indication range:

- Restorations made of hybrid ceramic (VITA ENAMIC)
- Restorations made of light-curing veneering material (e.g., VITA VM LC)
- Restorations made of CAD/CAM composites (e.g., VITA CAD-Temp)
- Prefabricated teeth (e.g., VITAPAN)
- Denture bases (e.g., VITA VIONIC BASE)

**Contraindication**

- In patients with allergies or sensitivities to the ingredients
- Do not use on occlusal contact points of restorations
- Not approved for intraoral use

**Product application/control**

**Workflow:**

**Hybrid ceramic:**

- Surface pretreatment
- Surface cleaning
- Silanization
- External characterization
- Glaze application
- Polymerization

**Veneering composite:**

- Surface pretreatment
- Application of modelling liquid (only for internal characteri-zation and if no inhibition layer is present).
- Internal characterization
- Application of veneering resin
- External characterization
- Glaze application
- Final polymerization

**CAD/CAM composites/prefabricated teeth/denture base resins:**

- Surface pretreatment
- Surface cleaning
- Silanization with highly filled CAD/CAM composites, such as Cerasmart (GC)
- External characterization
- Glaze application
- Polymerization

Please note the recommended processing time of approx. three minutes and clean the paint brushes with VITA AKZENT LC CLEANER after use.

**Polymerization**

The varnish can be polymerized with many standard dental light-curing units with a spectral range of 350 nm – 500 nm. You can find the latest information on the recommended poly-merization devices at: www.vita-zahnfabrik.com/akzentlc

**Recommended product combinations**

VITA AKZENT LC can be mixed with VITA VM LC flow at a max. ratio of 1:10.

**Storage/disposal**

Store at 4 °C - 25 °C (39 °F - 77 °F). Storage in a refrigerator is recommended.

Do not use after the expiry date!

The products labelled with a pictogram for hazardous sub-stances are to be disposed of as hazardous waste. All other products are to be disposed of as mixed industrial waste.

| Chemical composition                                  | Wt%     |
|---|---------|
| Methyl methacrylate and multifunctional methacrylates | 30 – 40 |
| Urethane(meth-)acrylates                              | 40 – 60 |
| Silicon dioxide                                       | 8 – 11  |
| Ethyl phenyl(2,4,6-trimethylbenzoyl)phosphinate       | 2 – 6   |
| Other   | < 1     |
| Pigments  | < 2     |

**General risks/reporting incidents**

Information on the general risks associated with dental treat-ments and on reporting serious incidents in connection with medical devices can be found at www.vita-zahnfabrik.com/service\_risks Safety data sheets can be downloaded at www.vita-zahnfabrik.com/vita-akzent-lc.

**Hazardous substances**

Danger

Contains methyl methacrylate, two propenoic acid reaction products with pentaerythritol, various acrylates
H317 May cause an allergic skin reaction
H 318 Causes serious eye damage
H 412 Harmful to aquatic life with long-lasting effects

**Disclaimer**

Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in con-junction with materials or equipment from other manufactur-ers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not a required component of the product. Publication of these work-ing instructions: 04.21

After the publication of these working instructions, any previ-ous versions become obsolete. The current version can be found under www.vita-zahnfabrik.com

VITA Zahnfabrik has been certified and the following products bear the mark:

**VITA AKZENT® LC**

**CE0124** **Rx only** **MD**