## Rinn®

### XCP®

## XCP-ORA®

## XCP-DS®

## XCP-DS FIT®

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**ENGLISH**

**Indications for use:** XCP® holders are indicated for holding film and phosphor plates for intraoral dental radiographs. XCP® and XCP-ORA® arms and rings are indicated for positioning and aligning dental intraoral radiographs. XCP-DS® and XCP-DS FIT® holders are indicated for holding digital sensors for intraoral dental radiographs.

**Contraindications:** None known.

**Warnings: Sterilize device before use.** DO NOT use phenol-based glutaraldehyde or ultrasonic cleaners. DO NOT chemiclave, dry-heat sterilize or cold sterilize. Only sterilize with steam sterilization according to the provided sterilization instructions. DO NOT exceed 273°F (134°C). It is the user’s responsibility to qualify alternative methods of disinfection and sterilization. Failure to follow properly validated sterilization processes and approved aseptic techniques may result in cross contamination. Any method of sterilization will shorten the life of plastic parts.

**Precautions:** Follow sterilization instructions carefully to prevent plastic damage during sterilization. When using with corded digital sensors, use a sensor cover. Do not pull on cord during use.

**Adverse Reactions:** None known.

**Residual Risks:** Failure to follow the validated reprocessing instructions provided may result in exposing the patient to cytotoxic residuals and/or cross-contamination.

**Step-by-Step Instructions:**

- Insert positioning arm into the appropriate color-coded channel of the aiming ring.
- Attach appropriate biteblock for area to be imaged. Correctly assembled, the imaging area should be centered when viewed through the aiming ring.
- Position holder in patient’s mouth and slide ring close to patient’s face.
- Align x-ray tube to ring.
- Take exposure.
- Adjust as needed for the next desired position.

**Reprocessing:**

Arm and Ring: Disassemble instrument(s) and separate all parts prior to reprocessing.

XCP-DS Fit Biteblocks: Separate the silicone band from the biteblock and retaining clip prongs, leaving the band loosely attached to the biteblock.

**Cleaning: Manual / Mechanical (Ultrasonic)**

- Rinse the disassembled device under potable running water while brushing with a soft-bristled brush to remove gross soil.
- Prepare an enzymatic instrument cleaning solution and immerse the device per the manufacturer’s instructions.
- Rinse the device under potable running water for at least 30 seconds and dry using a disposable lint-free cloth.
- Visually inspect the device for visible soil. If soil is seen, repeat the cleaning procedure.

**Cleaning and Disinfection: Automated – Instrument Washer-Disinfector**

Use only properly maintained, calibrated, and approved washer-disinfectors according to ISO 15883.

- Rinse the disassembled device under potable running water while brushing with a soft-bristled brush to remove gross soil.
- Follow manufacturer instructions for instrument washer solution, observing concentration rates and contact times as well as the instrument washer’s recommended cleaning and/or thermal disinfection cycle (1 minute at 90°C).
- Visually inspect the devices for visible soil and repeat the cleaning and disinfection procedure if necessary.

**Sterilization:** Ensure that the device is thoroughly dry prior to sterilization. Enclose the device into an appropriately sized sterilization pouch intended for steam sterilization. Place the pouched device into the sterilizer per the manufacturer’s instructions. Sterilize the device according to one of the recommended methods listed below.

**Sterilize metal and plastic components in separate pouches to avoid damage to the components. DO NOT exceed 273°F (134°C).**

Steam Sterilization Cycle	Minimum Temperature	Minimum Full Cycle Exposure Time	Minimum Dry Time
Pre-vacuum	270°F (132°C)	3 min	16 min
Steam-Flush Pressure-Pulse	270°F (132°C)	3 min	16 min
Gravity Steam Sterilization	270°F (132°C)	10 min	15 min

**Alternate Method:** Place the non-pouched device into the sterilizer and sterilize according to one of methods listed above. Non-pouched device must be used immediately following sterilization.

**Storage:** To maintain sterility, the device should remain in the sterilization pouch until ready for use.

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the Competent Authority of the Member State.

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