SDI	AURA EASYFLOW	E N G L I S H
	Ultra Universal	
	Restorative Material	

INSTRUCTIONS FOR USE

Aura Easyflow flowable composite, is the ideal versatile radiopaque, flowable, light cured material. The composite is directly injected into the cavity preparation for maximizing adaptation to the preparation.

COMPOSITION:
44% wt (68% vol) multifunctional methacrylic ester
56% wt (32% vol.) inorganic filler (0.2 - 1 micron)

INDICATIONS:

- Conservative Class I (including occlusal surfaces), II, III, IV restorations
- Class V restorations (cervical abrasion, root erosion, abfraction)
- Repair of defects in enamel and porcelain restorations
- Sealants (pit & fissure, implant)
- Radiopaque cavity liner under direct restorations
- Blocking out of undercuts
- Minor core build ups
- Sealing of temporary crowns
- Cementation of porcelain, ceramic, composite veneers
- Spinting (mobile teeth, fibre bridges)
- Cover stains

CONTRAINDICATIONS:

- Pulp capping
- Do not use in conjunction with any eugenol containing materials
- Anyone with allergy to acrylics
- Do not use if a dry field cannot be established

INSTRUCTIONS:

Clean and isolate tooth. A conservative cavity should be prepared, employing standard techniques and instruments, to form a slightly rounded internal form. Margins of the cavity preparation should end in sound and supported enamel with no levels in stress bearing locations. If beveling is desired in a low stress location it should extend no further than 1mm at an angle of no greater than 45°. Pre-wedging is also recommended to ensure that the restored tooth will have an adequate contact point. Prophy all surfaces to be etched including surfaces adjacent to the cavity with an oil free non-fluoride containing paste or a slurry of pumice and water. Rinse thoroughly with water.

Shade Selection

- Do not over dry the tooth before shade matching. Only use clean dry tooth for colour matching.
- Use the Aura Easy shade guide to select your shade before placement of rubber dam.
- When the adjacent teeth are significantly different in shade, try to match the lighter shade.
- For visual confirmation of shade selected, place small amounts of composite on a surrounding tooth surface and light cure before bonding.

Approximate equivalence to VITA® classical A1-B4* shades.

A1/B1 ≈ ae1
A2/B2 ≈ ae2
A3/B3 ≈ ae3
A3.5/B4 ≈ ae4

Isolation

Isolation techniques must be used to prevent contamination. Rubber dam is the preferred mode of isolation.

Pulp Protection

For deep cavities an appropriate liner or cement should be placed at the deepest point of the cavity.

- Acid Etching

Thoroughly dry the surface to be etched with dry, oil-free air. Etch tooth surface with 37% phosphoric acid

(a) Enamel only: Etch surface for at least 20 seconds.

(b) Dentin and enamel:

Using the "total etch" technique, etch the surface including any glass ionomer for at least 20 seconds.

Enamel subjected to fluoridation should be etched for 90 to 120 seconds.

Etching Precautions: Ensure that the dispensing tip hub is firmly attached to the syringe by twisting the hub securely onto the syringe. Avoid acid contact with oral tissues, eyes and skin. If accidental contact occurs wash thoroughly with water. In the case of eye contact, wash eye for 15 minutes and seek medical attention. Use matrix strips to protect adjacent tooth surfaces during etching. On first usage or after prolonged storage, extrude a small amount onto a mixing pad for familiarity with the etchant's viscosity and rate of extrusion.

 - Wash thoroughly with water.
 - Remove excess water. Keep moist. Avoid contamination e.g. saliva.
 - Apply bonding agent to saturate all internal surfaces according to manufacturer's instructions.
 - Directly inject the flowable composite in increments of 2mm or less in:
 - Class V restorations,
 - Pit and fissure sealants,
 - Conservative Class I, II, III and IV restorations
 - Other indications as required.

CAUTION: ensure that the dispensing tip is **firmly** attached to the syringe by twisting the tip hub securely onto the syringe.

WARNING: Apply restorative at normal room temperature (23°C/74°F).
- Light cure the composite using a suitable high power LED curing light (460-480nm wavelength) for a minimum of 20 seconds in increments of 2mm.
- Polishing of composite

(a) Remove excess composite and contour desired shape using a fine diamond or a 12-fluted carbide bur.

Tips for polishing using a flexible disc system:

 - The polishing motion **should be constant and unidirectional**.
 - A back and forth movement over the composite - enamel margin is not recommended.
 - Keep the surface and polishing disc dry while polishing. A dry surface will produce a smoother, more uniform finish.
 - Do not use a handpiece with a speed greater than 35,000 rpm.
 - Avoid touching the composite with the mandrel or disc eyelet.

(b) For gross reduction, set the speed of the handpiece to approximately 10,000 rpm. Use the coarsest grit disc of the range and polish the restoration starting at the gingival and **moving outward** over the restoration. Replace the disc if cutting efficiency impaired. Wash and dry.

(c) For final contouring, use a medium coarse grit at a speed of approximately 10,000 rpm. Wash and dry.

(d) For finishing, use a medium fine grit disc at a speed of approximately 30,000 rpm.

Wash and dry.

(e) Finally, use the finest grit disc at 30,000 rpm. Wash and dry. Note: After step 7, polishing paste may be used to give the composite a lustrous finish, using the following steps:

 - Apply polishing paste to a rubber cup or disc.
 - Smear a thin layer of paste onto the restoration.
 - Add a small amount of water to both the tooth and cup. Polish for 30 seconds at low speed and light pressure.
 - Wash and dry.

CEMENTATION OF PORCELAIN, CERAMIC AND COMPOSITE VENEERS

Prior to cementation with Aura Easyflow, the restoration must be pre-treated according to the manufacturer's instructions for use.

Aura Easyflow is only suitable for the cementation of restorations that have translucency and are of appropriate thickness to allow polymerization light to pass through, to ensure that the light reaches Aura Easyflow, so that it can cure completely.

PRECAUTIONS:

Avoid prolonged contact of bonding agents and composite with the skin or oral soft tissue, as it may cause inflammation of the oral soft tissues or skin sensitization.

Keep out of reach of children.

Do not take internally.

Store at temperatures between 2° and 25°C (35.6°-77°F). When not in use replace the cap tightly. Use at room temperature. Do not use after expiry date. Do not expose material to direct light. Caution: Federal Law restricts this device to sale by on the order of a dentist.

FIRST AID

- Eye (contact):** Wash thoroughly with water and seek medical attention.
- Skin (contact):** Remove using a cloth or sponge soaked in alcohol. Wash thoroughly with water.
- Ingestion:** Rinse mouth thoroughly. Drink plenty of water/ milk. Seek medical attention if symptoms persist.
- Inhalation:** No symptoms expected.

<div></div> <p>Made in Australia by SDI Limited 3-15 Brunson Street, Bayswater, Victoria 3153 Australia</p> <p>Australia 1 800 337 003 Austria 00800 0225 5734 Brazil 0800 770 1735 France 00800 0225 5734 Germany 0800 1005759 Italy 00800 0225 5734 New Zealand 0800 734 034 Spain 00800 0225 5734 United Kingdom 00800 0225 5734 USA & Canada 1 800 228 5166 www.sdi.com.au</p>	<div>EC REP</div> <p>SDI Germany GmbH Hansestrasse 85 51149 Cologne Germany T+49 2203 9255 0</p>
	<div></div> <div></div>

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