EN Instructions for use MD EU Medical device

Material-related characteristics

Chemical composition [wt. %]

Table with 2 columns: Component and value. Rows: Polymethylmethacrylate (PMMA) > 99, Coloring pigments < 1

Physical specifications

Table with 4 columns: Property, Value 1, Value 2, Value 3. Rows: Flexural strength, Fracture toughness, Density, Water sorption, Solubility, Residual monomer content

* determined by 3-point flexural test following ISO 10477

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CediTEC DB is a pre-colored dental milling blank made of impact resistant PMMA for the manufacture of denture bases for removable dentures.

2. Contraindications

- Intolerance to the contained components
Crowns, bridges, model casting with clasps, bars, narrow sublingual bars and transversal connectors

3. Safety information

Please pay attention to the information in the material safety data sheet in its current version. Avoid the inhalation of dust particles during processing.

4. Handling and storage

Only use clean, defect-free material from undamaged packaging and make sure that the blanks are stored at room temperature, in their original carton, dry and protected from light.

5. Construction

This medical device may only be processed by trained personnel. The denture base should be positioned centrally in the milling blank and should not be less than a minimum wall thickness of 1 mm.

6. Milling and further processing

Milling must be carried out with a milling strategy adapted to PMMA or milling parameters for CediTEC DB.

Only use one-edged milling tools in the milling machine or wet grinding procedures!

Further processing:

Carbide cutters with FSQ toothing and black-ring cutters for titanium are suitable for roughing. The separation of the denture bases from the milling blank and the fine finishing can be performed with simple, cross-toothed carbide millers.

7. Preparation of bonding surfaces

Check the fit of the prosthetic teeth before bonding. In case of interfering contacts, make the correction on the tooth. For an optimum bond, roughen the bonding surfaces of the restoration with aluminum oxide (1 - 2 bar / 50 to 125 μm).

8. Individualisation, relining and repair

Individualisations, relinings and repairs of CediTEC DB can be carried out with conventional auto-polymerizing denture base resin / light-curing composites.

9. Cleaning

- Ultrasonic bath (max. 40 °C / 104 °F) with mild cleaning agent

Do not use the following methods for cleaning:

- Steam blasting
Autoclave method (e.g. sterilization / disinfection)
Organic solvents and acids (e.g. acetone, hydrochloric acid)
Highly concentrated alcohols
Chemical denture cleaners
Abrasive cleaning agents

10. Disposal

Remains of PMMA blanks can be disposed of in the regular household garbage.

Explanation of symbols

Grid of symbols and their meanings: Manufacturer, Date of manufacture, Use-by date, Batch Code, Catalog number, Medical Device, Keep dry, Keep away from sunlight, Recommended storage temperature range, Consult instructions for use

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