**DC DENTAL CENTRAL**

# DC Ätzgel / DC Ätzgel Maxi

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Stand 08 / 2019

## GB Instruction for use

36 % phosphoric acid etching gel with optimum thixotropic flow characteristics. It serves as an etching of tooth substance for enamel and total etch technic. Taking the safety and application instructions into consideration, there are no known side effects to date. Could remove under filling material by prolonged usage of etch gel.

### Application:

#### Filling up of empty applicator syringe (only DC Etching Gel Maxi):

Empty applicator syringes are single used and must not be reprocessed. Only use empty unused applicator syringes for filling up!

The empty applicator syringe should be locked onto the adapter of the supply syringe (180° clockwise). Fill it up by carefully and slowly pressing the plunger. Eventually pull simultaneously the plunger and the syringe, paying attention not to press the plunger off the syringe. Lastly turn the adapter to remove it from the applicator syringe and the 50 ml syringe should be tightly closed. Fix the canule onto the applicator syringe.

### Preparation of the applicator syringe:

Take the applicator syringe out of the packaging and remove the cap. An application canule should be screwed onto the syringe (clockwise). It is now ready to be used by carefully pressing the syringe.

### Preparation:

The area of the tooth substance to be etched should be clearly isolated.

- 1) The use of a coffer dam is recommended. Any contact with soft tissues (skin, oral mucosa, eyes, etc.) is to be avoided. In case of contact with soft tissues, immediately rinse with plenty of water and seek medical advise.
- 2) Before treatment, clean tooth with fluoride and grease free polishing paste with the help of a rubber cup.
- 3) Remove cleaning residues with water and dry the surface with air pressure.
- 4) Prepare the respective cavity, chamfer edges to achieve permanent filling edges. Rinse the cavity with water and dry with air pressure. Avoid overdrying.
- 5) Protect the dental pulp in the cavity (e.g. cover with a thin layer of underfilling). The cavity should be dry and free from impurities but avoid overdrying the dentins.

### Doctors process:

- 6) Apply the Etching gel on enamel edges of the cavity and when necessary on the dentin.

#### Recommended resting time:

15 - 30 sec. on the enamel

10 - 15 sec. on the dentin

50 - 80 sec. on milk tooth and tooth with high fluoride content

- 7) Minimum 20 secs. rinsing with water.
- 8) Dry with oil and water free air pressure. The etched surface is not to be disturbed until the next process or contaminated by saliva. In case of contamination, repeat the etching process (approx. 10 sec.).

The etched enamel should be chalk white!

After etching it can be proceeded applying adhesives etc., according to the manufacturers instructions for use.

### Safety Instructions:

H314: **Causes severe skin burns and eye damage. Protective gloves and eye protection should be worn.**

P301+P330+P331: IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P302+P352: IF ON SKIN: Wash with soap and water.

P305+P351+P338: IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. P310: Immediately call a POISON CENTER or doctor/physician.

P404: Store in a closed container.

Disposal must be made according to official regulations.

Content: 36 % ortho-phosphoric acid

- Only for use in Dentistry!
- Not to be used after expiry date
- Syringe must be immediately closed after use.
- Keep out of reach of children.
- Store at max. temperature at 25° C.