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Light-curing resin based on (meth)acrylate, for the generative fabrication of flexible gingiva masks for DLP printers with UV-LED 385 nm

Suitable for the following DLP-Printer

Asiga 385 nm	MiiCraft 385 nm	Rapidshape 385 nm
MAX UV Pro2 UV Pico2 UV	MiiCraft y-Serie	DII-Serie

Printers may only be operated using material parameters authorised by DETAX.

Important notes

Only for dental use by qualified personnel.

Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description on page 2).
- ▶ Homogenize/roll the material prior to processing, i.e. with a roller mixer.
- ▶ Processing temperature 23 °C ± 2 °C.

Safety

- ▶ Please follow the instructions on the safety data sheet!
- ▶ Be sure to use personal protective equipment (protective gloves and protective glasses) during processing.
- ▶ Avoid direct contact with the liquid material and the components prior to post-curing. Irritating to eyes and skin (sensitisation is possible).
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.

Storage

- ▶ **FREEPRINT® gingiva** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

General

- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.

Contraindication

Contains (meth)acrylics and phosphine oxide.

Some ingredients of **FREEPRINT® gingiva** may cause allergic reactions in predisposed persons. In such cases refrain from using the product.

Adverse effects

Product may cause allergic reactions.


Indication:

Gingival masks for dental 3D models

Processing:

at 23 °C ± 2 °C

Storage

Ordering information:

FREEPRINT® gingiva 385
500 g bottle, gingiva **02820**
1.000 g bottle, gingiva **02843**

FREEPRINT® model 385
1.000 g bottle
sand **03778**
ivory **03780**
grey **03782**
caramel **02850**

FREEPRINT® model 2.0 385
1.000 g bottle
sand **02128**
light grey **02099**
grey **02177**
caramel **02850**
white **02148**

Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath.

Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath. Prior to post-curing, check for residues.

Post-curing

Post-curing is performed with a xenon photoflash unit (e.g. Otofash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

FREEPRINT® gingiva was developed for use in the dental field and must be used in accordance with the instructions for processing and safety. DETAX will not be responsible for damages caused by faulty or improper use of system and materials.