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Light-curing resin based on (meth)acrylate, biocompatible, for the generative fabrication of temporary restorations of anterior and posterior teeth for DLP printers with UV-LED 385 nm

Suitable for the following DLP-Printer

Asiga 385 nm	MiiCraft 385 nm	Rapidshape 385 nm	W2P 385 nm
MAX UV, Pro2 UV, Pico2 UV	MiiCraft y-Serie	DII-Serie	Solflex Serie

Printers may only be operated using material parameters authorised by DETAX.

Important notes

This is a medical device, only to be used by trained specialist personnel.

Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe process description on page 2).
- ▶ Homogenize/roll the material prior to processing, i.e. with a roller mixer.
- ▶ **Maximum curing depth* at direct post-exposure:** 3 mm
*In case of large objects and exposure on both sides, the material thickness can be up to 6 mm (Example FREEPRINT® temp – with a curing depth of 3 mm).
Do not use heat-based methods for disinfection or sterilisation. This could possibly deform the workpiece.
- ▶ **FREEPRINT® temp** is applied for the production of temporary c & b in the anterior and posterior area and can be used for single crowns and bridges up to 4 elements containing one pontic. The occlusal minimum thickness amounts to 1.5 mm, circularly up to 0.8 mm. The connector cross-section for pontics in the anterior area is 12 mm² minimum, in the posterior area 15 mm². Any non-compliance of these requirements can lead to an undesirable outcome.
- ▶ 25 and 50 µm are the layer thickness for the printing process.
- ▶ **FREEPRINT® temp** can be shaded individually by means of our smartrepair® System,
- ▶ Polish the surface mechanically. Prepolish by means of rotating brushes and prepolishing paste, highgloss by means of buffing wheels and highshine polishing paste for resins.
- ▶ The temporaries can be cemented by means of usual temporary cements, such as tempolink® (eugenol-free). Temporary cements containing eugenol must be removed completely after the removal of the provisional as otherwise the setting of subsequently used permanent cements could be affected.
- ▶ Processing temperature 23 °C ± 2 °C.

Safety

- ▶ Please follow the instructions on the safety data sheet!
- ▶ Be sure to use personal protective equipment (protective gloves and protective glasses) during processing.
- ▶ Avoid direct contact with the liquid material and the components prior to post-curing. Irritating to eyes and skin (sensitisation is possible).
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.
- ▶ Biocompatibility is only guaranteed with complete polymerisation.

Storage

- ▶ **FREEPRINT® temp** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

General

- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.

Contraindication

Contains (meth)acrylics and phosphine oxides.

Some ingredients of **FREEPRINT® temp** may cause allergic reactions in predisposed persons. In such cases refrain from using the product. **FREEPRINT® temp** only insert intraorally in completely polymerised state.

Adverse effects

Product may cause allergic reactions.

Indication:

Temporary crowns and bridges

Processing:

at 23 °C ± 2 °C

Storage

15 °C
59 °F

28 °C
82 °F

Ordering information:

FREEPRINT® temp 385

500 g bottle
A1 **04058**
A2 **04059**
A3 **04060**

1.000 g bottle
A1 **04062**
A2 **04063**
A3 **04064**

FREEPRINT® denture 385

rose transparent
500 g bottle **02060**
1.000 g bottle **02040**

FREEPRINT® tray 385

1.000 g bottle, green **04086**

FREEPRINT® ortho 385

1.000 g bottle, clear **03989**

FREEPRINT® splint 2.0 385

500 g bottle, clear **02080**
1.000 g bottle, clear **02076**

FREEPRINT® IBT 385

500 g bottle, clear **04248**
1.000 g bottle, clear **04249**

Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath.

Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath. Prior to post-exposure, check the openings and additional bore holes for residues.

Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otofash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

Surface processing

Polish surface mechanically.

Cementation

Temporary cementation.

FREEPRINT® temp was developed for use in the dental field and must be used in accordance with the instructions for processing and safety. DETAX will not be responsible for damages caused by faulty or improper use of system and materials.