||BOND® Etch 20 Gel||BOND® Etch 35 Gel||BOND® Etch 35 Fluid

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GB Instructions for use

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SE Bruksanvisning

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NO Bruksanvisning

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iBond® Etch 20 Gel, iBond® Etch 35 Gel, iBond® Etch 35 Fluid

Product description

iBond® Etch Gel and Fluid are etching agents based on orthophosphoric acid. They are used to etch enamel and condition dentine in the context of adhesive filling therapy.

Composition

iBond® Etch 20 Gel contains 20 % phosphoric acid
iBond® Etch 35 Gel/Fluid contain 35 % phosphoric acid
The products further contain Thickening agents, Pigments, Water

Indications

Enamel etching and enamel/dentine etching in the context of the total etch technique prior to the adhesive bonding of:

• Composite and compomer restorations • Laboratory-fabricated, indirect restorations (e.g. inlays, crowns, bridges, veneers) • Sealants

Contraindications

• Application to dentine close to the pulp (less than 1 mm) • Do not use in case of known hypersensitivity to iBond® Etch

Side effects

This product or one of its components may, in particular cases, cause allergic reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

Interactions

Liner and base materials, such as calcium hydroxide, may be dissolved by the acid, depending on their type and composition.

Warnings/precautions

iBond® Etch 20 Gel contains 20 % phosphoric acid and is an irritant. iBond Etch 35 Gel and Fluid 35 contain 35 % phosphoric acid and can cause chemical burns. May irritate eyes. Avoid contact with gingiva, skin, eyes or neighbouring teeth. In case of contact with oral mucosa, rinse thoroughly with plenty of water. If swallowed, thoroughly rinse mouth with water and then drink plenty of water. Consult a physician. In case of contact with skin, thoroughly rinse with soap and water. In the event of contact with eyes, rinse thoroughly with water and seek medical advice. Wear suitable protective gloves. Avoid skin contact. Wear protective clothing. Wear safety glasses/face protection. Immediately remove clothing contaminated with iBond® Etch 35 Gel or iBond® Etch 35 Fluid.

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Application

Prior to treatment, carefully clean tooth with fluoride-free polishing paste. The use of a rubber dam is recommended.

1. After preparation, rinse the cavity in accordance with adhesive filling therapy and dry with oil-free air.

2. Starting from the bevelled enamel areas, apply iBond® Etch to the entire cavity surface, including the dentine. Ensure that the enamel is etched for at least 20 s and no more than 30 s, and that the dentine is conditioned for at least 15 s and no more than 20 s. When using iBond® Etch 35 Gel/Fluid, pay special attention to complying with the dwell time on the dentine. Reducing the time to 15 s in these areas is permissible to avoid overetching. If desired, only apply iBond® Etch to the enamel and let set for a maximum of 30 s. Etch unprepared enamel, e.g. in case of fissure sealing, for at least 30 s and no more than 60 s. In case of fluoridated enamel, etching for no more than 60 s may also be required.

Caution: In case of deep cavities, we recommend covering the dentine near the pulp with a calcium hydroxide preparation and a base (e.g. glass ionomer cement) prior to etching. 3. Use a water sprayer to completely rinse off iBond Etch for 20 s. The cavity must be completely free of any remnants.

4. Dry the cavity with oil-free air (do not dry out dentine!). The etched enamel surface should appear chalky white. In case of contamination of the tooth surface (e.g. with saliva), the etching must be repeated (for a maximum of 10 s dwell time).

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