




# neodisher® MediZym



pH-neutral detergent for the reprocessing of thermostable and thermolabile instruments



## Liquid concentrate

### Main fields of application:

- Automated cleaning of surgical instruments, anaesthesia equipment, flexible endoscopes, containers and other medical equipment

### Performance spectrum:

- Removes residues of dried and denatured blood
- Stainless steel, instrument steel, light metals, glass, the usual plastics and the materials of anaesthesia equipment and flexible endoscopes are compatible with working solutions of neodisher MediZym

### Special properties:

- High degree of material protection
- pH neutral basis with enzymes and other cleaning adjuvants

### Application and dosage:

neodisher MediZym can be used in washer disinfectors as well as in immersion and ultrasonic baths. The dosing amount depends among others on the respective field of application and the individual degree of contamination of the instruments to be reprocessed. The following parameters are recommended when using neodisher MediZym.

Automated cleaning of thermostable and thermolabile instruments	2 - 10 ml/l (0.2 - 1.0 %), 40 - 50 °C, 5 - 10 min
---	--

Automated cleaning of flexible endoscopes endoscopes	5 ml/l (0.5 %), 5 - 10 min
Manual cleaning of thermostable and thermolabile instruments in immersion baths	5 - 20 ml/l (0.5 - 2.0 %), max. 45 °C, 20 - 30 min
Manual cleaning of thermostable and thermolabile instruments in ultrasonic baths	5 - 20 ml/l (0.5 - 2 %), max. 45 °C, 5 - 10 min

The working solution of neodisher MediZym has to be replaced at least daily, immediately in case of visible soiling.

The neodisher MediZym working solution has to be rinsed off completely (preferably with deionised water).

To avoid water stains, the use of fully demineralised water in the final rinse is recommended.

### Main fields of application:

- For professional use only.
- Do not mix with other products
- Rinse out dosing system including suction hose with water before changing product.
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes.
- Please observe the medical device manufacturer's recommendations for use according to the requirements of DIN EN ISO 17664.



# neodisher® MediZym

- The instructions given by the manufacturer of the washer disinfectant are to be observed.

MB 4040/3-2  
Date of issue: 05/2020

## Technical data:

pH range	7.6 – 7.7 (2 - 20 ml/l, determined in deionised water, 20 °C)
Density	approx. 1.0 g/cm <sup>3</sup> (20 °C)
Viscosity	< 10 mPas (concentrate, 20 °C)

## Ingredients:

Ingredients for detergents according to EC detergent regulation 648/2004.

< 5% anionic surfactants

also: enzymes, preservatives (2-phenoxy-ethanol, alkyl 4-hydroxybenzoate)

## CE-marks:


neodisher MediZym complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

## Storage information:

Store in a cool, frost-free place. Always store at a temperature between 0 °C and 25 °C.

Keep away from sunlight. Usable for 2 years when stored as recommended.

For expiry date refer to the stamp mark on the label behind the hourglass symbol .

## Hazard and precautionary statements:

For safety information see Safety Data Sheets.

These are available e.g. at [www.drweigert.com](http://www.drweigert.com) under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.