STATMATIC smart 1, 2, 3

Cleaning and Maintenance System

Instructions For Use

EN/DE

Brugsvejledning DK
Manuel d'instructions FR
Istruzioni per l'uso IT
Gebruiksaanwijzing NL
Instrucciones de uso ES

via: www.scican.uk.com/downloads
via: www.scican.uk.com/downloads
via: www.scican.uk.com/downloads
via: www.scican.uk.com/downloads
via: www.scican.uk.com/downloads

SciCan

Your Infection Control Specialist™

www.scican.com

STATMATIC smart 1, 2, 3
Copyright 2017 SciCan. All rights reserved.

Doc. No. 2.001.8462 – Rev. 2017-11

Table of Contents

1	Important Information.....	3	"Special instrument" program.....	15
	Type description.....	3	4 Maintenance.....	16
	Scope of delivery.....	3	Checking adapter fixation.....	16
	General.....	3	Checking the O-rings.....	16
	Intended use.....	4	Functional testing of flow rate using the testing	
	Disposal.....	4	adapter.....	16
	Qualifications of the operator.....	5	Changing the adapter.....	16
	Safety Notes SciCan "Dental Oil" and SciCan		Refilling of the media containers.....	17
	"Cleaner".....	5	Replacing the absorption pads.....	18
	Used symbols.....	5	Empty the collecting tray.....	18
	Disclaimer.....	6	Inside cleaning and disinfection of the maintenance	
	Device overview.....	7	area.....	18
2	Commissioning.....	9	Outside cleaning and outside disinfection of the	
	Environmental conditions.....	9	device.....	18
	Electrical Requirements.....	9	Recommendation.....	18
	Media connections.....	9	Note for users in countries with a legal validation	
	Filling the media system.....	10	requirement.....	18
	Setting the language.....	10	Data backup of the log file.....	19
	Setting the date/time.....	11	Deletion of the log file on the STATMATIC smart.....	19
	Assemble and check the adapter.....	11	Emptying the media system.....	19
	Functional testing of flow rate using the testing		Warning message: Servicing required soon.....	20
	adapter.....	11	Error message: Servicing required.....	20
3	Operation.....	12	Transport.....	20
	Manual pre-cleaning of instruments.....	12	5 Troubleshooting.....	21
	"Cleaning + Maintenance" program - for regular and		Replacing fuses.....	22
	contra-angle handpieces and turbines.....	12	6 Information on Electromagnetic Compatibility (EMC).....	23
	"Collet care" program.....	12	7 Accessories.....	25
	"Maintenance only" program - for regular and		8 Technical Data.....	26
	contra-angle handpieces and turbines.....	13	Type plate.....	27
	"Clean only" program - for regular and contra-angle		9 Documentation, Log File.....	28
	handpieces and turbines.....	13	10 Guarantee.....	29

STATMATIC PLUS is a trademark; OPTIM and STATIM are registered trademarks of SciCan Ltd. All other trademarks referred to in this manual are the property of their respective owners.

SciCan GmbH

Wangener Strasse 78
88299 Leutkirch, GERMANY
Phone: +49-7561-98343-0
Fax: +49-7561-98343-699

SciCan Ltd.

1440 Don Mills Road
Toronto, ON M3B 3P9, CANADA
Phone: +1-416-445-1600
Fax: +1-416-445-2727
Toll free: +1-800-667-7733
E-mail: techservice.ca@scican.com

SciCan Medtech AG

Alpenstrasse 14
6300 Zug, SWITZERLAND
Phone: +41-41-727-70-27
Fax: +41-41-727-70-29

SciCan Inc.

701 Technology Drive
Canonsburg, PA 15317, USA
Phone: +1-724-820-1600
Fax: +1-724-820-1479

Manufactured by:

SycoTec GmbH & Co. KG

Wangener Strasse 78
88299 Leutkirch
GERMANY
Phone: +49 7561 86-0
Fax: +49 7561 86-266
E-mail: info@sycotec.eu

1 Important Information

Congratulations on purchasing a *STATMATIC* smart. We are convinced that you have acquired the best product of its kind. The *STATMATIC* smart is a compact cleaning and maintenance system that can satisfy all your needs and expectations regarding the cleaning and maintenance of instruments thanks to its diverse programs and capabilities. This instruction manual contains all details on the installation, operation and maintenance of your *STATMATIC* smart system. In order to guarantee many years of safe and problem-free operation, you should read the instructions prior to commissioning the device and store them for future reference.

The instructions for the operation, maintenance and replacement of replacement parts must be followed so that the device can function properly. The content of this instruction manual may change without notice in order to reflect changes and improvements to the *STATMATIC* smart system.

Type description

The *STATMATIC* smart is manufactured in 3 different versions.

STATMATIC smart Type 1 is equipped with a DIN EN ISO 3964 coupling, which is operated by a motor as well as 2 neutral couplings that can be individually equipped with the required turbine adapters (see "Accessories" section). *STATMATIC* smart Type 2 is equipped with two DIN EN ISO 3964 couplings, which are operated by a motor as well as 1 neutral coupling that can be individually equipped with the required turbine adapter (see "Accessories" section). *STATMATIC* smart Type 3 is equipped with three DIN EN ISO 3964 couplings, which are operated by a motor.

Scope of delivery

Base *STATMATIC* smart unit
Power supply connection cable
Instructions for use
Absorption pad (set)
Compressed air hose
SciCan "Dental Oil" (250 ml)
SciCan "Cleaner" (500 ml)
O-Ring Set
2 seal plugs (red) for media container
USB stick (2 GB)

General

The *STATMATIC* smart has been developed and manufactured in accordance with the most recent quality and safety standards.

However, some precautions should be taken during installation and operation.

- Disconnect the *STATMATIC* smart system from the power supply when leaving the practice for a longer period of time.
- Protect the *STATMATIC* smart from temperatures above 50°C.
- Any changes or modifications made by a third party will invalidate the warranty as well as operating certification.
- Only use original replacement parts for operation or repair, as these have been tested extensively for safety, function and specific suitability. SciCan therefore also assumes any product responsibility for these parts and components.
- Observe all warnings on the device and in the instructions for use.
- The device may only be supplied with the voltage indicated on the type plate.
- The power cable must not be replaced by a weak or underpowered power cable.
- Operators are not allowed to perform maintenance on the unit in other respects than those stated in the instruction manual.
- Medical gloves have to be worn during commissioning, operation and maintenance of the device.
- The device may only be operated with original SciCan "Cleaner", SciCan "Dental Oil" and accessories.
- Prior to transport or before prolonged periods of standstill (> 30 days), the containers for cleaner and dental oil must be emptied.
- The operator must make sure that the device is in the proper operating condition prior to using the device as well as confirm its functional safety.
- A damaged device or a device with a damaged power cable has to be disconnected from the mains immediately. Operating the device is not permitted in this case.
- The device may only be operated inside buildings and areas that are not subject to explosive hazards.
- Never open the protective housing of the device or insert objects through the openings of the protective housing.

1 Important Information

Intended use

- The *STATMATIC smart* is a cleaning and maintenance system used for the machine-based treatment (internal cleaning and maintenance) of dental turbines as well as regular and contra-angle handpieces.
- Refer to the manufacturer instructions for further information on the suitability of certain instruments for cleaning and maintenance purposes.

! The *STATMATIC smart* treatment is restricted to the internal cleaning of the instruments. The outer surfaces of the instruments have to be cleaned manually or using machines prior to any internal cleaning. Following treatment with the cleaning and maintenance system, continue treating the instrument immediately afterwards and in accordance with the manufacturer's instructions.

NOTE: For all preceding and subsequent processing steps, the RKI guidelines as well as the Medical Products Law must be observed (applicable only for Germany).

Misuse

All applications not listed as part of the intend use constitute a case of misuse and may lead to damage to the device and the medical products to be reprocessed and/or a potential danger to the operator himself. SciCan shall accept no liability in such an event.

Disposal

Consumables

Any waste resulting from use of the device have to be disposed of or recycled without presenting a danger to human beings or the environment while existing national regulations must also be observed in this regard.

Devices and accessories after use

Based on EU Directive (WEEE 2012/19/EU) on waste electrical and electronic equipment we hereby inform you that this product is subject to the aforementioned directive and must be disposed of through special channels.

Contact the manufacturer or dental supply shop for further information.

This product must be disposed of or recycled without presenting a danger to human beings or the environment while existing national regulations must also be observed in this regard.

1 Important Information

Qualifications of the operator

The STATMATIC smart system may only be used by trained personnel, including dentists and their assistants who have been instructed in the use of this product.

Safety Notes SciCan "Dental Oil" and SciCan "Cleaner"

First-Aid-Measures - Dental oil

 General information	▶ Never give any oral medication to unconscious persons or persons suffering from cramps. Move the person to safety
When inhaled	▶ Bring the person to the fresh air and provide unobstructed breathing
Contact with skin	▶ Wash with plenty of water. Immediately changed out of contaminated clothes. Medical treatment required
Contact with eyes	▶ Immediately, carefully and thoroughly rinse with eye-bath or water. Consult an eye specialist in case of eye irritation.
When swallowed	▶ Vomiting may pose a danger of aspiration. Immediately rinse the mouth and drink plenty of water. Induce vomiting if the person is conscious. Medical treatment required.

Fire-fighting measures

The product is not flammable.

First-Aid-Measures - Cleaner

 General information	▶ Never give any oral medication to unconscious persons or persons suffering from cramps. Move the person to safety
When inhaled	▶ Provide fresh air.
Contact with skin	▶ Wash with plenty of water. Change out of contaminated clothing.
Contact with eyes	▶ Immediately, carefully and thoroughly rinse with eye-bath or water.
When swallowed	▶ Immediately rinse the mouth and drink plenty of water.

Fire-fighting measures


The product is not flammable.

Used symbols

In the manual:

	A potential danger to the operator
	A situation that may lead to a mechanical defect
	Important information
	Instructions for disposal
	Wear medical gloves

On the device:

	LED ready to operate		CE marking
	LED error		Fuse
	Manufacturer		Follow instruction manual
	VDE certification mark		MET certification mark
	Program selection / selection up		Menu / select / save
	Start program / selection down		Collet care / back
	USB port		

1 Important Information

Disclaimer


The installation of replacement parts, maintenance and upkeep may only be carried out by authorized personnel. SciCan assumes no liability whatsoever for coincidental, special or consequential damage resulting from upkeep or maintenance of the *STATMATIC* smart carried out by third parties or for the use of equipment or individual parts manufactured by third parties, including loss of profits, commercial losses of any kind, financial losses or losses due to the injury of persons.

Never remove the rear panel of the device and never insert any objects through the holes or openings on the housing. Noncompliance may result in damage to the device and/or danger to the operator.


1 Important Information

Device overview


- (1) Cover
- (2) Membrane keypad
- (3) Plug-in position for turbine adapter (2x Type 1 / 1x Type 2) or INTRAmatic ISO 3964 maintenance couplings (3x Type 3 / 2x Type 2 / 1x Type 1)
- (4) Collecting tray
- (5) Flap
- (6) Absorption pad (2x)
- (7) Collet care nozzle
- (8) Display

- (9) USB stick
- (10) Filling opening for dental oil
- (11) Filling opening for cleaner
- (12) Holding bracket for testing adapter
- (21) USB port
- (22) Media aeration


- (13) Power switch
- (14) Mains fuse drawer
- (15) Power line input
- (16) Drain hose for dental oil and cleaner
- (17) Pressure indicator
- (18) Pressure setting
- (19) Compressed air input
- (20) Type plate

1 Important Information

Membrane keypad

(L1) LED stand-by

(L2) LED error

(S1) Program selection / selection up

(S2) **Start program** / selection down

(D1) Display

(S3) Collet care / back


(S4) Menu / select / save

Display

Date / Time

Status lines

Slots

1

Slot number

2

Slot is finished

3

Slot is occupied

3

Frame is flashing:
Slot is being processed

2 Commissioning

Environmental conditions

There are a number of factors that can influence the performance of your STATMATIC smart system. Please consider these factors carefully and choose a suitable location for installing the device.

- **Temperature and air humidity**

Do not set up your STATMATIC smart system at a location where it is exposed to direct sunlight or a heat source (for example a heater). The recommended operating temperature lies between 5 and 35°C with air humidity between 25 and 85%.

- **Environment**

The STATMATIC smart should be operated in a clean, dust-free and well-ventilated environment.

- **Work surface**

The STATMATIC smart should be installed on a flat, level and water-resistant as well as oil-resistant surface. Do not install or operate your unit on a sloped surface.

- **Electromagnetic environment**

The STATMATIC smart has been tested for electromagnetic emissions and meets the decisive standards. While the product itself emits no radiation, it can be influenced by other radiating devices. We therefore recommend that you keep the device away from potential sources of interference.

Electrical Requirements

DANGER! Mortal danger with improper handling of power sources! Danger of electrical shock

- For connection to the mains, use a properly grounded and secure socket providing the stated voltage (see type plate on the back of the STATMATIC smart).
- Avoid using multiple sockets. When using power bars with surge protection, always connect only a single device.
- The mains input connector on the back of the device is to be used as a separator from the power supply.

Media connections

- ▶ Plug the mains plug into the socket.
- ▶ Attach the compressed air hose to the coupling and tighten.
The set pressure should be between **4.5 bar** and 4.3 bar!
In the event of a deviating display, pull out the grip for regulating pressure and turn it to the right or the left until the manometer shows 4.5 bar.
- ▶ Push the grip back into the pressure regulator and fixate it.
In order to ensure a safe and trouble-free operation, only use compressed air in accordance with the information in the "Technical Data" section.

2 Commissioning

Commissioning

- ▶ Turn on the mains switch.
- ▶ Select the desired language with the arrow keys next to the display (Page 10).
- ▶ Set the date and time (Page 11).
- ▶ "WARNING Fill with dental oil" and/or "WARNING Fill with cleaner" will be displayed.
- ▶ Pull off top cover with both hands approx. 15 mm (overcome magnetic resistance) and pull off upwards.
- ▶ Refill the media containers (also see chapter 4 "Maintenance").

Observe the color code on the media container! Blue = dental oil / green = cleaner
Only use original SciCan "Dental Oil" and SciCan "Cleaner".

Use of other media may put patients and users at risk and damage the device and the instruments.


- ▶ Turn the metal flaps to the back in order to open the media container.
- ▶ Fill up the container to the filling level marking. Upon exceeding the minimum filling level, a 1-second beep will sound. If the container is full, a 2-second beep will sound 2 times. The beep is always followed by a visual sign on the container.
- ▶ Close the media container.
- ▶ Attach the top cover and close until it stops.
- ▶ Turn the device off using the power switch and then back on again.
- ▶ The message "Fill media system" will be displayed after filling the dental oil and cleaner containers. The device will automatically start this function and run for about 3 minutes.

Filling the media system

This is only necessary if the function did not start automatically or if a container has been completely empty for more than 3 days. The containers must be filled already for this function.

- ▶ Turn on the device
- ▶ Press the  button. "Export log data" will appear on the display.
- ▶ Press the  button once until "Fill media system" is highlighted in black on the display.
- ▶ Using the  button, the "Fill media system" menu is selected.
The media system will be filled automatically. During this time, a clock will run backwards until 00:00.
- ▶ Press the  button to return to the "Cleaning + Maintenance" main menu

Setting the language

- ▶ Turn on the device
- ▶ Press the  button. "Export log data" will appear on the display.
- ▶ Press the  button 4x until "Language setting" is highlighted in black on the display.
- ▶ The "Language setting" menu is selected using the  button.
- ▶ The desired language is selected using the  button.
- ▶ Save the language selection using the  button.
- ▶ Return to the "Language setting" menu using the  button.
- ▶ Press the  button to return to the "Cleaning + Maintenance" main menu.

2 Commissioning

Setting the date/time


- ▶ Turn on the device
- ▶ Press the  button. "Export log data" will appear on the display.
- ▶ Press the  button 3x until "Date/time" is highlighted in black on the display.
- ▶ The "Date/time" menu is selected using the  button.
- ▶ Set the desired value using the  and  buttons
- ▶ Using the  button, the cursor is moved to the next possible setting. The values are set in the following order: Year, month, day, hour, minute.
- ▶ After final actuation using the , the program jumps back again.
- ▶ Press the  button to return to the "Cleaning + Maintenance" main menu.

Assemble and check the adapter

- ▶ Assemble the adapter, if required
- ▶ Check the MULTIflex receptacle (and/or other turbine adapter) for firm seating.
 - ➔ The clamping nut must be tight.
 - Retighten clamping nut, if required.

! Please note that only white *STATMATIC* smart / *STATMATIC* PLUS adapters can be used. Adapters (black) of the *STATMATIC* maintenance unit are not suitable.


- ▶ Check the output quantity on all cleaning slots using a testing adapter intended for the coupling.

Functional testing of flow rate using the testing adapter

- ▶ Attach the matching testing adapter to the corresponding instrument couplings (such as INTRAmatic, DIN ISO EN coupling).
- ▶ Close the front cover.
- ▶ Press the START button.
- ▶ After completion of the program, check the filling quantity in each of the testing adapters.
 - ➔ The fluid must be above the mark ("Min.") !
 - If the filling quantity is below the mark, please consult the SciCan customer service.

! Please note that only *STATMATIC* smart testing adapters can be used. *STATMATIC* PLUS testing devices are not suitable.

After this process has completed, the *STATMATIC* smart system is operational.

3 Operation




Manual pre-cleaning of instruments

- Clean the instrument primarily in the head area with a medium-hard toothbrush for 20 sec. under running water (non strong water jet) until the instrument looks fully clean. (drinking water, temperature 30 °C ± 5 °C).
- Leave the instrument for 5 minutes.

"Cleaning + Maintenance" program - for regular and contra-angle handpieces and turbines

- ▶ Turn on the mains switch.
LED  (Stand-by) lights up green. "Cleaning + Maintenance" will appear on the display.
- ▶ Open the front door.
Place regular and contra-angle handpieces and turbines on the adapters provided for this purpose (instruments must snap into place). A frame will appear on the display around the station number per occupied slot.
- ▶ Close the front cover.
- ▶ Push the  (Program Start) button on the control membrane keyboard.
Program runs automatically. During the "Cleaning + Maintenance" program, a clock will count backwards until 00:00 on the display.
All attached regular and contra-angle handpieces and turbines are cleaned and maintained for inside in turn. After a successful process, the program slots will be displayed inversely. The display shows: "Cleaning + Maintenance successful".
- ▶ Open the front door.
- ▶ Remove regular and contra-angle handpieces and turbines by pressing the locking button (A). Pull off turbines downwards (press support lever down).
The frame around the slot number will disappear from the display.


Wipe off the instruments using a damp cloth soaked with disinfectant. For Germany we recommend using disinfectants listed by the VAH (such as SciCan OPT/IM BLUE).

- Immediately following this treatment, the instrument must be treated further in accordance with the manufacturer's instructions on reprocessing (such as disinfection/sterilization in the STAT/IM).
- Always perform a function test before using the instruments on patients. Some liquid may come out.

"Collet care" program

- ▶ Turn on the mains switch.
LED  (Stand-by) lights up green. "Cleaning + Maintenance" will appear on the display.
- ▶ Open the front door.
- ▶ Push instrument with chuck opening into the collet care nozzle.
- ▶ Activate  (Collet care) key on the control membrane keyboard.
Program runs automatically. A clock on the display runs backwards until 00:00
The following is then shown on the display: " Collet care measure successful".
- ▶ Remove instrument from the nozzle.


- Wipe off the instruments using a damp cloth soaked with disinfectant. For Germany we recommend using disinfectants listed by the VAH (such as SciCan OPT/IM BLUE).
- Immediately following this treatment, the instrument must be treated further in accordance with the manufacturer's instructions on reprocessing (such as disinfection/sterilization in the STAT/IM).
- Always perform a function test before using the instruments on patients. Some liquid may come out.


3 Operation

The processes "Clean only" and "Maintenance only" can also be carried out individually.

"Maintenance only" program - for regular and contra-angle handpieces and turbines


- ▶ Turn on the mains switch.
LED  (Stand-by) lights up green. "Cleaning + Maintenance" will appear on the display.
- ▶ Open the front door.
Place regular and contra-angle handpieces and turbines on the adapters provided for this purpose (instruments must snap into place). A frame will appear on the display around the station number per occupied slot.
- ▶ Close the front door.
 - ▶ Press the  button 1x. "Maintenance only" will appear on the display.
 - ▶ Push the  (Program Start) button on the control membrane keyboard.
Program runs automatically. During the maintenance program, a clock will count backwards until 00:00 on the display.
All attached regular and contra-angle handpieces and turbines are maintained for inside in turn. After a successful process, the program slots will be displayed inversely. The display shows: "Maintenance successful".
- ▶ Open the front door.
- ▶ Remove regular and contra-angle handpieces and turbines by pressing the locking button.
Pull off turbines downwards (press support lever).
The frame around the maintenance slot number will disappear from the display.

- ! • Wipe off the instruments using a damp cloth soaked with disinfectant. For Germany we recommend using disinfectants listed by the VAH (such as SciCan OPT/IM BLUE).
- • Immediately following this treatment, the instrument must be treated further in accordance with the manufacturer's instructions on reprocessing (such as disinfection/sterilization in the STAT/IM).
- Always perform a function test before using the instruments on patients. Some liquid may come out.


"Clean only" program - for regular and contra-angle handpieces and turbines

- ▶ Turn on the mains switch.
LED  (Stand-by) lights up green. "Cleaning + Maintenance" will appear on the display.
- ▶ Open the front door.
Place regular and contra-angle handpieces and turbines on the adapters provided for this purpose (instruments must snap into place). A frame will appear on the display around the station number per occupied slot.
- ▶ Close the front cover.
- ▶ Press the  button 2x. "Clean only" will appear on the display.
- ▶ Push the  (Program Start) button on the control membrane keyboard.
Program runs automatically. During the cleaning program, a clock will count backwards until 00:00 on the display.

3 Operation

All attached regular and contra-angle handpieces and turbines are cleaned inside in turn. After a successful process, the program slots will be displayed inversely. The display shows: "Cleaning successful".

- ▶ Open the front door.
 - ▶ Remove regular and contra-angle handpieces by pressing the locking button. Pull off turbines downwards (press support lever). The frame around the maintenance slot number will disappear from the display.
- !
- Wipe off the instruments using a damp cloth soaked with disinfectant. For Germany we recommend using disinfectants listed by the VAH (such as SciCan OPT/M BLUE).
 - Immediately following this treatment, the instrument must be treated further in accordance with the manufacturer's instructions on reprocessing (such as disinfection/sterilization in the STAT/M).
 - Always perform a function test before using the instruments on patients. Some liquid may come out.


3 Operation

"Special instrument" program

Special instruments can be maintained and cleaned using this program (such as prophylaxis). The times require for maintenance and cleaning can be individually adjusted by a service technician.

- ▶ Turn on the mains switch.
LED  (Stand-by) lights up green. "Cleaning + Maintenance" will appear on the display
- ▶ Open the front door.
Place the instruments on the adapters provided for this purpose (instruments must snap into place). A frame will appear on the display around the station number per occupied slot.
- ▶ Close the front cover.
- ▶ Press the  button 3x. "Special instrument" will appear on the display.
- ▶ Push the  (Program Start) button on the control membrane keyboard.
Program runs automatically. During the "Cleaning" program, a clock will count backwards until 00:00 on the display
All attached instruments are cleaned and maintained for inside in turn. After a successful process, the program slots will be displayed inversely. The display shows: "Special instrument successful".
- ▶ Open the front door.
- ▶ Remove regular and contra-angle handpieces by pressing the locking button.
Pull off turbines downwards (press support lever).
The frame around the slot number will disappear from the display.

- !
 - Wipe off the instruments using a damp cloth soaked with disinfectant. For Germany we recommend using disinfectants listed by the VAH (such as SciCan OPT/M/BLUE).
 - Immediately following this treatment, the instrument must be treated further in accordance with the manufacturer's instructions on reprocessing (such as disinfection/sterilization in the STAT/M).
 - Always perform a function test before using the instruments on patients. Some liquid may come out.


4 Maintenance

Checking adapter fixation

Recommended:

- After installation
 - During operation: 1x week
- ▶ Check the MULTiflex receptacle (and/or other turbine adapter) for firm seating.
 - ➔ The clamping nut must be tight.
 - Retighten clamping nut, if required.
 - ▶ Check the output quantity on all cleaning slots using a testing adapter intended for the coupling.

! Please note that only *STATMATIC* smart testing adapters can be used. *STATMATIC PLUS* testing devices are not suitable.

Type: 2

Checking the O-rings

Recommended:


- During operation: 1x week
- ▶ Check all O-rings on the adapters for damage
 - ▶ Replace damaged O-rings

Functional testing of flow rate using the testing adapter

Recommended:

- Upon initial commissioning
 - After changing an instrument coupling
 - During operation: 1x day before use
- ▶ Attach the matching testing adapter to the corresponding instrument couplings (such as INTRAmatic, DIN ISO EN Coupling).
 - ▶ Close the front cover.
 - ▶ Press the START button.
 - ▶ After completion of the program, check the filling quantity in each of the testing adapters.
 - ➔ The fluid must be above the mark ("Min.") !
 - If the filling quantity is below the mark, please consult the SciCan customer service.

! Please note that only *STATMATIC* smart testing adapters can be used. *STATMATIC PLUS* testing devices are not suitable.

Type: 3

Pressure test:

Recommended:

- Upon initial commissioning
 - during operation: 1x day, preferably before starting work
- Read off the "operating pressure" using the pressure gage.
- ➔ Pressure should be between **4.5 bar** and 4.3 bar!

Changing the adapter

- ▶ The loosening of the clamping nut allows for the exchangeable adapter to be removed. Insert the desired adapter and tighten it again using the clamping nut.

! Please note that only white *STATMATIC* smart / *STATMATIC PLUS* adapters can be used. Adapters (black) of the *STATMATIC* maintenance unit are not suitable.

- ▶ Check the output quantity on the adapter using a testing adapter intended for the coupling.


! Please note that only *STATMATIC* smart testing adapters can be used. *STATMATIC PLUS* testing devices are not suitable.

4 Maintenance

Refilling of the media containers

Check the filling levels of the containers every day prior to commencing work



- ▶ Turn on the mains switch
- ▶ The following will appear on the display: "WARNING Fill Dental oil" or "WARNING Fill Cleaner"
- ▶ Pull off top cover with both hands approx. 15 mm (overcome magnetic resistance) and pull off upwards.

! Observe the color code on the media container!

Blue = dental oil (SciCan "Dental Oil")
Green = cleaner (SciCan "Cleaner")

- ▶ Turn the metal flaps to the back in order to open the media container.


- ▶ Fill up the container to the filling level marking.



Upon exceeding the minimum filling level, a 1-second beep will sound.

If the container is full, a 2-second beep will sound 2 times.

The beep is always followed by a visual sign on the container.

- ▶ Attach the top cover and close until it stops.

! Only use original SciCan "Dental Oil" and SciCan "Cleaner"!

Use of other media may put patients and users at risk and damage the device and the instruments.


4 Maintenance

Replacing the absorption pads

Replace both pads after 35 instrument cleaning processes or at least 1x week (more often if needed).

- ▶ Close the front door.
- ▶ Pull the pads from the pins.
- ▶ Remove pads upwards from fastening clips.
- ▶ Dispose of pads in household waste.
- ▶ Hook new fleeces into the fastening clips.
- ▶ Fixate fleeces on the pins.
- ▶ Close the front cover.


Only use original SciCan absorption fleeces.

Empty the collecting tray

Empty the collecting tray every day and clean/disinfect it

- ▶ Open the front door.
- ▶ Pull out tray forwards (horizontal) (overcome slight catch).
Empty the tray and clean/disinfect it with a disinfectant wipe.
For Germany we recommend using disinfectants listed by the VAH (such as SciCan OPT/IM BLUE).
- ▶ Push tray into the side guides until it locks.

The rear fleece must not be pinched in, but rather has to hang loosely inside the tray.

- ▶ Close the front cover.

Inside cleaning and disinfection of the maintenance area

Clean the maintenance area and door with disinfectant wipes every day. For Germany we recommend using disinfectants listed by the VAH (such as SciCan OPT/IM BLUE).

Outside cleaning and outside disinfection of the device

The surface and control panel of the STATMATIC smart system can be cleaned and disinfected using disinfectant wipes. For Germany we recommend using disinfectants listed by the VAH (such as SciCan OPT/IM BLUE).

Recommendation

SciCan recommends an annual maintenance measure to be performed by an authorized dealer or the manufacturer.

NOTE: For users in countries with a legal validation requirement

- An annual maintenance by the manufacturer or an authorized specialist company is recommended.
- Please observe the national specifications for validation. An authorized specialist company must carry out a first validation with performance qualification (PQ) in accordance with the legal requirements if validation is required. Please observe the validation instructions for the device.
- SciCan recommends that a revalidation be carried out at intervals of one year (acc. to EN ISO 15883-1), unless shorter cycles are required by local authorities or the validator. These revalidations should be carried out during or following maintenance work.


4 Maintenance

Data backup of the log file

- ▶ Switch off the device
- ▶ Insert the USB stick on the back of the device
- ▶ Turn on the device
- ▶ Press the  button. "Export log data" will appear on the display.
- ▶ Use the  button to start exporting the log files to the USB stick.

When the process is finished, the display shows the following: Log data export complete.


- ▶ Use the  button to return to the menu.
- ▶ Switch off the device
- ▶ Remove the USB stick.

- ! • If the USB port is to be extended using a USB cable, this cable must be equipped with a ferrite core. The cable must be inserted with the ferrite core facing the device side.
- ! • Please note that the data on the USB stick must be copied to a suitable permanent storage device. Please check that the data on the storage medium can also be read correctly.

Deletion of the log file on the STATMATIC smart

! Deletion of the log file causes the data to be permanently lost. Please create a backup of the log file first and check whether the data on the USB stick can also be read correctly.

- ▶ Turn on the device
- ▶ Press the  button. "Export log data" will appear on the display.
- ▶ Press the  button 6x until "Delete log file" is highlighted in black on the display.
- ▶ "Deletion" is selected using the  button.
- ▶ Press the  button once more. The deletion of the log data will be performed.
- ▶ Use the  button to return to the menu.
- ▶ Press the  button to return back to the "Cleaning + Maintenance" main menu

Emptying the media system

This is only required if the system is to be emptied completely - for example for transport purposes or longer standstill periods (>30 days).

First, empty the containers completely using drain hoses. Then, tilt the device to the rear by approx. 45 degrees in order to also empty the residual liquids located at the bottom of the container using the drain hoses. Media can only be re-used again if they have been kept in their original bottles (mind the best before dates on the bottom or on the sides of the original bottles).

- ▶ Turn on the device
- ▶ Press the  button. "Export log data" will appear on the display.
- ▶ Press the  button 2x until "Empty media system" is highlighted in black on the display.
- ▶ Using the  button, the "Empty media system" menu is selected.
The media system will be automatically emptied. During this time, a clock will run backwards until 00:00.
- ▶ Press the  button to return to the "Cleaning + Maintenance" main menu

4 Maintenance

Warning message: Servicing required soon

This warning message appears if more than 25,000 regular and contra-angle handpieces and turbines have been cleaned and maintained for with the STATMATIC smart. Please contact the service department. You can clean a further 3,300 regular and contra-angle handpieces and turbines.

Error message: Servicing required


This error message appears if more than 28,300 regular and contra-angle handpieces and turbines have been cleaned and maintained for with the STATMATIC smart. Please contact the service department.

Transport

- If the device was already filled with dental oil and cleaner, the media must be emptied out fully through the drain hoses.
- Empty the media system (see p. 20)
- Unscrew the safety screw of the media aeration and blow compressed air into the valve several times. Tighten the screw again.
- Place the accompanying seal plugs in the filling funnel of the media container (see fig.)
- Close the flaps.
- If possible, send in the original box.

5 Troubleshooting


Error	Display	Acoustic signal	Cause	Correction
LEDs don't light up	No display	-	No power supply	Activate mains switch Plug in the power supply plug
			Defective fuse	Replace the fuse (see Chapter 6 "Troubleshooting")
Compressed air	ATTENTION Compressed air error Fault LED (L2) flashes	Interval – Beep 0.5 seconds on, 8 seconds off; 10 times	Compressed air hose not connected	Connect compressed air hose
			Compressed air hose pinched	Remove flow disruption
			Mains supply defective	Create supply
			Pressure reducer defective	Please consult the SciCan customer service.
			Pressure switch defective	Please consult the SciCan customer service.
Top cover open	WARNING Close the top cover Fault LED (L2) flashes	-	Top cover not closed correctly	Close top cover tightly
			Magnet on top cover is missing	Please consult the SciCan customer service.
Program doesn't start	ATTENTION Close the door ATTENTION No slots occupied	-	Front doors are not correctly closed	Close front doors tightly
			No instruments in the device	Attach the instruments to the couplings
Process interrupted	ATTENTION Process interrupted	-	Switching flap does not recognize instrument	Turn the instrument until a frame is displayed around the station number
			Front doors are not correctly closed	Close front doors tightly.
Undefined fault	ATTENTION Please consult the customer service.	Interval – Beep 2 seconds on, 4 seconds off	Instrument is not attached to the coupling correctly	Push the instrument firmly onto the coupling
			Internal device fault	Please consult the SciCan customer service.
Memory full	ATTENTION Not enough memory available	-	Internal memory full	Export log data to USB stick

5 Troubleshooting

Error	Display	Acoustic signal	Cause	Correction
Cleaner empty	WARNING Fill the cleaner Fault LED (L2) flashes	Beeping sound 2 seconds	Cleaner container empty	Fill the cleaner container with SciCan "Cleaner"
Oil empty	WARNING Fill in dental oil Fault LED (L2) flashes	Beeping sound 2 seconds	Oil container empty	Fill the oil container with SciCan "Dental Oil"
WARNING Fill Dental oil or Cleaner, even though the respective container is full.	WARNING Fill Dental oil or Cleaner Fault LED (L2) flashes	Beeping sound 2 seconds	Air in the fill measurement system	Fully empty out the affected container through the drain hose. Tilt the device back 45°. Unscrew the safety screw of the media aeration and blow compressed air into the valve several times. Tighten the screw again and fill the container.
Export data	ATTENTION Storage medium not readable Fault LED (L2) flashes	-	USB stick not recognized	Turn off the device, insert the USB stick, turn on the device
			USB stick not formatted correctly	- Perform NTFS formatting - Use another USB stick
Liquids (dental oil and cleaner) have been switched around	-	-		Empty media container: <ul style="list-style-type: none"> ▶ Drain liquids from the container using a drain hose. Draw any residual liquid from the corresponding media container using a 50 ml syringe. ▶ Completely remove residue of the "false "medium from the media container with a highly absorbent textile cloth that is free of fibers. ▶ Turn on the device Select "Empty media system" on the display (see Chapter 4 "Maintenance") ▶ Refill dental oil and cleaner (see Chapter 4 "Maintenance") <ul style="list-style-type: none"> - Open the front door - Select the "Fill media system" menu item in the display - Close the front door ▶ After this process has completed, the device is operational.

Replacing fuses

- ▶ Remove the STATMATIC smart device plug.
- ▶ Remove fuse tray from the power input module (push both tabs on the fuse tray inward).
- ▶ Replace fuses.
- ▶ Insert fuse tray again.

Only use the fuse type indicated in the technical data.

6 Information on Electromagnetic Compatibility (EMC)

The STATMATIC smart is designed for operation in an environment as described below. The customer or the user of the STATMATIC smart should ensure that it is operated in an environment of this kind.			
Electromagnetic interference measurements	Compliance	Electromagnetic environment – guidelines	
HF emissions according to CISPR 11	Group 1	The STATMATIC smart utilizes HF energy exclusively for its internal function. Its HF emissions are therefore very low and it is unlikely that neighboring electronic devices will be disrupted.	
HF emissions according to CISPR 11	Class B	The STATMATIC smart is suitable for usage in all facilities, including those located in living areas and those directly connected to a public power supply also supplying buildings used for residential purposes.	
Emissions from harmonics according to IEC 61000-3-2	Class A		
Emissions from voltage fluctuations / flickers according to IEC 61000-3-3	Conforms		
Susceptibility tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Discharge of static electricity (ESD) according to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Flooring should be of wood or concrete, or be covered with ceramic tiles. When the floor is covered with synthetic material, the relative air humidity must be at least 30 %.
Rapid, transient electrical disturbances/ bursts according to IEC 61000-4-4	± 2 kV for power cables ± 2 kV for input and output lines	± 2 kV, not applicable for power cables Not applicable	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Surges according to IEC 61000-4-5	± 1 kV normal mode voltage ± 2 kV common mode voltage	± 1 kV normal mode voltage ± 2 kV common mode voltage	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Voltage drops, short-term interruptions and fluctuations of the supply voltage according to IEC 61000-4-11	0 % UT for 0.5 periods 40 % UT for 5 periods 70 % UT for 25 periods 0 % UT for 250 periods	0 % UT for 0.5 periods 40 % UT for 5 periods 70 % UT for 25 periods 0 % UT for 250 periods	The quality of the supply voltage should correspond to that of a typical business or hospital environment. When the user of the STATMATIC smart also requires continuing function with the occurrence of interruptions of the power supply, it is recommended that the STATMATIC smart draw power from an uninterruptable power supply.
Magnetic field of the supply frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	Not applicable	Not applicable
Conducted HF disturbance according to IEC 61000-4-6	3 Veff 150 kHz to 80 MHz outside of the ISM bands ^a	10 Veff	Portable and mobile radio devices should be used at a distance from the STATMATIC, including the cable, no less than the recommended protective distance calculated according to the equation applicable to the transmission frequency. Recommended protective distance: $d = 1.17 \sqrt{P}$ $d = 0.35 \sqrt{P}$ for 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ for 800 MHz to 2.5 GHz with P as the maximum power rating of the transmitter in watts (W) in accordance with the information of the transmitter manufacturer and d as the recommended protective distance in meters (m). ^b The field strength of stationary radio transmitters should be less than the compliance level for all frequencies according to experiments carried out on location. ^d Interference is possible in the vicinity of devices bearing the icon shown below.
Radiated HF disturbances according to IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	10 V/m	
			

6 Information on Electromagnetic Compatibility (EMC)

NOTE U_T : is the alternating supply voltage prior to application of the test levels.

NOTE 1: The higher frequency range applies for 80 MHz and 800 MHz.

NOTE 2: These guidelines may not be applicable in all cases. The spread of electromagnetic disturbances is influenced by absorption and reflection of the building, objects and people.

^a The ISM frequency bands (for industrial, scientific and medical applications) between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz and 40.66 MHz to 40.70 MHz.

^b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range between 80 MHz and 2.5 GHz are intended to reduce the probability of mobile/portable communication devices causing any interference whenever they are inadvertently brought into patient areas. For this reason, the additional factor of 10/3 is used for the calculation of the recommended protective distances in these frequency ranges.

^c The field strengths of fixed transmitters, such as base stations for radio telephones and mobile radios, amateur radio stations, AM and FM radio broadcasts as well as TV broadcasts cannot be predicted in a fully accurate manner. In order to determine the electromagnetic environment with respect to the stationary transmitter, a study of the location should be considered. When the measured field strengths at the location at which the *STATMATIC* smart is used exceeds the compliance level above, the *STATMATIC* smart should be observed in order to demonstrate correct functioning. When unusual performance features are observed, additional measures may be necessary, such as a changed orientation or another location for the *STATMATIC* smart.

^d Above the frequency range of 150 kHz to 80 MHz, the field strength should be less than 3 V_{eff} V/m.

Recommended protective distances between portable and mobile HF telecommunications devices and the *STATMATIC* smart.

The *STATMATIC* smart is designed for operation in an electromagnetic environment in which the HF disturbances are controlled. The customer or the user of the *STATMATIC* smart can help avoid electromagnetic disturbances by observing the minimum distance between portable and mobile HF telecommunications devices (transmitters) and the *STATMATIC* smart, depending upon the output power of the communications device, as indicated below.

Rated power of the transmitter / W	Protective distance depending upon the transmission frequency in m		
	150 kHz to 80 MHz $d=0.35 \sqrt{P}$	80 MHz to 800 MHz $d=0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d=0.70 \sqrt{P}$
0.01	0.04	0.04	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.70
10	1.11	1.11	2.21
100	3.50	3.50	7.00

For transmitters with a maximum rated power not shown in the table above, the recommended protective distance d can be determined in meters (m) using the equation from the respective column, whereby P is the maximum rated power of the transmitter in watts (W) in accordance with the information provided by the transmitter manufacturer.

NOTE 1: The higher frequency range applies for 80 MHz and 800 MHz.

NOTE 2: These guidelines may not be applicable in all cases. The spread of electromagnetic disturbances is influenced by absorption and reflections of the building, objects and people.

7 Accessories

Please order accessories only from your authorized dealer.

Article no.	Description	Article no.	Description
 S32102	SciCan "Cleaner" (6x 500 ml)	 S30146	Chuck tip
 S32101	SciCan "Dental Oil" (2x 500 ml)	 S32117	Fuse, T 2.0 A (H)
 S32100	Absorption pad (set) (10 pieces)	 S30108	Rubber base
 S32114	Compressed air hose (line)	 S30122	Maintenance coupling gasket
 S30116	MULTiflex* O-Ring set	 S33120	Front door seal
 S32118	INTRAMatic O-Ring set	 S32201-E	Testing adapter INTRAMatic DIN EN ISO 3964- Coupling
 S32101-CX	Maintenance coupling MULTiflex*/ CONNEX	 S32201-CX	Testing adapter MULTiflex*
 S32101-ST1T	Maintenance coupling Sirona T1 turbine	 S32201-ST1T	Testing adapter Sirona T1 turbine
 S32101-NSK	Maintenance coupling NSK turbine	 S32201-NSK	Testing adapter NSK turbine
 S32101-W&H	Maintenance coupling W&H turbine	 S32201-W&H	Testing adapter W&H turbine
 S32101-BA	Maintenance coupling BienAir turbine	 S32201-BA	Testing adapter BienAir turbine
 S33126	Seal plugs	 S32202	Spare test tube STATMATIC smart

*MULTiflex is a registered trademark of Kaltenbach & Voigt

8 Technical Data

Device:

Device dimensions:	Width:	300 mm (11.8 inches)
	Depth:	280 mm (11.0 inches)
	Height:	345 mm (13.6 inches)
Weight:		7.3 kg (15.8 lbs)
Power ratings:	Rated voltage:	100 – 230 V \pm 10%
	Rated frequency:	50 / 60 Hz
	Rated power:	45 VA max.
Protection class:		I
Protection:		covered
Installation category:		Cat. II
Contamination level:		P 2
Ambient temperature:		5 – 35 °C (41 – 95 °F)
Noise level:		60 (65) dB(A)
Relative air humidity:		85 % max.
Max. operating altitude:		2000 m above sea level
Operating location:		only indoors
Compressed air:		4.5 – 6 bar (65 – 87 psi)
Air consumption:		40 NI/min
Compressed air quality:	Residual oil content:	max. 0.1 mg/m ³
	Residual dust:	max. 1 μ m; max. 1 mg/m ³
	Residual water:	max. 0.1 g/m ³ at -40 °C (-40 °F)
Power fuses:	Rated current:	2 A
	Rated voltage:	250 V
	Triggering characteristic:	Large (T)
	Breaking capacity:	H
Number of ISO 3964 integrated couplings:	Type 1	1
	Type 2	2
	Type 3	3

Device:

	Storage conditions	Transport conditions
Ambient temperature:	5 – 35 °C C (41 – 95 °F)	-30 – 50 °C (-22 – 122 °F)
Relative air humidity:	5 – 95 %	5 - 95 %
Compressed air:	700 - 1060 hPa	700 - 1060 hPa
Empty the media containers prior to transport! (see Chapter 4 "Maintenance")		


Media:

	Storage conditions	Transport conditions
Ambient temperature – dental oil:	5 – 50 °C (41 – 122 °F)	-30 – 50 °C (-22 – 122 °F)
Ambient temperature – cleaner:	5 – 35 °C (41 – 95 °F)	-30 – 50 °C (-22 – 122 °F)
Relative air humidity:	5 - 95 %	5 - 95 %
Compressed air:	800 - 2500 hPa	800 - 2500 hPa


8 Technical Data

Type plate


STATMATIC smart 1

STATMATIC smart 2

STATMATIC smart 3

The symbols used are described in Chapter 1 "Important instructions"

9 Documentation, Log File

STATMATIC smart : SN 000050 : Log-File																				
22.09.2014 15:56																				

Software-Version :2.00																				

Start time	Program	Intervals cleaning Spray	Intervals cleaning Drive	Intervals lubricating	Operating pressure	Start Place 1	Finish Place 1	Start Place 2	Finish Place 2	Start Place 3	Finish Place 3	Magnetic valve 10	Magnetic valve 11	Frontdoor Start	Frontdoor End	Finish time	Result	SYS-INFO	Valid	
07.08.2014 09:36	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	07.08.2014 09:36	OK	2	0	
07.08.2014 09:41	1	6	6	2	1	1	1	1	1	0	0	1	1	1	1	07.08.2014 09:48	OK	0	1	
07.08.2014 12:02	1	6	6	2	1	1	1	1	1	0	0	1	1	1	1	07.08.2014 12:10	OK	0	1	
07.08.2014 13:53	1	6	6	2	1	1	1	1	1	0	0	1	1	1	1	07.08.2014 14:00	OK	0	1	
07.08.2014 13:54	1	6	6	2	1	1	1	1	1	0	0	1	1	1	1	07.08.2014 14:01	OK	0	1	
07.08.2014 13:56	1	6	6	2	1	1	1	1	1	0	0	1	1	1	1	07.08.2014 14:04	OK	0	1	
07.08.2014 14:03	1	6	6	2	1	1	1	1	1	0	0	1	1	1	1	07.08.2014 14:10	OK	0	1	
07.08.2014 14:16	1	6	6	2	1	1	1	1	1	0	0	1	1	1	1	07.08.2014 14:23	OK	0	1	
11.08.2014 15:40	1	6	6	2	1	1	1	1	1	1	1	1	1	1	1	11.08.2014 15:50	OK	1	1	
13.08.2014 08:46	1	6	6	2	1	1	1	1	1	1	1	1	1	1	1	13.08.2014 08:56	OK	0	1	
14.08.2014 14:06	1	6	6	2	1	0	0	1	1	1	1	1	1	1	1	14.08.2014 14:13	OK	0	1	
14.08.2014 14:49	1	6	6	2	1	0	0	1	1	0	0	1	1	1	1	14.08.2014 14:53	OK	0	1	
14.08.2014 14:52	1	6	6	2	1	1	1	1	1	0	0	1	1	1	1	14.08.2014 14:59	OK	0	1	
22.09.2014 10:14	1	6	6	2	1	1	0	0	0	1	1	1	1	1	1	22.09.2014 10:21	Fail	0	0	
22.09.2014 10:16	1	6	6	2	1	0	0	0	0	1	1	1	1	1	1	22.09.2014 10:19	OK	0	1	
22.09.2014 10:58	1	6	6	2	1	0	0	0	0	1	1	1	1	1	1	22.09.2014 11:02	OK	0	1	
22.09.2014 12:03	1	6	6	2	1	1	1	1	1	1	1	1	1	1	1	22.09.2014 12:13	OK	0	1	
22.09.2014 12:04	2	0	0	2	1	0	0	0	0	1	1	1	1	1	1	22.09.2014 12:08	OK	0	0	
22.09.2014 12:05	3	6	6	0	1	0	0	0	0	1	1	1	1	1	1	22.09.2014 12:09	OK	0	0	
22.09.2014 12:07	5	0	0	2	1	0	0	0	0	0	0	1	1	0	0	22.09.2014 12:08	OK	0	0	
22.09.2014 12:07	4	1	1	1	1	0	0	1	1	0	0	1	1	1	1	22.09.2014 12:09	OK	0	0	
22.09.2014 14:47	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	22.09.2014 14:47	OK	2	0	
*Program 1 was started on 08/11/2014 at 3:40 p.m. with 3 instruments and successfully ended at 3:50 p.m. with tested operation data.																				
○ Fault in station 1; instrument slipped during operation from the coupling -> Fail																				
Title	= Explanations regarding the headline															Program no.				
Start time	= Starting time DD.MM.YYYY HH:MM															0	Log data deleted			
Program	= Program no. (see separate table)															1	Cleaning + Care			
Intervals cleaning Spray	= Number of cleaning intervals of the spray channel															2	Care only			
Intervals cleaning Drive	= Number of cleaning intervals of the drive channel															3	Clean only			
Intervals lubricating	= Number of lubrication intervals															4	Special instrument			
Operating pressure	1 = Compressed air available; 0 = No compressed air															5	Collet care (chuck)			
Start Place 1	1 = Slot 1 occupied upon start; 0 = Not occupied															6	Continuous operation (service)			
Finish Place 1	1 = Slot 1 occupied at the end; 0 = Not occupied															7	Filling the media system			
Start Place 2	1 = Slot 2 occupied upon start; 0 = Not occupied															8	Emptying the media system			
Finish Place 2	1 = Slot 2 occupied at the end; 0 = Not occupied																			
Start Place 3	1 = Slot 3 occupied upon start; 0 = Not occupied																			
Finish Place 3	1 = Slot 3 occupied at the end; 0 = Not occupied																			
Magnetic valve 10	1 = Magnet valve 10 (cleaner) activated; 0 = Not activated																			
Magnetic valve 11	1 = Magnet valve 11 (oil) activated; 0 = Not activated																			
Frontdoor Start	1 = Door closed at the end; 0 = Open																			
Frontdoor End	1 = Door closed at the start; 0 = Open																			
Finish time	= End time DD.MM.YYYY HH:MM																			
Result	Fail = faulty program execution; 0 = OK																			
SYS-INFO	= System information (see separate table)																			
Valid	1 = Program 1 with tested operating data completed successfully																			
																	6	Cleaner refilled		
																	7	Power supply interrupted		
																	8	Cleaner & oil empty		
																	9	Data backup is enforced		

10 Guarantee

Limited warranty

SciCan extends a guarantee of two years (from the date of purchase) against any malfunction of the **STATMATIC smart** due to material or processing errors not resulting from obvious misuse, false use or accidents, with the foregoing applying to devices in a new and unused condition under normal conditions of use.

The two-year guarantee applies to the functional capability of all parts of the device, provided that the device is used and maintained in accordance with the description included in the user manual.

In the event of any malfunction that can be traced back to such defects arising during the guarantee period, **SciCan** is obligated only to either carry out repairs at no charge or to replace any defective parts (with the exception of O-rings and seals) at no charge at our own discretion, on condition that we are notified in writing within 30 days of the occurrence of such a malfunction and that the defective parts are also returned freight paid.

This guarantee will be considered to have been agreed upon when the product is accompanied by the original invoice issued by our authorized **SciCan** dealer and when that invoice identifies the item with an item number and clearly indicates the purchase date. No other type of receipt is acceptable. After the expiration of two years, all of the guarantees provided by **SciCan**, as well as other obligations with respect to the quality of the product are considered to be conclusively fulfilled. Any liability is thus terminated and no legal proceedings may thereafter be initiated against **SciCan** on account of a breach of any such guarantee or obligation.

We exclude and reject any specific guarantee not indicated here and any implied guarantee or representation with reference to performance capability, as well as any legal remedy in the event of a contractual violation, which, apart from this provision, could arise due to implication, legal interpretation, commercial practice or custom of trade, including any implied guarantee with respect to saleability or suitability for a particular purpose with reference to one or all products manufactured by **SciCan**. To find out more about **SciCan** products and services, please visit our website at www.scican.com.

(DE = original)

STATMATIC smart 1, 2, 3

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DE

Brugsvejledning DK
Manuel d'instructions FR
Istruzioni per l'uso IT
Gebruiksaanwijzing NL
Instrucciones de uso ES

via: www.scican.uk.com/downloads
via: www.scican.uk.com/downloads
via: www.scican.uk.com/downloads
via: www.scican.uk.com/downloads
via: www.scican.uk.com/downloads

SciCan

Your Infection Control Specialist™

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