

Buffalo Dental Arbor Bands Operating and Instruction Guide

Silicone carbide abrasive bands designed for aggressive trimming and finishing of acrylics, silicone relines, and other dental materials.

-Available in coarse (60 grit), medium (80 grit), and fine (120 grit)

-No-lap seams for smooth operation and trouble-free service.


-Economical bulk packaging promotes convenient, disposable use.

Always wear approved safety glasses to avoid eye injuries when grinding.

Always use a dust mask to avoid inhalation of grinding dust.

Always grind in such a manner as to direct particles away from the body and face.

Always follow the manufacturer's instructions and safety procedures when using this product with power tools.



Item No.	Description	Item No.	Description	Item No.	Description
02721	3/4" Fine, Pkg. of 100	02722	3/4" Fine, Pkg. of 500	02723	3/4" Fine, Pkg. of 1000
02711	3/4" Medium, Pkg. of 100	02712	3/4" Medium, Pkg. of 500	02713	3/4" Medium, Pkg. of 1000
02701	3/4" Coarse, Pkg. of 100	02702	3/4" Coarse, Pkg. of 500	02703	3/4" Coarse, Pkg. of 1000
02751	1/2" Fine, Pkg. of 100	02752	1/2" Fine, Pkg. of 500	02753	1/2" Fine, Pkg. of 1000
02741	1/2" Medium, Pkg. of 100	02742	1/2" Medium, Pkg. of 500	02743	1/2" Medium, Pkg. of 1000
02731	1/2" Coarse, Pkg. of 100	02732	1/2" Coarse, Pkg. of 500	02733	1/2" Coarse, Pkg. of 1000
02771	3/8" Fine, Pkg. of 100	02772	3/8" Fine, Pkg. of 500	02773	3/8" Fine, Pkg. of 1000
02761	3/8" Medium, Pkg. of 100	02762	3/8" Medium, Pkg. of 500	02763	3/8" Medium, Pkg. of 1000
02811	1/4" Fine, Pkg. of 100	02812	1/4" Fine, Pkg. of 500	02813	1/4" Fine, Pkg. of 1000
02801	1/4" Medium, Pkg. of 100	02802	1/4" Medium, Pkg. of 500	02803	1/4" Medium, Pkg. of 1000
02791	1/4" Coarse, Pkg. of 100	02792	1/4" Coarse, Pkg. of 500	02793	1/4" Coarse, Pkg. of 1000

Buffalo Dental Manufacturing Co., Inc.: PO Box 678, 159 Lafayette Dr., Syosset, NY 11791-0678 USA
 Tel: 516.496.7200 or 800.828.0203 Fax: 516.496.7751
www.buffalodental.com Sales@buffalodental.com

102721 05/18