(f)((in)() video.wh.com Now at your dealer or wh.com



DIY-Service

Alegra Spare Parts

In addition to the W&H Service you also have the opportunity to maintain your Alegra turbines and contra-angle handpieces yourself and to carry out minor repairs quickly and professionally on-site. Regular on-site maintenance extends the service life of your turbines and contra-angle handpieces and reduces downtimes.

Do it yourself

Take a look at our Alegra spare parts overview and follow our video tutorials on our W&H Video Channel.



alegra



Coupling seal REF 01000700 REF 02207300



Short nozzle cleaner for turbine, contra-angle and straight handpieces REF 02015101





Rotor with hexagon wrench for TE-98 LQ, TE-98, TE-98 BC, TE-98 RM REF 06787500 for TE-97 LQ, TE-97, TE-97 BC, TE-97 RM REF 07234100 for TE-95 RM REF 07495400



Water filter with resuction stop for BC connections REF 07092500 for RM connections REF 07095500

PLEASE NOTE!

Tools: For your own safety and the safety of your patients, please always order the compatible tools to carry out the repair in good quality.

Original Spare parts:



O-ring set for Roto Quick couplings REF 07508900

Generator

for RQ-53

for RQ-54

REF 06840300

REF 06793000

Torque wrench for replacing the rotor REF 07508800

Push-button

REF 06641900

REF 07548000

for TE-95 RM, TE-95 VC,

TE-98, TE-98 RM, TE-98 BC

for TE-97, TE-97 RM, TE-97 BC



Tool kit for replacing of contra-angle heads REF 07550300



Contra-angle head for WE-56 REF 07549900 for WE-57 REF 07550000 for WE-66 REF 07550100



The medical device loses its CE conformity if repairs are carried out with a non-original spare part. The use of non-original spare parts does not meet the high W&H quality standards.

Safe use and a long product life

can no longer be guaranteed.