

Robust. Powerful. Efficient.

KaVo **expert** series

KaVo EXPERTtorque LUX E680 L

- **Max. working efficiency**, high torque and optimal power of up to 18 W
- **Optimal service life** due to its excellent head-protect system and stainless steel heads
- **Proven configuration efficiency** due to customer-replaceable water filter and robust chuck system

KaVo EXPERTtorque Mini LUX E677 L

- **Perfect vision** due to smaller head size and 100°/19° head/elbow angle combination
- **Efficiency** due to robust KaVo chuck system for standard and short-stem drills
- **Optimal service life** due to high chuck retention force of up to 30 N



KaVo EXPERTtorque
LUX E680 L

KaVo EXPERTtorque
Mini LUX E677 L



KaVo EXPERTtorque Turbines

KaVo EXPERTtorque LUX E680 L

Power:

- 18 watts power
- Stable speed

Improved vision:

- Head height reduced to 13.6 mm
- 100°/19° head/elbow angle combination (KaVo patent)
- Integrated light

Efficiency and infection protection:

- Four-nozzle spray
- Head-protect system: New improved back-suction effect while the turbine comes to a stop
- Chuck system: Low trigger force

Low noise level:

- Up to 62 dB(A)

Better hygiene conditions:

- Spray air/spray water guided on the interior
- New customer-replaceable water filter



KaVo EXPERTtorque Mini LUX E677 L

Improved vision:

- 12.2 mm head height (miniature head)
- 100°/19° head/elbow angle combination (KaVo patent)
- Integrated light

Power:

- 15 watts power
- Stable speed

Efficiency and infection protection:

- Three-nozzle spray
- Head-protect system: New improved back-suction effect while the turbine comes to a stop
- Chuck system: Low trigger force, use of standard and short-stem drills

Low noise level:

- Up to 65 dB(A)

Better hygiene conditions:

- Spray air/spray water guided on the interior
- New customer-replaceable water filter



Special features of the EXPERTtorque:

- EXPERTgrip surface for a safe grip
- Customer-replaceable water filter for consistently good spray
- Chuck retention force of up to 30 N
- Carbide guide sleeve for perfectly running drill
- Original KaVo ceramic ball bearing
- 18-month warranty
-  Made in Germany

Handpiece servicing:

KaVo QUATTROcare PLUS 2124 A

- Optimal servicing quality with QUATTROcare Plus Spray
- Simple and comfortable handling and rapid programme flow
- Low initial purchase and consumable costs



KaVo expert series

Article	Description	Material number
E680 L	EXPERTtorque LUX	1.006.8700
E680 C	EXPERTtorque (without light)	1.006.9000
E680 LS	EXPERTtorque LUX, for Sirona couplings	1.006.4600
E680 LN	EXPERTtorque LUX, for NSK couplings	1.006.4700
E677 L	EXPERTtorque Mini LUX	1.007.3600
E677 LN	EXPERTtorque Mini LUX, for NSK couplings	1.006.0100
MULTIflex 465 LED	MULTIflex LUX coupling 465 LED	1.007.3197
MULTIflex 465 RN	MULTIflex coupling 465 RN – without light	0.553.1600
QUATTROcare PLUS 2124 A	Handpiece servicing system	1.008.3805



KaVo. Dental Excellence.