

# MICRO-TOUCH® Nitrile Accelerator-Free

High Touch

# Nitrile, powder-free examination glove with advanced allergy prevention

This non-sterile accelerator-free offers a solution for those who suffer from chemical Type IV allergies. A comfortable fit and feel with micro-textured fingers ensure a secure and safe grip in wet or dry conditions.



- Provides a high-level of barrier protection against punctures as well as excellent resistance from a broad range of commonly used chemicals.
- · Chlorinated inner surface for easy donning
- Protects from latex Type I allergy
- 100% thiuram-free formulation significantly reduces the risk of Type IV allergies.

WARNING: No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.


## Recommended For

- · Tested for use with chemotherapy drugs
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients
- Standard patient examinations

## FEATURED TECHNOLOGIES







# MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.



# MICRO-TOUCH® Nitrile Accelerator-Free

High Touch

PRODUCT INFORMATION				
Material	Nitrile	Not Made From Natural Rubber Latex	Yes	
Color	Blue	Cuff Length	Standard	
Powder Content	Powder-Free	External Glove Surface	Textured Fingers	
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.08 / 3.1	
Finger Thickness (mm/mil)	0.12 / 4.7	Allergy Prevention	Latex (Type I), Chemical (Type IV)	
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)	
Sterile	No	Antistatic	Not Tested	
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	High Touch	

ORDERING INFORMATION						
Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	
Product Code	700101	700102	700103	700104	700105	

PACKAGING AND STORAGE		
Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case	
Shelf Life	3 Years	
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.	

## STANDARDS AND CERTIFICATIONS

CE 0493, Category III, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 9001, ISO 10993-10:2002, ISO 13485

#### **Vulcanization Chemical Accelerators**

**Accelerator Free** 

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America US Tel: 800 952 9916 CA Tel: 1-844-494-7854 Email: insidesalesus@ansell.com

Latin America & Caribbean Tel: +52(442) 296 20 50 Email: cslac@ansell.com Brazil Tel: +55-11-3356-3100 Brazil Email: Iuvas.medicas@ansell.com Europe, Middle East & Asia Pa Africa Tel: +6 Tel: +32 (0) 2 528 74 00 Email: Email: info@ansell.eu apac.m

Asia Pacific Tel: +603 8310 6688 Email: apac.medical@ansell.com

Australia & New Zealand Tel: +61 3 9270 7270 Email: protection@ap.ansell.com