Professional prophylaxis







Professional tooth cleaning is an effective method of preventing plaque deposits from forming, giving your patients a beaming smile.

With Proxeo polishing handpieces, W&H offers a compact, easy to handle solution, which makes it even easier to give your patients what they want.



Technical data

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.

Proxeo Universal Handpieces 18mm







Type:	WP-64 TU	HP-44 T
Max. motor speed:	10,000 rpm	10,000 rpm
Transition ratio:	4:1	4:1
Recommended speed for cups and brushes:	1,250 – 1,500 rpm	1,250 – 1,500 rpm
Maximum speed for cups and brushes:	2,500 rpm	2,500 rpm
Attachments that can be used:	All common screw-in cups and brushes (1-72 UNF 2A screw thread) as well as snap-on cups and brushes (adaptor enclosed)	Disposable contra-angle handpiece with plastic shaft, Doriot system
Coupling system:	For motors with ISO connection in accordance with ISO 3964/DIN 13940	For motors with ISO connection in accordance with ISO 3964/DIN 13940
Diameter of sheath:	18 mm	18 mm
Length:	89 mm	64 mm
Weight:	40 g	29 g
Sterilizable:	In the autoclave up to 135°C	In the autoclave up to 135°C

W&H Dentalwerk Bürmoos GmbH

Ignaz-Glaser-Straße 53, Postfach 1 5111 Bürmoos, **Austria** t + 43 6274 6236-0 f + 43 6274 6236-55 office@wh.com

wh.com